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
Pyro Light
Pyro Light XL

DE

EN

Enjoy mobility.

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Index

| | |
|--|-----------|
| 1. Preview | 53 |
| 1.1. Signs and Symbols | 54 |
| 2. Important safety instructions | 55 |
| 2.1. Common safety instructions | 55 |
| 2.2. Tipping hazard | 58 |
| 2.3. Safety in vehicles | 59 |
| 2.4. Participation in road traffic | 59 |
| 2.5. Brakes | 60 |
| 2.6. Indications | 60 |
| 2.7. Contraindications | 61 |
| 2.8. Declaration of conformity | 61 |
| 2.9. Responsibility | 61 |
| 2.10. Useful life | 61 |
| 3. Purpose | 62 |
| 4. Product and Delivery overview | 62 |
| 4.1. Checking the delivery | 62 |
| 4.2. Scope of delivery | 62 |
| 4.3. Overview | 63 |
| 4.4. Type plate and Serial number | 64 |
| 5. Assembly / Adaptations | 65 |
| 5.1. Setting up the wheelchair | 65 |
| 5.2. Folding of the wheelchair | 65 |
| 5.3. Foot plates | 66 |
| 5.4. Legrests | 66 |
| 5.4.1. Horizontal adjustable legrests (optional) | 66 |
| 5.5. Amputation leg support (optional) | 67 |
| 5.6. Tipping aid | 67 |
| 5.7. Side parts | 68 |
| 5.8. Armrests | 68 |
| 5.9. Seat | 69 |
| 5.9.1. Seat depth | 69 |
| 5.9.2. Seat height | 69 |
| 5.10. Seat angle | 69 |
| 5.11. Wheels | 70 |
| 5.11.1. Drive wheel height | 70 |
| 5.12. Height of the front wheel | 71 |

| | |
|--|-----------|
| 5.13. Wheelbase extension (optional) | 71 |
| 5.14. Parking brakes | 71 |
| 5.15. One-hand brake (optional) | 72 |
| 5.16. Operation of the parking brake with brake lever extension (optional) | 73 |
| 5.17. Drum brake (optional) | 74 |
| 5.18. Backrest | 75 |
| 5.18.1. Backrest cover | 75 |
| 5.19. Push handles | 75 |
| 5.20. Thru axles | 76 |
| 5.21. One-hand operation (optional) | 76 |
| 5.22. Passive lighting | 76 |
| 5.23. Accessories | 77 |
| 5.24. Anti-tip wheels | 77 |
| 5.25. Lap belt | 78 |
| 5.26. Therapy table (optional) | 78 |
| 5.27. Togo (optional) | 78 |
| 5.28. Cane holder (optional) | 79 |
| 5.29. Spoke protector (optional) | 79 |
| 5.30. Handrim cover (optional) | 79 |
| 5.31. Infusion / Oxygen bottle holder (optional) | 79 |
| 5.32. Carry-on bag (optional) | 80 |
| 5.33. Headrest (optional) | 80 |
| 6. Use of the wheelchair | 81 |
| 6.1. Check before setting off | 81 |
| 6.1.1. Checking the brakes | 81 |
| 6.1.2. Checking the tyres | 81 |
| 6.2. Side entry and exit | 81 |
| 6.3. Entering and exiting from front | 82 |
| 6.4. Overcoming stairs and high steps | 83 |
| 6.5. Driving with the wheelchair | 83 |
| 7. Technical data | 84 |
| 8. Care instructions for the user | 86 |
| 8.1. Cleaning and Disinfection | 86 |
| 8.2. Cleaning | 86 |
| 8.3. Disinfection | 87 |

| | | |
|-------|--|----|
| 9. | Maintenance instructions | 88 |
| 9.1. | Maintenance advice for the specialist dealer | 88 |
| 10. | Repair | 89 |
| 11. | Transport | 89 |
| 11.1. | Transport of the wheelchair | 89 |
| 11.2. | Transport in motor vehicles | 89 |
| 12. | Passing on the wheelchair | 90 |
| 13. | Storage / Shipping | 91 |
| 14. | Disposal | 91 |
| 15. | Warranty | 92 |

1. Preview

Dear user,
You have chosen a wheelchair from Bischoff & Bischoff. We thank you for your trust. The wheelchair has been designed to offer you many advantages and to meet your needs. The Pyro Light can be adjusted to your individual body dimensions due to the many adjustment options and the optimal adaptability with regard to seat height, seat angle and lower leg length.

The Pyro Light is designed for both indoor and outdoor use.

Please read and observe the operating instructions carefully before using your new wheelchair for the first time. This is an integral and necessary part of the wheelchair. Keep the operating manual handy and give it with the wheelchair when passing it on. For users with visual impairments, this document is accessible as a PDF file on our website www.bischoff-bischoff.com. Repair and adjustment work requires special technical training and may therefore only be carried out by Bischoff & Bischoff authorised dealers.

Notice!

Misprints, product, colour changes and technical specifications are reserved. Illustrations may contain accessories. The colours of the products shown may vary slightly in reality.

The latest version of the operating instructions can be found in the download area on our website: www.bischoff-bischoff.com

1.1. Signs and Symbols

Attention! Indicates particularly safety-relevant information. Observe operating instructions!

Important! Indicates particularly useful information in the respective context.

2. Important safety instructions

2.1. Common safety instructions

To avoid falls and dangerous situations, you should first practise using your new wheelchair on level, manageable terrain. An accompanying person is recommended in this case.

- Only drive on level, solid ground.
- Be aware of the risk of entrapment between moving parts when adjusting, using and maintaining the wheelchair.
- The accompanying person should be physically and mentally capable of driving a wheelchair.
- The driver and the accompanying person must not be under the influence of alcohol or drugs.
- Before using the wheelchair, please check that all attached parts are properly fastened.
- Only lift the wheelchair by firmly mounted parts (the armrests or the legrests).
- Make sure that the quick-release axles are correctly seated in the axle sockets, otherwise there is a risk of injury.
- Keep clothing and body parts away from gaps between moving parts!
- Please note that the risk of tipping over is increased if there is a shift in balance (e.g. due to strong upper body movements or when overcoming obstacles).
- When negotiating obstacles (steps etc.), an attendant must be consulted. Use the tipping aid (see chapter 5.24) to tilt the wheelchair accordingly (fig. A).
- Stairs and larger obstacles may only be climbed with the help of at least two accompanying persons. Use the lifting points for this (see chapter 6.4). If facilities such as ramps, lifts or stair-climbing aids are available, use them.

- When driving on a slope/gradient, it must be noted that depending on the set axle position, there is a risk of tipping over from just a few degrees. The use of an anti-tipper is recommended (Fig. A).
- Avoid collisions of any kind! Avoid driving into an obstacle (step, kerb) without braking or jumping down from steps (Fig. B).
- Replace the tyres as soon as the tread is worn down.
- When getting in or out of the wheelchair, apply both parking brakes.
- Use the parking brake only for parking, it must not be used to slow down the drive.
- Avoid driving on uneven, loose surfaces.
- In road traffic, the Road Traffic Act (StVZO) must be observed.
- Check the passive lighting (reflectors) of your wheelchair for integrity and visibility to other road users before each journey.
- The seat and upholstery materials meet the requirements of resistance to flammability according to EN 1021- 2, nevertheless do not expose upholstery or other parts to fire sources - such as cigarettes.
- The wheelchair may only be used to transport one person.
- The wheelchair must not be used to transport loads or more than one person.
- When braking from a fast speed or on long downhill stretches, using the hand rims, fingers and palms heat up. Caution, danger of burns! (Fig. C).
- We recommend using leather gloves for outdoor riding. These increase grip and at the same time protect your fingers and palms from dirt and injuries.
- Note that upholstered parts that are exposed to direct sunlight can heat up and cause skin injuries when touched. Therefore, cover these parts or protect the wheelchair from sunlight.

Protect the wheelchair from sunlight to avoid burns from heated components.

- Also note that the frame and upholstered parts can cool down considerably in winter. If possible, do not park the wheelchair outside in cold weather.
- Observe the storage and operating conditions.

Attention: Due to the dimensions of the wheelchair, it is necessary to check whether your wheelchair can also be moved in the planned environment. This applies in particular to the use of escape routes.

When transferring between wheelchair and bed, please always keep in mind:

- Apply both parking brakes.
- Swing the side part of the entry/exit side backwards.
- Swing the foot plates up or away (when getting in/out).

- Reporting of incidents:

If you, as the operator, user or their relatives, discover possible defects or functional limitations, please contact an authorised specialist dealer immediately. Operators, users or their relatives should inform their specialist dealer who provided them with the product of any suspected serious incidents that have had or may have had, directly or indirectly, any of the following consequences and that affect them.

You can also report this to the competent higher federal authority.

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of the state of health of a patient, user or other person,
- a serious threat to public health.

2.2. Tipping hazard

If available, the anti-tip wheels on both sides (see chapter 5.) prevent the wheelchair from tipping backwards as far as possible. Nevertheless, please note that driving with anti-tip wheels is only possible to a limited extent on stepped terrain, lifting platforms and ramps. Ensure sufficient distance upwards, downwards and to the sides.

2.3. Safety in vehicles

The Pyro Light XL is not suitable as a vehicle seat. It must not be used as such under any circumstances. It does not meet the requirements of ISO 7176-19.

Failure to do so can result in serious injury and death in the event of an accident (Fig. A).

The Pyro Light XL wheelchair is not crash-tested and must therefore never be used as a seat in motor vehicles for the transport of mobility-impaired persons (KMP)!

You can tell if your wheelchair is not approved as a seat in motor vehicles for the transport of mobility-impaired persons by the means of the symbol on the type plate, (chapter 4.4).

| | |
|----------------------|---|
| Pyro Light XL | Must not be used as a vehicle seat (Fig. A) |
| Pyro Light | May be used as a vehicle seat (Fig. B) |

2.4. Participation in road traffic

The wheelchair is designed for indoor and outdoor use. Please note that you are participating in public road traffic and must comply with the regulations of the Road Traffic Act. Do not endanger other participants by driving recklessly, this applies especially on pavements.

2.5. Brakes

Operating brake:

When driving, brake the wheelchair by using the handrims. Please note that your hands may get hot.

Parking brake/parking:

As soon as the wheelchair is stationary, push both levers of the knee lever brake forward (see chapter 5.14). The wheelchair is safely braked when it can no longer be pushed away and both levers are pushed forward as far as they will go.

Pull the levers back towards you to release the brake (see chap. 5.14). The locking levers located on the brake levers for the drum brake can also be used to lock the drum brake for parking (chap. 5.17).

Drum brake:

The wheelchair can be optionally equipped with a drum brake, which can only be operated by an attendant via the brake levers on the push handles.

The brake levers on the push handles only act on one drive wheel at a time! Always use both brake levers to brake the wheelchair.

2.6. Indications

Inability to walk or severe walking disability in adults due to

- Paralysis
- Limb loss
- Limb defect / deformity
- Joint contractures / joint damage (not on both arms)
- Other diseases

Use with reinforced wheelchairs is indicated when the body weight of the disabled person is 125 kg or more and a standard wheelchair no longer provides the necessary safety.

2.7. Contraindications

The use of the wheelchair by adults is unsuitable in the case of

- Perceptual disorders
- Severe balance disorders
- Loss of limbs on both arms
- Joint contractures / joint damage in both arms
- Incapacity to sit
- Reduced or insufficient vision

2.8. Declaration of conformity

We Bischoff & Bischoff GmbH declare under our sole responsibility that the Pyro Light / Pyro Light XL wheelchair complies with all requirements of MDR 2017/745 that are applicable.

2.9. Responsibility

We shall only provide a warranty if

- the product is used under the specified conditions and for the intended purpose,
- modifications, extensions, repairs and maintenance work are only carried out by persons who have been authorised by us to do so, and
- the wheelchair is used in compliance with all instructions for use.

2.10. Useful life

The expected service life is up to five years if used as intended and the safety, care and maintenance instructions are observed.

3. Purpose

The wheelchair is designed exclusively for mobility enhancement and the transport of adult persons with walking disabilities in accordance with stated indications. The maximum net weight is:

Pyro Light 125 kg

Pyro Light XL 170 kg

Please note:

We can only provide a warranty if the product is used under the specified conditions and for the intended purposes.

4. Product and Delivery overview

4.1. Checking the delivery

The Pyro Light wheelchair is delivered fully assembled in a cardboard box. After unpacking, please keep the carton if possible. It can be used for later storage or return of the wheelchair. The delivery and instruction of the wheelchair is carried out by a qualified and authorised medical supply retailer.

In the case of shipment by rail or forwarding agent, the goods must be checked immediately for transport damage in the presence of the deliverer.

Check the contents for completeness (see chapter 4.2) and intactness. In case of irregularities or damage, contact our customer service. The contact details can be found on the back cover of this manual.

4.2. Scope of delivery

The content consists of the following main components:

- Wheelchair, pre-assembled
- 1 pair of leg rests
- 1 Instruction manual
- 1 back cushion


4.3. Overview

Illustration of the most important components (Fig. 01)

4.4. Type plate and Serial number

The type plate (see fig. 02) and the serial number are located on the cross brace.

- A** Manufacturer logo
- B** Manufacturer information
- C** Medical product
- D** CE-mark
- E** Fastening in cars
- F** Attention! - Observe instructions for use
- G** Important - Observe operating instructions
- H** Model designation
- I** Model number
- J** Max. load
- K** Production date
- L** UDI number
- M** Serial number

5. Assembly / Adaptations

5.1. Setting up the wheelchair

To unfold your wheelchair after unpacking, stand in front of the wheelchair.

- Grasp the armrest or seat tube on one side of the wheelchair with one hand and tilt it towards you (Fig. 3).
- With the other hand, push the cross brace apart until the seat tube is fully seated in the plastic holder (Fig. 4).
- The seat cushion must be fully spread out.
- Operate the parking brakes on the right and left of the drive wheels.
- Hang up the leg supports. (see chap. 5.4.)
- Now you can sit down.

5.2. Folding of the wheelchair

To fold your wheelchair, first swing up the foot plates (Fig. 5).

- Stand next to the wheelchair
- Grasp the seat cushion at the front and rear and pull it up (Fig. 6).
 To stow the wheelchair in a passenger car, you can easily remove the rear wheels via the quick-release axle.

PINCHING! Make sure that you do not embrace the seat tube with your hands!

5.3. Foot plates

- The lower leg length can be adjusted in five steps (Fig. 7).
- To do this, loosen both hexagon socket screws, inside and outside.
- Set the desired height.
- Tighten both screws again.
- Proceed in the same way on the other side.

5.4. Legrests

- To make it as easy as possible for you to get in and out, you can fold the foot plates upwards and swing the leg supports outwards or remove them completely.
- To swing the legrest outwards, open the lock (Fig. 8).
- In the „outside position“ the legrest can be completely removed (Fig. 9, Fig. 10).
- To assemble the legrest, proceed in reverse order.
- Swing the legrest forward again until you hear it click into place.

5.4.1. Horizontal adjustable legrests (optional)

Assembly, disassembly and swivelling down are carried out in the same way as for the standard leg supports.

To adjust, loosen the clamping lever, then adjust the legrest to the desired length and tighten the clamping lever.

The horizontal adjustment works as follows:

- Push the locking lever forward.
- By raising or lowering the footplate support tube, bring the legrest into the desired position.
- Release the locking lever; the legrest locks automatically.

Never climb onto the foot plate and watch your fingers when swivelling, there is a danger of getting caught.

The legrest must be swivelled away or removed in order to sit down. Never sit on the legrest. Danger of tipping over!

5.5. Amputation leg support (optional)

Assembly, disassembly and swivelling down are carried out as with the standard legrests, adjustment as with the horizontally adjustable legrest.

The legrest must be swivelled away or removed in order to sit down. Never sit on the legrest. Danger of tipping over!

5.6. Tipping aid

- The tipping aid makes it easier for the attendant to tilt your wheelchair, e.g. to drive over a step.
- To do this, press the tipping aid down with the foot (Fig. 11).
- You keep your balance by using the push handles.
- The tipping aid protrudes to the right and left of the frame and is fitted with rubber profiles to ensure sufficient grip.
- Optionally, anti-tip wheels with integrated tipping aid can be fitted instead of the rubber profiles.

5.7. Side parts

- To make it easier to get in/out of the chair, the two side sections with the armrests can be swivelled upwards.
- Operate the locking lever with light finger pressure and at the same time swivel the side part upwards (Fig.12 / 13).
- To lock the side part, fold it forwards and let it engage with a slight pressure on the armrest.
- The side panels can be removed completely.
- Fold up the side parts as described. Then pull out the knob at the back of the fastening (Fig. 14) and at the same time pull the side part on the armrest upwards away.

When swivelling the side parts, make sure that no items of clothing or other utensils can be trapped. Check that the side panels are firmly in place before putting them into operation!

5.8. Armrests

If necessary, the armrests can be angled in the area of the wrists so that your hands rest comfortably or the armrest is shorter.

- Press the push button on the outside, the front part of the support moves inwards towards the seat surface (Fig 15).
 - Now you can fold the front part down until it snaps back into place and can no longer be moved (Fig 16). You can also push the entire front part inwards and fold it down.
- You can adjust the height of the armrests:
- Press the button and at the same time pull the armrest up to the desired position (Fig 17).
 - Release the button at the desired height.
 - Check that the armrest is locked in place so that it cannot be moved.

5.9. Seat

5.9.1. Seat depth

You can extend the seat depth (see chapter 7. Technical data).

- Swivel the side parts upwards (see chapter 5.7).
- Loosen the five screws of the seat cover on each side with a Phillips screwdriver.
- Pull out the front tube on the seat by the distance of the threaded bushes, approx. two cm (Fig 18).
- Also move the seat cover forward by the distance of the threaded bushes (approx. 2 cm).
- Secure the seat cover by tightening the screws in new front threaded bushes.

5.9.2. Seat height

The seat height can be changed by adjusting the drive wheel and the front wheel at the same time (see chapter 5.11). Please note that the drive and front wheels must be adjusted in the same direction and by the same amount.

5.10. Seat angle

If you tilt the seat of the wheelchair backwards, you will sit lower and more firmly in the wheelchair. To tilt the seat backwards, you have to attach the quick-release axle (see chap. 5.11.) in a position further up. If you want to tilt the seat forwards, attach the quick-release axle in a position further down. With the front wheels it is the other way round. If you attach the front wheel further down in the fork, the seat will tilt backwards. If you fix the front wheel further up in the fork, the result is a negative seat tilt towards the front.

5.11. Wheels

5.11.1. Drive wheel height

To facilitate the work, remove the drive wheel by pressing on the quick-release axle (Fig. 19 /Chapter 5.20). To adjust the height of the drive wheel and thus raise or lower the entire wheelchair, you must first remove the drive wheel bracket.

For Pyrolight


Using a spanner on the inside (Fig. 20) and an Allen key (Fig. 21) on the outside, you can loosen the fastening nuts (Fig. 20) of the axle bracket. Now you can remove the bracket and attach it at the desired height. Tighten the nut again and attach the drive wheel. Make sure that both drive wheels are mounted at the same height.

For Pyrolight XL

Use two spanners to loosen the nut (Fig. 22). Now you can remove the bracket and attach it at the desired height. Tighten the nut again and attach the drive wheel. Make sure that both drive wheels are mounted at the same height.


If the drive wheel height is changed, the knee lever brake must be readjusted.

5.12. Height of the front wheel

The height of the front wheel can be adjusted in four steps. Hold the hexagon socket screw (Fig. 23), loosen the nut on the opposite side and pull out the wheel axle (screw). Place the front wheel at the desired height and tighten the screw and nut again. Note that both front wheels are mounted at the same height

5.13. Wheelbase extension (optional)

The wheelbase extension improves stability and thus reduces the risk of tipping over backwards. However, this makes it more difficult to overcome obstacles or to tip over.

5.14. Parking brakes

The Pyro Light is equipped with a knee lever brake that can be operated by the occupant. This is only to be used for parking. It is not suitable for braking the ride! The handrims are used for this purpose and light pressure is applied to them with the palms of the hands. When the wheelchair has been braked in this way, apply both knee lever brakes by pushing the brake levers forward. The wheelchair is now safely braked (Fig. 24).

- To ensure a safe braking effect at all times, the knee lever brake can be adjusted by changing the distance (Fig. 25).
- To adjust the distance, loosen the two hexagon socket screws on the brake bracket (Fig. 26).
- Then slide the brake system along the bracket into the correct position.
- The brake bolt must be aligned horizontally, i.e. parallel to the ground.
- Finally, tighten both hexagon socket screws again.

- Make sure that the brake fully engages when the lever is operated and does not snap back.

The adjustment on the brakes may only be carried out by trained specialists. Bear in mind that incorrectly adjusted brakes can have life threatening consequences.


Never use the brake lever as a support - e.g. when transferring or standing up - otherwise there is a risk of tipping and falling! The lever could break off!

The parking brake must not be used as a operating brake, as there is a risk to life due to the fall risk

5.15. One-hand brake (optional)

With the one-hand brake you have a parking brake where you brake both drive wheels by operating only one brake lever.


5.16. Operation of the parking brake with brake lever extension (optional)

The operation of the parking brake with brake lever extension is the same as the operation of the parking brake.

Fold in brake lever extension:

- The brake lever can be folded in while driving or when transferring sideways.
- Pull the plastic handle out of the lever (Fig. 27).
- Fold down the plastic handle (Fig. 28).
- Before driving off, the brake must be returned to the operating position in reverse order.

Never use the brake lever as a support - e.g. when transferring or standing up - otherwise there is a risk of tipping and falling! The lever could break off!

The parking brake - with and without brake lever extension - only acts on one drive wheel and must therefore always be used on both sides!

5.17. Drum brake (optional)

The wheelchair can be optionally equipped with a drum brake that can only be operated by an attendant via the brake levers on the push handles (Fig. 29/1).

The drum brake can be locked for parking by means of the locks on both brake levers (Fig. 29/2).



Adjusting the drum brake

The basic adjustment of the drum brake is made at the cable clamping screw. The fine adjustment is made via the adjusting nut on the brake cable. If the adjusting nut is turned up, the braking effect increases (Fig. 30).

- The drum brakes only act on one drive wheel at a time!
- Drum brakes are also subject to natural wear. If the braking effect decreases, the wheelchair must no longer be operated. To avoid this, for your own safety, you should check the brake system for proper functioning before every journey!
- After completing the adjustment process, tighten all screws and carry out a brake test!

The adjustment on the brakes may only be carried out by trained specialists. Bear in mind that incorrectly adjusted brakes can have lifethreatening consequences.

5.18. Backrest

To enable you to sit comfortably, the hardness of the backrest upholstery can be adjusted.

5.18.1. Backrest cover

The tighter you adjust the backrest upholstery, the more harder the backrest feels. First remove the backrest padding attached with a Velcro Fastener.

You can then tighten or loosen the individual straps (Fig. 31).

To relieve certain regions of your back, tighten the straps there a little more loosely. Now you have to fasten the backrest upholstery again with the Velcro.

5.19. Push handles

To allow an accompanying person to push comfortably, the push handles are adjustable in height in three steps. Swivel the side parts upwards. Loosen the hexagonal screw on the inside of the push handle (Fig. 32). Remove the screw completely. Pull the tube with the push handle to the desired height (Fig 33). Then reinsert the screw and tighten again with the lock nut.

5.20. Thru axles

The spoked wheels are attached to the chassis with a quick-release axle (Fig. 34).

- Press the locking button of the quick-release axle with your finger and pull off the wheel (Fig. 34).
- When reattaching the wheel, simply put it back on the holder. The locking button must also be pressed in.

5.21. One-hand operation (optional)

With the one-hand operation you have two handrims on one wheel, so you can drive both drive wheels with one hand.

The option can be mounted either on the left or on the right.

5.22. Passive lighting

Your wheelchair is equipped with reflectors on the back and in the spokes as standard. These must not be covered.

5.23. Accessories

5.23.1. Anti-tip wheels (optional)

Always use the anti-tip wheels in pairs.

The anti-tip wheels prevent the wheelchair from tipping backwards. The anti-tip castors must have a distance of approx. 3-5 cm from the floor (Fig. 35).

The anti-tip wheels are particularly recommended for inexperienced wheelchair users or for wheelchair users with leg amputations.

- The anti-tip wheels are attached as a pair to the left and right of the frame.
- First pull the rubber profiles of the tipping aid off the frame.
- Slide the round profile of the anti-tip wheels onto the side frame tube.
- Now screw the anti-tip wheels into the holes provided on the frame.
- Adjust the height by pulling the washer of the holder towards you and then snap the tube with the anti-tip wheels into place at the appropriate height (Fig. 36).
- The anti-tip wheels can also be used as a tipping aid when the anti-tip wheels are raised.

5.23.2. Lap belt (optional)

A lap belt can be ordered as an option and secures users who do not have the necessary support in the wheelchair (Fig. 37).

- The Pyro Light can be fitted with a lap belt on the back tube.
- The belt is opened and closed at the locking tab.
- You can adjust the strap length using the sliding strap tab on the strap.

5.23.3. Therapy table (optional)

For side sections with long armrests, the table is simply pushed onto the armrests from the front using the guides provided. To remove the therapy table, simply pull it out again. For „desk“ side sections with arm pads, the table is placed with the enclosed guides and fastened with star grip screws. To remove the therapy table, first loosen the star screws and then pull the table out again.

5.23.4. Togo (optional)

The Togo push aid can be fitted to your wheelchair.

All relevant information on assembly and operation can be found in the operating instructions for the Togo push aid.

5.23.5. Cane holder (optional)

Place the cane in the cane holder and then fix to the backrest.
Never climb on the pole holder and never use it as a tipping aid.

5.23.6. Spoke protector (optional)

The spoke protector is attached to the spokes with cable ties.

5.23.7. Handrim cover (optional)

Put the handrim cover over the handrim to get more grip on the handrim,
for example if you have limited hand strength.

5.23.8. Infusion / Oxygen bottle holder (optional)

The infusion / oxygen bottle holder is attached to the lower frame
tube of your wheelchair.

The risk of the wheelchair tipping over backwards increases with the
increased load of oxygen cylinders. The use of anti-tip castors is
strongly recommended.

5.23.9. Carry-on bag (optional)

If you want to attach a bag to your wheelchair, hang the loops of your
bag over the back tubes of your wheelchair.

Please note that the reflectors must not be covered.

Caution: Increasing the load of the bag also
increases the risk of your wheelchair tipping over
backwards. The use of anti-tip wheels is strongly
recommended.

5.23.10. Headrest (optional)

The headrest is screwed onto the push handles.
After loosening the clamping lever or star screw, you can adjust the
position of the headrest. Then tighten the screws firmly.

6. Use of the wheelchair

6.1. Check before setting off

6.1.1. Checking the brakes

Before each journey, check that the brake system is functioning properly.

When the parking brakes are applied, the wheelchair cannot be pushed away on a dry surface. When the brakes are released, the wheelchair moves without grinding noises and with good directional stability.

If this behaviour does not apply, inform your specialist dealer immediately and discontinue use until repairs have been carried out.

6.1.2. Checking the tyres

Always check that the tyres are undamaged before setting off. If the tyre is damaged, have it repaired by an authorised specialist workshop.

6.2. Side entry and exit

- Bring the wheelchair and the current or future seat as close together as possible to the side.
- If you want to transfer from another wheelchair or indoor undercarriage, secure it by applying the parking brakes.
- Secure your wheelchair against unintentional rolling away by operating the brake levers on both sides of the wheelchair.
- Swing the side parts away to the rear.
- Fold the foot plates to the side or remove the leg supports.

- Now slide sideways onto the other seat.
- Make sure you sit as far back in the seat as possible.
- Finally, hang the leg supports back in the starting position.

Never step on the footplates! Danger of tipping

6.3. Entering and exiting from front

- First swivel the foot plates upwards or remove the leg supports.
- Bring the Pyro Light and the current or future seating as close to each other as possible from the front.
- If you want to transfer from another wheelchair or indoor undercarriage, secure it by applying the parking brakes.
- Secure your wheelchair against unintentional rolling away by operating the brake levers on both sides of the wheelchair.
- Now slide onto the seat by turning your body.
- Finally, hang the leg supports back in the starting position.

6.4. Overcoming stairs and high steps

For high stairs and stairs with more than 3 steps, 2 accompanying persons must help.

- The helpers may only grasp at firmly mounted frame parts. Do not lift the wheelchair by the legrests or armrests.
- When going up stairs, the rear assistant pulls the wheelchair over the step using the permanently mounted handles (Fig. 38).
- The lower assistant grips the front frame tubes (Fig. 39), NOT the legrests, and stabilises the position. In doing so, he presses the rear wheels of the wheelchair against the steps.
- When descending, the lower helper brakes by pushing the wheelchair straight against the steps. The upper helper holds the wheelchair by the push handles, secures it and keeps it in the correct position. The wheelchair should roll down one step at a time.

The lifting points for the transport of the wheel chair are exclusively the frame and the push handles. It is not allowed to use the armrest for lifting the wheelchair!

6.5. Driving with the wheelchair

Sit in your Bischoff & Bischoff wheelchair. Make sure that you sit as far back as possible. Locomotion and braking is done via the handrims on the drive wheels (Fig. 40). Carry out the first driving attempts carefully until you have become accustomed to the wheelchair and its driving behaviour. Wheelchairs have only limited tipping and sliding safety. Take special care when braking, starting or turning on slopes.

7. Technical data

| Dimensions | |
|--|---|
| Seat width: | 380, 410, 430, 450, 480, 520 mm / 510/560 mm* |
| Seat depth: | 420, 440 mm / 440, 460 mm* |
| Seat height front (without cushion): | 460, 475, 490 mm / 500 mm* |
| Rear seat height back (without cushion): | 440, 460 / 480* |
| Seat angle: | 0 - 4,5° |
| Total width: | Sitzbreite +200 mm |
| Total width folded: | 290 mm |
| Total length: | |
| with leg supports: | 1040 mm / 1070 mm* |
| without leg supports: | 810 mm / 840 mm* |
| Total height: | 945 mm |
| Backrest height (without cushion): | 400, 420, 440 mm |
| Lower leg length: | 390 - 470 mm, 410 - 490 mm* |
| Backrest angle to the seat surface: | 100° |
| Distance between armrest and seat: | 190, 200, 210, 220, 230, 240, 250, 260 mm |
| Length of armrest: | 350 mm |
| Width of armrest: | 50 mm |
| Push handle height: | 945 mm |
| Dimension folded: | 805/830* x 290 x 940 mm (L/B/H) |
| Weights | |
| Max. load capacity: | 125 kg, 170 kg* |

| | |
|---|--------------------------------------|
| Unload weight: | 18,6 kg |
| Heaviest part: | 10,2 kg |
| Driving characteristics | |
| steering range: | 1120 mm |
| Stability | |
| Static uphill / downhill: | 10° |
| Static sideways: | 10° |
| Parking brake up to inclination: | 7° |
| Environmental conditions | |
| Temperature (°C): | -10 - 50 |
| Humidity (%): | 20 - 80 |
| Storage conditions | |
| Temperature (°C): | -20 - 65 |
| Humidity (%): | 20 - 60 |
| Wheels | |
| Drive wheels: | 24" x 1 3/8" PU (pannensicher) |
| Steering wheels: | 6" x 1 3/4" Softgummi (pannensicher) |
| Materials | |
| Frame: | Powder coated steel |
| Armrest cushions: | PU foamed |
| Seat / back cushion: | Nylon and foam |
| Flame retardand material: | yes (EN 1021-2) |
| *Model Pyro Light XL | |
| Accessories: anti-tip wheels, pelvic belt, therapy table, Togo, cane holder, spoke protector, handrim cover, infusion/oxygen bottle holder, headrest, one-hand operation drum brake | |
| The dimensions given are subject to manufacturing tolerances of +/- 10 mm | |

8. Care instructions for the user

8.1. Cleaning

The seat and back materials can be cleaned with warm water using a sponge and/or a soft brush. In case of stubborn dirt, a commercial cleaning agent can be added to the water.

Note: Never use aggressive cleaning agents such as solvents / scouring agents or hard brushes.

Plastic parts are best cleaned with a commercially available plastic cleaner. Please observe the product information of the cleaning agent.

- Plastic claddings are attacked by non-ionic surfactants and solvents, especially alcohols!
- Wipe the frame parts of the wheelchair with a damp cloth. For heavier soiling, also use a mild cleaning agent.
- Clean the wheels with a damp brush with plastic bristles (do not use a wire brush!).
- Do not use a high-pressure cleaner.
- Wash the push handles, arm pads, seat and back pads with a mild detergent.
- Upholstery parts can be machine washed at 40°C. Do not machine dry, only dry at room temperature.

8.2. Disinfection

- Before disinfection, the rollator must be thoroughly cleaned. The wheelchair can be disinfected with a household disinfectant. Use approved spray disinfectants or alcohol-based wipe-on disinfectants.
- In the case of multiple product users, a commercial disinfectant is prescribed.
- Only use the wipes approved by the Robert Koch Institute (RKI) and the Association for Applied Hygiene e.V. disinfectants (e.g. Bacillol AF or Kohrsolin FF); see manufacturer's recommendations for exposure

time and concentration.

- Caution: Damaged upholstery parts cannot be sufficiently disinfected. These should be replaced.
- Under certain circumstances, the surfaces can be affected by disinfectants. Longer-term restrictions of the functionality can be the result. Please note our enclosed overview.

Attention: Danger of infection! If there are several users, the wheelchair must be cleaned after each use.

9. Maintenance instructions

We recommend that the wheelchair is inspected by a specialist dealer at least once a year.

In the event of faults or defects on the wheelchair, it must be handed over immediately to the specialist dealer for repair. The repair and replacement of parts is carried out by the specialist dealer.

The necessary information and documents for repair and maintenance will be made available to the specialist dealer on request.

If you have any questions, please contact the Bischoff & Bischoff customer service.

You must check your wheelchair regularly according to the maintenance schedule and find it to be in perfect working order.

9.1. Maintenance advice for the specialist dealer

Maintenance must be carried out in accordance with the Bischoff & Bischoff maintenance plan. A maintenance schedule will be provided upon request. A training programme for the specialised trade is also available.

The wheelchair has to be checked according to the maintenance schedule and be in perfect working order.

10. Repair

Repairs and maintenance may only be carried out by authorised partners (medical supply stores or specialist dealers). Please contact our customer service by telephone to find a partner in your area. You will find the contact details on the back of this instruction manual.

11. Transport

11.1. Transport of the wheelchair

To make the wheelchair as handy as possible for transport, remove all detachable parts (leg supports, drive wheels, armrests, insertable accessories) and fold the wheelchair (see chapter 5.2).

Your wheelchair should only be transported when folded (Fig. 42). When folded, the wheelchair should be secured to fixed frame parts in the vehicle.


11.2. Transport in motor vehicles

You can tell whether or not your wheelchair is approved as a seat in a motor vehicle for transporting mobility-impaired persons by the symbols on the type plate on the product (see chapter 4.4):

Wheelchair not suitable as a seat in a motor vehicle for transporting mobility-impaired persons (Fig. 43)

Wheelchair suitable as a seat in a motor vehicle for transporting persons with reduced mobility (Fig. 44)

All products approved by Bischoff & Bischoff as seats in vehicles are tested according to ISO 7176-19.

| | |
|----------------------|-------------------------------------|
| Pyro Light | May be used as a vehicle seat. |
| Pyro Light XL | Must not be used as a vehicle seat. |

If you have any queries regarding the securing of your product in a motor vehicle for the transport of persons with reduced mobility, please contact the Bischoff & Bischoff customer service. You will find the contact details on the back of these operating instructions.

If your wheelchair is not crash-tested, it must never be used as a seat in a motor vehicle for transporting mobility-impaired persons. There is a danger to life!

12. Passing on the wheelchair

The wheelchair is suitable for reuse. Before it is passed on to another user, it must be serviced and hygienically prepared by a specialist dealer. When passing on the wheelchair, please remember to give this user manual and the address of your specialist dealer to the new user.

13. Storage / Shipping

If the wheelchair is to be stored or shipped, all insertable and unattached parts, except the drive wheels, must be removed and packed individually in suitable cardboard boxes. The individually packed parts can then be packed together in a larger box. It is recommended to keep the original packaging and store it for these purposes so that it is available when needed. Then your wheelchair is optimally protected against environmental influences during storage or transport.

For the storage of your wheelchair, please observe the information in the technical data.

To avoid mould and damage to the upholstery parts, the storage location should be as dry as possible and not exposed to direct sunlight.

14. Disposal

If your wheelchair is no longer used and is to be disposed of, contact your specialist dealer.

If you want to dispose of the waste yourself, ask local recycling companies about the regulations in your area.

15. Warranty

1. For products supplied by Bischoff & Bischoff, the period of warranty is 24 months from purchase date. The products are delivered free from production and material defects. In case of evidence of material or manufacturing defect, defective parts will be replaced free of charge.
2. This warranty excludes claims, if a repair or replacement of a product or a part is required for the following reasons:
 - a) normal wear and tear, which include but is not limited to the following parts where fitted; Batteries, motor brushes, handles, armrests, upholstery, tyres, brakes, caps, etc.
 - b) overloading the product such as, for example, by exceeding the maximum user weight or load.
 - c) the product or part has not been maintained or serviced in accordance with the Instruction manual, the care and hygiene instructions or the manufacturer's recommendations as shown in the maintenance instructions.
 - d) Accessories have been used which are not original accessories.
 - e) The product or part has been damaged through negligence, accident or incorrect use.
 - f) changes/modifications have been made to the product or parts, which deviate from the manufacturer's specifications.
 - g) Repairs have been carried out, before our customer service has been informed of the circumstances.

3. To enforce the warranty, please contact your dealer with the exact details of the nature of the difficulty. The product must be repaired by a dealer authorised by Bischoff & Bischoff.
 4. For parts, which have been repaired or exchanged within the scope of the warranty, the warranty will be extended to the remaining warranty period for the product in accordance with point 1).
 5. For original spare parts which have been fitted at the customer's expense, these will have a 12 months guarantee, following the fitting.
 6. If the repair fails after a reasonable time, the customer can, at his discretion, request a reduction, compensation or the termination of the purchase contract.
 7. The guarantee is subject to the law of the country in which the product was purchased from Bischoff & Bischoff & Bischoff.
- Please also observe the Bischoff & Bischoff care, warranty, hygiene and maintenance instructions. Your dealer will be happy to provide these to you.
- For replacement devices and repairs, the warranty period is one year.

Non-observance of the operating instructions as well as improper use, maintenance jobs and technical changes or adaptations without the permission of the Bischoff & Bischoff company, lead to guarantee loss and the loss of product liability in general.

Warranty certificate

Produkt: Pyro Light / XL

Model number:

Serial number:*

* (to be entered by the dealer)

Specialist retailer:

Date and stamp:

(Pay attention to our general terms and conditions)

*Notice signs are located on the cross brace

Note!

Print errors, product and colour changes and technical information are excepted. The colours of the products shown may be slightly different in real life. The images may show accessories. The current version of each Instruction manual is ready for download on our website:

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