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---

**Content**

<b>1. Introduction</b>	<b>64</b>		
1.1. Signs and symbols	65		
<b>2. Important safety instructions</b>	<b>65</b>		
2.1. General safety instructions	65		
2.2. Safety in vehicles	67		
2.3. Participation in road traffic	67		
2.4. Notes on EMC interference	68		
2.5. Brakes	68		
2.6. Indications	69		
2.7. Contraindications	70		
2.8. Declaration of conformity	70		
2.9. Responsibility	70		
2.10. Information about range	70		
2.11. Service life	71		
<b>3. Purpose</b>	<b>71</b>		
<b>4. Description of product and contents</b>	<b>72</b>		
4.1. Inspection on delivery	72		
4.2. Contents	72		
4.3. Type plate and serial number	73		
4.4. Summary	74		
<b>5. Settings</b>	<b>75</b>		
5.1. Basic settings	75		
5.1.1. Mechanical footrests	75		
5.1.2. Electronic footrests	76		
5.1.3. Seat	77		
5.1.4. Electronic adjustment functions	79		
5.1.5. Sag of the back rest	79		
5.1.6. Seat suspension	80		
5.1.7. Armrests	80		
5.1.8. Side parts	81		
5.1.9. Control unit	81		
5.1.10. Lap belt	82		
5.2. Accessories	82		
<b>6. Using the wheelchair</b>	<b>84</b>		
6.1. General functions	84		
6.2. Description of functions of VR-2 control	86		
6.3. Functional description of the VR2 control unit for electronic adjustment options (optional)	88		
6.4. Functional description of the VR2 attendant control unit	90		
6.5. Locking and unlocking of the wheelchair	93		
6.6. Description of R-Net control	94		
6.7. Functional description of the R-NET control unit for electrical adjustment options (optional)	95		
6.8. Functional description of the R-NET attendant control unit	96		
6.9. Locking and unlocking of the wheelchair	98		
6.10. Adjustments to the Control unit	99		
6.11. Getting in and out of the chair	99		
6.12. Leg supports	100		
6.13. Check before use	100		
6.13.1. Brake check	100		
6.13.2. Checking lights, indicators and horn	101		
6.13.3. Checking the tyres	101		
6.14. Batteries	101		
6.14.1. Charging	102		
6.14.2. Replacing batteries	103		
6.15. Driving the wheelchair	104		
6.16. Pushing the wheelchair	106		
6.17. Head rest (optional)	107		
6.18. Kerb climber (optional)	107		
6.19. Mechanical parking brake	108		
<b>7. Technical data</b>	<b>110</b>		
<b>8. Maintenance tips for users</b>	<b>112</b>		
8.1. Cleaning and disinfection	112		
8.2. Disinfection	112		
<b>9. Maintenance note</b>	<b>113</b>		
9.2.1. Maintenance advice for the specialist trade	114		
<b>10. Repair</b>	<b>114</b>		
<b>11. Parking</b>	<b>114</b>		
<b>12. Lifting and holding points</b>	<b>114</b>		
<b>13. Transport</b>	<b>115</b>		
14. Passing on the wheelchair	115		
<b>15. Storage and dispatch</b>	<b>115</b>		
<b>16. Disposal</b>	<b>116</b>		
<b>17. Warranty conditions</b>	<b>116</b>		

## 1. Introduction

Dear user,

You have decided on a high-quality Bischoff & Bischoff electric wheelchair, and we would like to thank you for your trust.

The electric wheelchair was designed to offer you many advantages and meet your needs. The wheelchair can be fitted to your individual body size thanks to a multitude of adjustable features and optimal adaptability such as the length of the lower leg, the recline of the backrest, the sag of the backrest as well as the swing away function of the control unit and the choice of mounting the control unit either to the left or the right side.

The Taiga can be supplied with two different control systems. In the standard configuration, the Taiga is delivered with the VR-2 controller, which offers you the option of adjusting up to two components electronically. The R-Net controller permits electrical adjustment of the leg rests and the seat and back angles, or when more than two electrical options are needed. Thus, the Taiga can easily be customised to your specific needs.

The electric wheelchair is designed for indoor and outdoor use.

Before using your new electric wheelchair for the first time, please read and follow the instructions in this manual carefully. The instruction manual is an integral and necessary part of the wheelchair. This is why you should keep it close by and include it when passing on the wheelchair to someone else. For visually impaired users, this document is available as a PDF file at our website [www.bischoff-bischoff.com](http://www.bischoff-bischoff.com).

Any repairs and adjustments require special technical training and may therefore only be carried out by authorized Bischoff & Bischoff dealers.

## 1.1. Signs and symbols

Important! This is especially useful information on a given subject.

Warning! These are especially safety-relevant instructions. Follow the directions in the manual!

## 2. Important safety instructions

### 2.1. General safety instructions

General safety instructions

- To avoid falls and dangerous situations, you should first practise using your new wheelchair on level, manageable terrain. An accompanying person is recommended in this case.
- Before using the wheelchair, please check that all attached parts are properly fastened.
- Before using the electric wheelchair for the first time, practise riding on level and manageable terrain. Familiarise yourself intensively with the braking and acceleration behaviour when driving straight ahead and cornering.
- Please note that the risk of tipping to the side and back may increase when shifting the balance due to body movements or loading the wheelchair.
- The wheelchair must not be misused to transport several people or loads. Observe the intended use.
- Never use the wheelchair while under the influence of alcohol or other substances that affect alertness or physical and mental receptiveness.
- Please be sure to use access ramps to overcome obstacles.

- Avoid driving into an obstacle (e.g. step, kerb) without braking.
- Observe the anti-tip protection on stepped terrain, lifting platforms and ramps. Ensure sufficient clearance above, below and to the sides.
- The control unit must always be switched off when you get in or out of the vehicle.
- The road traffic regulations must be observed in road traffic.
- The rear lights on the wheelchair must always be visible to other road users.
- If possible, wear light-coloured, conspicuous clothing. This will make you more visible to other road users.
- Please note that upholstered parts that are exposed to direct sunlight can heat up and cause skin injuries if touched. Therefore, cover these parts. Be aware of the risk of trapping between moving parts when adjusting, using and maintaining the wheelchair.
- The seat and upholstery materials fulfil the requirements for resistance to flammability in accordance with EN 1021-2/ ISO 8191-2 Nevertheless, do not place any upholstery or other parts on fire sources- such as cigarettes.
- Please note the increased risk of entrapment for the user, the occupant and third parties when using the electric seat adjustment.
- Always adjust the seat and backrest to an upright position when travelling with the wheelchair.
- Travelling on slopes is only permitted if the seat and backrest are in an upright position.
- At low ambient temperatures, metal parts of the wheelchair can cool down considerably, so avoid skin contact.
- If you notice any unusual behaviour or damage to your power wheelchair, stop using it and contact your specialist dealer.

- Reporting of incidents.

If you as the operator, user or their relatives discover possible defects or functional limitations, please contact an authorised immediately. Operators, users or their relatives should inform their authorised dealer who provided them with the product about suspected serious incidents that have had or could have directly or indirectly one of the following consequences and which affect them. You can also report this to the competent higher federal authority:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration in the state of health of a patient, user or other person,
- a serious threat to public health.

## 2.2. Safety in vehicles

The Taiga is not suitable for use as a vehicle seat. It must not be used as such under any circumstances. It does not fulfil the requirements of ISO 7176-19 and failure to do so may result in serious injury or death in the event of an accident (Fig. A).

## 2.3. Participation in road traffic

The wheelchair is designed for indoor and outdoor use. Please note that you must take part in public road traffic and comply with the road traffic regulations.

Do not endanger other participants by driving recklessly, especially on pavements

## 2.4. Notes on EMC interference

Despite compliance with all applicable EMC guidelines and standards, it is possible that the power wheelchair may be disturbed by or interfere with other electrical devices (e.g. mobile phones). If you notice such behaviour, either switch off the external device or move your power wheelchair out of the range of interference.

## 2.5. Brakes

Service brake:

In driving mode, the motors act as a service brake. To brake the wheelchair slowly, move the steering and drive lever (joystick) slowly back to the neutral position (starting position).

To brake harder, release the steering and drive lever.

The braking distance is significantly longer on a downhill gradient than on a flat road. Bear this in mind when descending ramps and brake in good time.

Make sure you brake your power wheelchair in good time. Especially in front of people and when travelling on slopes.

Parking brake/parking:

Switch off the wheelchair using the on/off switch (see chapter 6.1). If the drives are engaged, the power wheelchair cannot be pushed or moved (see chapter 6.16). The parking brake function complies with the provisions of the German Road Traffic Licensing Regulations (StVZO).

Sliding mode:

To be able to push the wheelchair, the unlocking levers must be actuated (see chapter 6.6). If the wheelchair is switched on, the Battery status display. You can switch off the wheelchair for pushing mode.

Emergency brake:

For an emergency stop while travelling, press the on/off switch. This results in very heavy braking. This can lead to injuries and should therefore only be used in the event of danger

An emergency stop may only be carried out in the event of danger, as this can lead to injuries.

## 2.6. Indications

Inability or severely impaired ability to walk due to

- paralysis
- loss of limb
- defect/deformity of limb
- joint contracture
- joint damage (not in both arms)
- other diseases

The provision of an electric wheelchair is indicated if the use of a hand-driven wheelchair is no longer possible but the proper operation of an electric drive is possible.

## 2.7. Contraindications

The provision of an electric wheelchair is unsuitable for persons

- with severe balance disorders
- with impaired or insufficient vision
- with strongly impaired cognitive abilities

## 2.8. Declaration of conformity

We Bischoff & Bischoff GmbH declare on our sole responsibility, that the wheelchair Taiga complies with all requirements of MDR2017/745 which are applicable.

## 2.9. Responsibility

We can only provide a warranty if the product is used under the specified conditions and for the intended purposes, if modifications, extensions, repairs and maintenance work are only carried out by persons authorised by us to do so and if the wheelchair is used in accordance with all instructions for use.

## 2.10. Information about range

The range of the electric wheelchair Taiga is 35 km. The range that can be achieved by electric wheelchairs depends on the following factors:

- Composition of the carriageway
- Way of driving
- Battery condition
- Vehicle load weight
- Ambient temperature

Driving on inclines has a negative effect on the range.

Please note that the range specified by us was determined under test conditions:

- Batteries were new and fully charged
- Ambient temperature of 21 °C
- User weight of 85 kg
- Level, firm surface
- Constant speed
- Deactivated lighting system

## 2.11. Service life

If used as intended, and the safety, care and maintenance instructions are observed, the expected service life is up to 5 years. Beyond this time the electric wheelchair can be used continuously, provided it is in a safe condition.

## 3. Purpose

Unless otherwise indicated, the instructions in this operating manual are intended for the occupant as the operator.

The Taiga electric wheelchair is able to handle long distances and surmount obstacles outdoors. It is classified under category B according to the EN 12184 standard. The maximum load capacity is 170 kg.

## 4. Description of product and contents

### 4.1. Inspection on delivery

All Bischoff & Bischoff products undergo a proper final inspection in our company and are labeled with the CE mark.

The wheelchair is delivered in a special carton. Keep the carton after unpacking, if possible. It is useful for later storage of the product or for sending it back. The Taiga wheelchair comes disassembled in its component parts. Assembly, adjustment to your body size and first instructions are taken care of by your qualified medical supply dealer.

If the wheelchair is shipped by rail or road, the goods should be checked immediately for transport damages (carton) in the presence of the bearer. Should there be any damage, please contact our customer service immediately.

Check the contents for completeness and make sure nothing is damaged. Should you find any irregularities or damage, please contact our customer service.


### 4.2. Contents

After receipt of the goods, please check the contents for completeness right away. The contents consist of:

- outer package
- electric wheelchair (ready for use incl. 2 batteries)
- charging unit
- instruction manual

### 4.3. Type plate and serial number

The type plate (fig. 01) and serial number are located under the seat.

- (A) Manufacturer logo
- (B) Manufacturer information
- (C) Caution! Follow instructions for use
- (D) CE mark
- (E) Not suitable for use as a seat in a vehicle
- (F) Caution! - Follow instructions for use
- (G) Important - Follow instructions for use
- (H) Model designation
- (I) Model number
- (J) Application class
- (K) Date of manufacture
- (L) Maximum weight capacity
- (M) Maximum speed
- (N) UDI number
- (O) Serial number

#### 4.4. Summary

#### 5. Settings

##### 5.1. Basic settings

It is necessary to adjust the wheelchair to your body size before first use. This is the only way to ensure a physiologically correct posture allowing for long drives without tiring. Below you find a description of the basic settings of the wheelchair. You should carry out these adjustments with the help of your medical supply dealer. You need a set of 5 mm to 8 mm Allen keys.

Make sure to always use two screws for secure fastening!

##### 5.1.1. Mechanical footrests

###### **Length of lower leg supports**

- The lower leg supports can be adjusted from 470 to 550 mm.
- To do so, loosen the two Allen screws (fig. 03).
- Adjust to the desired height.
- Tighten the screws again. Use medium-strength Loctite for fixation.

###### **Angle of footrest**


- The angles of the footrests are adjustable in 25° steps.
- Loosen the screwing at the side of the footrest.
- Adjust to the desired angle.
- Tighten the screw again. Use medium-strength Loctite for fixation.

### 5.1.2. Electronic footrests


#### Length of lower leg supports

- The lower leg supports can be adjusted from 420 - 580 mm.
- To do so, loosen the two Allen screws.
- Adjust to the desired height (fig. 05).
- Tighten the screws again.
- Use medium-strength Loctite to lock screws in place.

#### Angle of footrest

- The angle of the footplate is fully adjustable.
- Undo the screw to the side of the footplate.
- Adjust the footplate to the desired angle (fig. 06).
- Tighten the screw again.

Make sure there is an appropriate distance between footrests and the ground.

#### Height of the calf pad

- The height of the calf pad is adjustable.
- Undo the screw to the side of the footrest.
- Adjust the calf pad to the desired height (fig. 07).
- Use medium-strength Loctite for fixation.


### 5.1.3. Seat

#### Seat width (mechanical)

- The seat width can be adjusted in three steps from 415 mm to 515 mm.
- Ensure that the wheelchair is switched off.
- Remove the attachment parts such as the leg rests, side parts, back pad and seat pad.
- Detach the Velcro straps of the back rest cover (mechanically).
- You first adjust one half of the seat (left or right) and then the other.
- For each side, you must loosen two fastening screws under the seat (fig. 08) and remove two Allen head screws from the back rest (fig. 09).
- Move the seat half to the desired width. Markings on the seat (fig. 10) serve as orientation aids. Shift the seat and back rest equally.
- Secure the screws to the back rest and tighten them and the fastening screws under the seat. Repeat the procedure with the other half of the seat.
- You can pull the Velcro straps of the back rest taut again and reattach all other attachment parts.

Ensure that both halves of the seat are adjusted equally.

**Seat depth**

- The default setting of the seat depth is 450 mm.
- You can reduce the seat depth to 410 mm . You should then use a newly fitted seat cushion.
- To adjust the seat depth, undo the fastenings of the backrest, slide it forward and tighten the screws again (fig. 11).

**Seat angle - mechanical**

This description only applies, if you have a configuration with a mechanically adjustable seat angle.

- The default setting of the seat angle is 1° backwards.
- The seat angle can be adjusted backward in five settings.
- To do this, undo the fastening screw on the telescopic section at the front under the seat.
- Remove the screws and adjust the seat to the desired angle (fig. 12).
- Place the screws into the nearest possible hole and tighten them again.

**Angle of the back - mechanical**

This description only applies, if you have a configuration with a mechanically adjustable back angle.

- The angle of the backrest of your electrical wheelchair can be easily adjusted. In addition, you also have the option of entirely folding down the backrest forwards.
- To fold down the backrest, remove the fastening pin on the backrest and then fold down the backrest. (fig. 13)
- The back rest angle can be adjusted from 90° to 120° in 10° increments.
- Pull the cable behind the back. While doing so, hold the back rest firmly with your other hand.
- The pins latch into the holes after the cable is released.

Make sure that both pins, on the right and left, are fully engaged.

**5.1.4. Electronic adjustment functions**

Depending on the configuration it may also be possible to adjust the angle of the seat, backrest and footrests electronically. For further details, see chapter 6, "Using the wheelchair".

Beware of the risk of to the user, carer and third parties when electronically adjusting the seat functions.

**5.1.5. Sag of the back rest**

A well-adjusted backrest makes it easier to sit continuously in a relaxed position and is effective at reducing the risk of pressure sores. Ensure that there is sufficient slack in the area of the pelvis.

- First completely remove the pad layer of the back cover.
- Undo the Velcro straps of the tensioning straps (fig. 14).
- The straps are adjusted one by one, starting from the bottom.
- Place the pad layer of the backrest over the tensioning straps.


### 5.1.6. Seat suspension

- You can adjust the seat suspension to your weight for comfortable driving (fig. 15)
- A special spanner is required in order to adjust the spring travel (please contact your specialist trade).
- Compression of the spring is equivalent to a low spring travel and a „harder“ setting.
- Releasing the tension of the spring increases the spring travel and ensures a „softer“ setting.

Ensure that both shock absorbers, right and left, are adjusted the same.

### 5.1.7. Armrests


The armrests can be completely taken off to be able to get in and out of the chair from the side. In addition, the height of the armrests can be adjusted. They can also be adjusted sideways by adjusting the side parts (see chap. side parts).

#### Removing the armrests

- To remove the armrest, loosen the wing bolt and pull out the armrest to the top (fig. 16).
- When installing the armrest, make sure it is inserted all the way and fastened again.

#### Adjusting the height

- First remove the arm rest from the bracket.
- Remove the check screw (fig. 17) and screw it into the desired position of the arm rest.
- The lowest hole corresponds to the max. height setting of the arm rest.

### 5.1.8. Side parts

You can slide the side parts sideways.

- Undo the fastening screw underneath the seat and slide the side part into the desired position (fig. 18).
- Tighten the screw again.

Ensure that the adjustment is restricted to the marked range!

### 5.1.9. Control unit

The control unit can be mounted either on the right or left-hand side underneath the armrest each. The control unit must be adjusted by a qualified specialist trade to warrant safe functionality of the wheelchair.

#### Distance from the armrest

- You can adjust the position of the control unit to your forearm length.
- Loosen the hand screw on the outer side of the side part.
- Move the control unit into the desired position, then tighten the screw again (fig. 19).

#### Adjusting the height

- The distance to the armrest can be adjusted.
- Loosen the screw connection on the inner side of the side part (fig. 20).
- Move the control unit into the desired position, then tighten the screw again.

Please note that the wires to the control unit must be re-routed. Have this done by a technician.

### 5.1.10. Lap belt

The lap belt is attached to the side of the backrest.

- Adjust the length of the lap belt by positioning the buckle accordingly (fig. 21).
- Fasten the seatbelt with the latch.
- Open the lap belt by pushing the red button at the latch.

Caution: The lap belt must always be fastened during the journey.

## 5.2. Accessories

### 5.2.1. Electrically adjustable backrest (optional)

It is possible to equip your wheelchair with an electrically adjustable backrest for comfortable adjustment of the backrest angle.

### 5.2.2. Headrest (optional)

The wheelchair can also be equipped with a height- and depth adjustable headrest.

### 5.2.3. Horizontal and width adjustable leg rest, mechanical or electric (optional)

For comfortable adjustment of the legrest, it is possible to equip your wheelchair with an electrically adjustable legrest that can be adjusted horizontally and in width.

### 5.2.4. VR-2 companion control (optional)

To be able to control your wheelchair from an accompanying person, a companion control is available.

### 5.2.5. R-Net control (optional)

An R-Net control system is available for up to 4 electrical adjustment options.

### 5.2.6. R-Net attendant control (optional)

If you would like to be able to control your wheelchair with the R-Net control unit from an accompanying person, it is possible to attach an attendant control unit.

### 5.2.7. Rear-view mirror (optional)

An optional rear-view mirror is available for better all-round visibility.

### 5.2.8. Kerb lifter (optional)

The kerb lifter makes it easier for you to negotiate steps of up to 100 mm.

### 5.2.9. Cane holder (optional)

There is an optional cane holder so that you can take your cane with you while travelling.

### 5.2.10. Comfort seat (optional)

The power wheelchair can also be equipped with more comfortable seat and back cushions for more comfortable sitting.

## 6. Using the wheelchair

After adjusting the wheelchair to your personal needs, familiarize yourself with the control functions.

### 6.1. General functions

#### Light

Use this button to switch on the active lighting system of your wheelchair when darkness falls or when travelling through dark areas such as tunnels. This will make you more visible to other traffic (fig. 22).

#### Warning indicators


Press this button to switch on the warning indicators in emergency situations, such as inability to move the wheelchair. All four indicator lights will flash and signal other traffic that your wheelchair is disabled (fig. 23).

#### Indicators (left/right)

Use this buttons (fig. 24) to switch on the left or right indicators of your wheelchair to signal other traffic that you wish to change your direction of travel to the left or right. Press the button again to switch off the indicator.


#### Battery status display

The battery status display (fig. 25) shows the charge status. If the display is in the green range, the battery is sufficiently charged. If the display is in the red range, the wheelchair can only be used for a short distance and the battery should be charged as soon as possible (see charging procedure).


#### ON/OFF button

Press this button to switch on the wheelchair (fig. 26).

#### Display of mode

The mode display shows the preselected mode. Settings from 1 to 5 LEDs are available. The modes indirectly define the maximum available speed when the joystick is moved, i.e. in mode 5 the highest maximum speed is available and in mode 1 the lowest (fig. 27).

#### Setting the mode

Press these buttons to change the mode, i.e. increase (right) or reduce (left) the maximum available speed. The speed is controlled by movement of the joystick (fig. 28).



#### Horn

Press this button to sound the horn (fig. 29).

## 6.2. Description of functions of VR-2 control

### Using for the first time

Sit on the seat of your Bischoff & Bischoff wheelchair and place your feet on the footplates.

Make sure that you are sitting as far back in the seat as possible. Make sure that loose clothing such as, for example, scarves, shirts, etc. do not hang out of the wheelchair to the side where they could get caught in the wheels.

Switch your wheelchair on using the on/off button (Fig. 26). You can now operate the control unit as described in this chapter.

### Joystick

The main function of the joystick is to control the speed and direction of the wheelchair. The further you push the joystick from the central position, the faster the wheelchair will go. When you move the joystick back, the brakes are applied automatically (see Chapter 2.5).

If the wheelchair is fitted with electronic adjustment options, the joystick can also be used for selection and adjustment, see Section 6.3 for further details.

### Control panel


All electronic functions for the wheelchair are operated using the control panel. The on/off button (Fig. 26) powers the electronics for the control system, which in turn powers the motors of the wheelchair. Do not use the on/off button to stop the wheelchair, unless it is an emergency. Doing this reduces the service life of the drive components. The battery display indicates that the wheelchair is switched on.

The functions of the control panel are described in Chapter 6.1.

### Charging and programming socket

Using the charging and programming socket is described in Chapter 6.14.

### 6.3. Functional description of the VR2 control unit for electronic adjustment options (optional)

#### **Controlling the adjustment options**

By pressing the adjustment function button (Fig. 33) the electronic adjustment options are activated. By pressing the button once the LED display below the button is activated.

#### **Wheelchairs with one adjustment option**

Pressing the adjustment function button (Fig. 33) once activates the electronic adjustment options. This is indicated by both LEDs below the adjustment function button lighting up.

By deflecting the joystick forwards or backwards, you can control the adjustment direction for this channel.

To return to drive mode, either press the adjustment function button (Fig. 33) or one of the two the speed level buttons (Fig. 28).

#### **Wheelchairs with two adjustment options**

By pressing the adjustment function button (Fig. 33) the electronic adjustment options are activated. By pressing the button once, the left LED lights up, and then by deflecting the joystick forwards or backwards you can control the adjustment direction for this channel.


Select between the two adjustment options by deflecting the joystick to the left or right. If the other adjustment option is selected, the relevant LED will light up. LEFT for adjustment option 1 and RIGHT for adjustment option 2.

To return to drive mode, press the adjustment function button (Fig. 33) again or press one of the two speed level buttons. (Fig. 28).

## 6.4. Functional description of the VR2 attendant control unit

As an option the electric wheelchair can be fitted with an attendant control unit.

- A** Toggle switch for operating mode
- B** Operating display
- C** Adjustment function button
- D** Display for activated adjustment option
- E** Speed level adjustment button
- F** Speed level display

### **Toggle switch for operating mode**

Use the toggle switch (A) to select the operating mode, whether the wheelchair will be controlled by the attendant control unit or by the user. The selected operating mode will be shown through the LEDs (B). The red LED lights up when the wheelchair is being controlled by the user. The green LED lights up when the wheelchair can be controlled by the attendant.

### **Adjustment function button**

All VR2 attendant control units are fitted with an adjustment function button (C) as standard. If the VR2 control unit is not programmed with any adjustment options, this button has no function.

However if the VR2 control unit is programmed with 1 or 2 adjustment options, the functions for this button are described in the relevant section below.

### **Wheelchairs with one adjustment option**

Pressing the adjustment function button (C) once will activate the electronic adjustment options.

This is indicated by both LEDs above the adjustment function button lighting up (D). By deflecting the joystick forwards or backwards, you can control the adjustment direction for this channel.

To return to drive mode, either press the adjustment function button (C) or the speed level button (E).

### **Wheelchairs with two adjustment options**

By pressing the adjustment function button (C) the electronic adjustment options are activated. By pressing the button once, the left LED lights up, and then by deflecting the joystick forwards or backwards you can control the adjustment direction for this channel.

Select between the two adjustment options by deflecting the joystick to the left or right. If the other adjustment option is selected, the relevant LED will light up. LEFT for adjustment option 1 and RIGHT for adjustment option 2.

To return to drive mode, press the adjustment function button (C) again or press one of the two speed level buttons. (Fig. 28).

### **Speed levels display**

Indicates the adjustment of the speed level for the wheelchair, if the attendant control unit is activated. There are five settings - Setting 1 is the lowest speed and setting 5 is the highest possible speed.

### **Button for speed level adjustment**

Use this button (E) to adjust the speed level of the wheelchair, if the attendant control unit is activated. By pressing this button, the maximum speed of the relevant drive speed will be increased by one level, until you reach level 5. Pressing again will reset the speed back to level 1. If the control unit is programmed for operation with drive profiles, only the speed level within the selected profile will be adjusted using the button for speed level adjustment of the positioning system.

## **6.5. Locking and unlocking of the wheelchair**

The VR2 control system can be locked to prevent unauthorised use. The locking method is by a sequence of button presses and joystick movements as described below.



To lock the wheelchair

- With the control system switched on, press and hold the on/off button.
- After 1 second, the control system will beep. Now release the on/off button.
- Move the joystick forward until the control system beeps.
- Move the joystick backwards until the control system beeps.
- Release the joystick, you will hear a long beep.
- The wheelchair is now locked.

To unlock the wheelchair;

- Switch on the control system with the on/off button. The maximum speed/profile display moves up and down.
- Move the joystick forward until the control system beeps.
- Move the joystick backwards until the control system beeps.
- Release the joystick, you will hear a long beep.
- The wheelchair is now unlocked.

### 6.6. Description of R-Net control

### 6.7. Functional description of the R-NET control unit for electrical adjustment options (optional)

The functions of the control panel are the same as for the VR2 control unit.



#### Commissioning

Sit in the seat of your Bischoff & Bischoff wheelchair and rest your legs on the footplates.

Make sure that you are sitting as far back in the seat as possible.

Make sure that no loose clothing such as scarves, shirts etc. is hanging out of the side of the wheelchair and could get caught in the wheels. Switch on your wheelchair by pressing the on/off switch (Fig. 26). You can now activate the control of the adjustment options by pressing the Mode button (Fig. 37).

#### Display panel for adjustment options

The user can activate the adjustment functions with the Mode button. The lamps for the relevant LED for an adjustment option (left footrest, right footrest, seat, backrest) show which adjustment function is active at that time.


### Adjustment functions

Once the adjustment functions have been activated using the Mode button, move the joystick (Fig. 35) left or right, to select between the possible adjustment functions. These are shown through the LEDs on the display panel (Fig. 36).

By deflecting the joystick forwards or backwards, you can control the relevant adjustment direction for the individual channels. In order to return to drive mode, press the Mode button again or one of the buttons for speed level adjustment.

## 6.8. Functional description of the R-NET attendant control unit

As an option, an attendant control unit is available so that an attendant is able to control the wheelchair.

- A** Toggle switch for operating mode
- B** Operating display
- C** Adjustment function button
- D** Display for active drive and speed level
- E** Speed level adjustment button

### Toggle switch for operating mode

Use the toggle switch (A) (Fig. 39) select the operating mode, whether the wheelchair will be controlled by the attendant control unit or by the user. The selected operating mode will be shown through the LEDs (B). The red LED lights up when the wheelchair is being controlled by the user. The green LED lights up when the wheelchair can be controlled by the attendant.

### Adjustment function button

All R-NET attendant control units are fitted with an adjustment function button as standard. If the R-NET control unit is not programmed with any adjustment options, this button has no function.

However if the R-NET control unit is programmed with 1 or up to 4 adjustment options, the functions for this button are described in the relevant section below.

### Wheelchairs with up to 4 adjustment options

Pressing the adjustment function button (C) once will activate the setup mode for the electronic adjustment options.

By deflecting the joystick forwards or backwards, you can control the relevant adjustment direction for the individual channels. Select between the adjustment options by deflecting the joystick to the left or right. When you change the selected adjustment option, this is shown on the control panel of the user's control unit.

To return to drive mode, either press the adjustment function button (C) or the speed level button (E).

### Speed levels display

Indicates the speed level adjustment for the wheelchair, if the attendant control unit is activated. There are five settings - setting 1 is the lowest speed and setting 5 is the highest speed. The speed level is adjusting using the (E) button.

**Button for speed level adjustment**

Use this button (E) to adjust the speed level of the wheelchair, if the attendant control unit is activated. By pressing this button, the maximum possible speed of the relevant drive speed will be increased by one level, until you reach level 5. Pressing again will reset the speed back to level 1. If the control unit is programmed for operation with drive profiles, only the speed level within the selected profile will be adjusted using the button for speed level adjustment of the positioning system.

**6.9. Locking and unlocking of the wheelchair**

The R-Net control system can be locked to prevent unauthorised use. The locking method is by a sequence of button presses and joystick movements as described below.

To lock the wheelchair

- With the control system switched on, press and hold the on/off button.
- After 1 second, the control system will beep. Now release the on/off button.
- Move the joystick forward until the control system beeps.
- Move the joystick backwards until the control system beeps.
- Release the joystick, you will hear a long beep.
- The wheelchair is now locked.

To unlock the wheelchair;

- Switch on the control system with the on/off button. The maximum speed/profile display moves up and down.
- Move the joystick forward until the control system beeps.
- Move the joystick backwards until the control system beeps.
- Release the joystick, you will hear a long beep.
- The wheelchair is now unlocked.

**6.10. Adjustments to the Control unit**

Changes to the programming of the control unit must not be carried out and voids the product conformity and road approval. Changes to the programming have an impact on speed, acceleration and brake delay. If carried out incorrectly, this can lead to fatal injuries to the user and/or a third party.

**6.11. Getting in and out of the chair****Entering and leaving from the side**

- Move the wheelchair as close as possible to the present or future seat so they are side by side.
- If you wish to move from another wheelchair or room chair, lock it by setting the parking brakes.
- Lock the Taiga to prevent it from moving accidentally by pressing the On/Off button to switch off the control (all displays off).
- You must also activate the mechanical locking brake when the drives are disengaged.
- Remove the arm rest and side section.
- Fold the foot plates to the side or remove the leg rests.
- Now slide sideways onto the other seat.
- Make sure that you are sitting as far back in the seat as possible.
- Move the arm rest and leg rests to the original position.

**Entering and leaving from the front**

- First fold the foot plates up or remove the leg rests.
- Move the Taiga as close as possible to the present or future seat front to front.
- If you wish to move from another wheelchair or room chair, lock it by setting the parking brakes.
- Now slide onto the seat by turning your body around.
- Then move the arm rest and leg rests to the original position.

### 6.12. Leg supports

To make it as easy as possible to enter and leave your wheelchair, the foot plates can be folded up, or the leg rests can be swung to the side out of the way or removed.

- Open the lock to swing the leg rest to the side (fig. 40).
- The leg rest in the swung-out position can be completely removed (fig. 41).
- For electrically adjustable leg rests, the connection cables under the seat must initially be disconnected (fig. 42).
- When you want to replace the leg rest, proceed in reverse order.
- Move the leg rest to the forward position until you hear it lock into place.

### 6.13. Check before use

#### 6.13.1. Brake check

Every time you use the wheelchair check the brake system. Move forward slowly and release the joystick. The wheelchair must come to a stop in the normal manner.

When it is at a standstill the electromagnetic brake must lock with an audible click. If you detect unusual brake behaviour, inform your mobility dealer immediately and do not use the wheelchair.

The motor brake operates correctly if you cannot push the wheelchair when it is switched off. However, if it can be moved, the electromagnetic brake is faulty. Do not use the wheelchair if this is the case. Contact your mobility dealer immediately.

#### 6.13.2. Checking lights, indicators and horn

Before moving off in the wheelchair always check the operation of the lighting system, the indicators and horn. If you press the button and there is no response from the equipment, inform your mobility dealer immediately and do not use the wheelchair.

#### 6.13.3. Checking the tyres

Before moving off in the wheelchair always check that the tyres are undamaged and the tyre pressure is adequate. If the pressure is too low, the tyres must be pumped up.

Check the tyre tread depth at least once a month. It should be at least 1 mm. Defective or worn tyres must be replaced by the mobility dealer.

### 6.14. Batteries

The wheelchair is equipped with maintenance-free lead acid battery whose fluid level does not need to be checked.

If you do not use your wheelchair for an extended period, the batteries should be charged every two months, otherwise they may be destroyed by deep discharge. Make a habit of checking the battery status display on the control every time before using your wheelchair.

If the battery status display shows only the red LEDs, do not travel long distances and connect the battery charger to the wheelchair as soon as possible. We recommend charging the battery after every time the wheelchair is used.

### 6.14.1. Charging

Only use approved chargers to charge the batteries (type in accordance with the technical data or as delivered by the Bischoff & Bischoff dealer Fig. 44).

Please also follow the operating instructions for the charger!

Do not charge frozen batteries.

The batteries must be charged when the temperature is between 0° and 40°C.

When using the electric wheelchair for the very first time, we recommend driving until just 3 bars are visible in the display (see Chapter 6.1) and recharge the batteries at that point.

In the event that this charge level is displayed, you should recharge the batteries as soon as possible in order to avoid deep discharge.

Recharging requires approx. 12 - 16 hours.


When using the product after that, you can recharge the batteries at any battery charge level.

You can interrupt the charging procedure even if the batteries are not fully charged, if you need to continue driving, for example.

But where possible, always charge to 100% so that you have the greatest range.

Use only by Bischoff & Bischoff approved battery chargers to charge the batteries (see Technical Data). Follow the instructions for use of the battery charger (pic. 43).

#### Charging batteries:

- First switch the electric wheelchair off.
- Before you start using the charger, check that the power supply is 230 V.
- Connect the charger with the charging socket of the electric wheelchair and then with the power socket (230 V).
- The charging socket is located on the side below the ignition switch (Fig. 45).

- Once the charging procedure is complete, unplug the mains plug and disconnect the electric wheelchair from the charger.
- When not in use, the battery should be charged once a month in order to avoid deep discharge.
- Deep discharge shortens the service life of the battery significantly, or the battery may also be damaged by the deep discharge.
- Batteries should be stored upright, in dry, frost-free conditions.
- As soon as the range with fully-charged batteries starts to decrease significantly, you should replace the batteries.

Never use your electric wheelchair with empty or almost empty batteries. You may endanger yourself and others by stopping unexpectedly.

### 6.14.2. Replacing batteries


When working on the batteries be very careful, do not use metal tools, do not touch the battery terminals and follow the directions of the battery manufacturer.

Caution! This work should only be done by a trained technician. Danger of electrocution!

Use only batteries approved for this wheelchair (see technical data). Make sure that the wheelchair is switched off before starting to replace batteries.

- First remove the battery cover (fig. 46).
- First loosen the mounting and pull the battery tray out (fig. 47).
- Disconnect the battery cable. First the red terminal and then the black terminal. This does not require any tools.

- Lift the batteries carefully out by the handles.
- Install the new batteries in their correct position and connect the terminals. Note the correct colour-coding (red and black).
- Secure the battery cover with the Velcro strap and push the battery compartment back into the initial position. It latches in again automatically.

When the batteries have reached the end of their life, i.e. they can no longer be adequately charged, they must not be disposed of with household rubbish. Old batteries are returned to the dealer when new batteries are purchased and correctly disposed of.

If you wish to dispose of it yourself, please contact a local recycling company for the disposal regulations of your jurisdiction.

### 6.15. Driving the wheelchair

Before moving check that the electromagnetic brake is correctly locked. The safest way is to make sure the electrical wheelchair is switched OFF!

- Make sure that you are comfortably and securely seated.
- Press the On/Off button to switch on the Taiga.
- All electrical driving functions are now available.

#### **Driving functions**

Your dealer can adjust the values for speed, acceleration and deceleration to your requirements subject to the legal regulations.

For a start select a low mode and increase it while moving if required.

Use the joystick to start moving. The joystick is used simultaneously to increase the speed and to change the direction of movement. When the joystick is in central position the electro-magnetical brake is active, unless the joystick is not moved the wheelchair stands still.

Moving forward: To move forward in a straight line move the joystick cautiously forward. The further forward the joystick is moved, the faster your wheelchair gets. If you release the joystick the wheelchair will automatically be braked to a stop.

Reversing: To move in reverse mode, move the joystick cautiously back. The further back the joystick is moved, the faster your wheelchair gets. If you release the joystick the wheelchair will automatically be braked to a stop.

Turning corners: To change direction of travel, tilt the joystick to the left or right. This can also be done when the joystick is tilted forwards or backwards, i.e. while driving backwards or forwards.

Caution! Driving on inclines is only permitted if the seat and back rest are in an upright position.

Caution! Always place the seat and back rest in an upright position when driving the wheelchair.

### 6.16. Pushing the wheelchair

If you wish to push the wheelchair in some situations, the brakes must be unlocked. The wheelchair is easy to push in freewheel mode.

Caution! The electromagnetic brake is disabled in this mode.

- With the yellow rotary switch the actuators can be unlocked and also coupled again (fig. 48).
- To be able to push the electric wheelchair manually, turn the switch downwards. The switch will be automatically engage.
- Repeat this procedure with the other actuator. Now you can use the wheelchair in push mode.
- To engage the actuators turn the lever upwards, it will click into place automatically. Also repeat this procedure on the other side.
- A slight jerk tells you that the drive mechanism is now engaged and you can use the driving function.

Please note that operating the wheelchair in push mode by uncoupling the drive mechanism will deactivate the electro-magnetic brakes. The wheelchair should therefore only be used on level ground in this mode. Non-observance may lead to life-threatening situations.

Caution! Always use both drives in the same position (accident hazard).

Caution! Danger of fatal injury! The drives must only be unlocked on level ground.

### 6.17. Head rest (optional)

The wheelchair can be equipped with an optional head rest. The settings for the version with the comfort seat are described below; the settings of the standard version are analogous.

- You can easily adapt the height of the head rest by pulling out the push button (fig. 49) and moving the head rest to the desired height. The push button latches in automatically and audibly.
- The position of the head rest can be adjusted and adapted to your needs by loosening the lever screws. Then retighten the lever screws.

### 6.18. Kerb climber (optional)

The kerb climber (fig. 50) makes it easier for you to negotiate steps of up to 100 mm.

Using the kerb climber requires a bit of skill and practice. Practice using the kerb climber initially on manageable terrain with an accompanying person.

Follow all safety instructions of the wheelchair when using the kerb climber.

Caution: risk of falls! Always drive directly at the kerb. Approaching it at an angle can lead to accidents.

It is preferable to use ramps to surmount obstacles. Only use the kerb climber if there are no ramps available to you.

- Drive up to the kerb carefully and ensure that you can clear the height of the kerb with the leg rests.
- Be sure to only drive directly at the kerb.
- The kerb should only be surmounted in the slowest driving mode.
- Carefully make contact with the kerb using the kerb climber and slowly accelerate until the front wheels have surmounted the obstacle.
- Then continue to drive slowly until the rear wheels have also surmounted the obstacle.

### 6.19. Mechanical parking brake

The mechanical parking brake allows you to park the electric wheelchair even when the drives are decoupled and secure it against accidentally rolling away.

#### **Using the parking brake**


To use the brake just push the brake lever forward until it snaps (fig. 51). To loosen the brake, pull the lever backwards.

#### **Adjustment of the parking brake**

If you should notice any unusual braking behaviour you have to do the following:

- Check the tread depth. Pay attention whether there is even abrasion. If the tread depth is > than 1 mm, the tyres must be replaced by a qualified mobility dealer.

- Now check the adjustments of the brake. It must be adjusted in the way that there is a distance between the running surface of the wheel and the pivot bolt of 15 mm when brakes are open. The pivot bolt must be horizontal e.g. parallel to the ground.
- After loosening the Allen screw with a hexagon socket (5 mm) on the brake mount (fig. 52), you can position the brake correctly.
- The parking brakes effect only one driving wheel at a time.
- Once the adjustment process is completed, firmly tighten the loosened screw up again and carry out a brake test!
- The brake system should be adjusted only by a qualified mobility dealer.

## 7. Technical data

### Dimensions and weights

Max. load capacity :	170 kg
Seat width:	415 - 515 mm
Effective seat width :	415 - 615 mm (via armrest adjustment)
Seat height :	500 mm (without seat cushion)
Seat depth :	410 - 560 mm
Seat angle :	-2° bis +10° (manually) -2° bis +20° (electrically)
Seat back height :	530 mm (without cushion, standard) 630 - 670 mm (without cushion, comfort seat)
Seat back angle :	90° bis 120° (in 4x10° steps, manually) 90° bis 135° (electrically)
Overall length :	900 mm (without supports) 1200 mm (with supports)
Height of the armrests :	210 - 270 mm
Length of the armrests :	410 mm
Leg rest length :	440 - 540 mm (manually) 420 - 580 mm (electrically)
Angle of the leg rest :	75°
Overall width :	670 mm
Unloaded weight :	121 kg (incl. batteries)
Weight of the batteries :	21,4 kg
Drive wheels :	14" (3.00 - 8mm), PU
Steering wheels :	10" (3.00 - 4mm), PU

### Driving characteristics

Speed :	10 km/h (forwards)
Braking distance at 6 km/h:	930 mm
Climbing capacity :	max. 10°
Obstacle height :	max. 50 mm (without kerb climber)
Range :	35 km
Turning circle :	810 mm

### Electrical system

Main fuse :	50 A
Battery :	maintenance-free lead acid battery 2 x 12V / 80Ah
Battery charger :	24V DC 8A
Motor:	2 x 24V / 600W / 4400 rpm
Controller:	PG Drive VR-2, R-Net 90A
Lights:	headlights and rear lights according to the German traffic regulations [StVZO]

### Ambient conditions

Temperature :	-25 - 50 °C
Air humidity :	20 - 80 %

### Storage conditions

Temperature :	-40 - 65 °C
Air humidity :	45 - 60 %

### Equipment

Frame:	Aluminium, powder-coated
Seat and back frame:	Steel, painted
Pad materials:	flame-resistant, tested according to EN1021-1/2
Brake:	automatic braking system Electromagnetic engine brake according to the German traffic regulations [StVZO]

### Operating forces

Joystick:	approx. 1 N
Buttons on the control panel:	approx. 1 N
Disconnecting the drive mechanisms:	> 60 N

## 8. Maintenance tips for users

### 8.1. Cleaning and disinfection

The seat and back materials can be cleaned with warm water using a sponge and/or a soft brush. In case of stubborn dirt, a commercial cleaning agent can be added to the water.

Note: Never use aggressive cleaning agents such as solvents / scouring agents or hard brushes.

Plastic parts are best cleaned with a commercially available plastic cleaner. Please observe the product information of the cleaning agent.

- Plastic claddings are attacked by non-ionic surfactants and solvents, especially alcohols!
- Wipe the frame parts of the wheelchair with a damp cloth. For heavier soiling, also use a mild cleaning agent.
- Clean the wheels with a damp brush with plastic bristles (do not use a wire brush!).
- Do not use a high-pressure cleaner.
- Wash the push handles, arm pads, seat and back pads with a mild detergent.
- Upholstery parts can be machine washed at 40°C. Do not machine dry, only dry at room temperature.

### 8.2. Disinfection

- Before disinfection, the rollator must be thoroughly cleaned. The wheelchair can be disinfected with a household disinfectant. Use approved spray disinfectants or alcohol-based wipe-on disinfectants.
- In the case of multiple product users, a commercial disinfectant is prescribed.
- Use only wipe-on disinfectants approved by the Robert Koch Institute (RKI) and the Association for Applied Hygiene (Verband für Angewandte Hygiene e.V.) (e.g. Bacillol AF or Kohrsolin FF);

see manufacturer's recommendations for exposure time and concentration.

- Caution: Damaged upholstery parts cannot be sufficiently disinfected. These should be replaced.
- Under certain circumstances, the surfaces can be affected by disinfectants. Longer-term restrictions of the functionality can be the result.

Please note our enclosed overview of the most important care, hygiene and maintenance tasks.

Attention: Risk of infection! If there are several users, the wheelchair must be cleaned after each use.

Caution! Do not use high pressure cleaners or steam cleaners.

## 9. Maintenance note

We recommend that the wheelchair is inspected by a specialist dealer at least once a year.

In the event of faults or defects on the wheelchair, it must be handed over immediately to the specialist dealer for repair.

The repair and replacement of parts is carried out by the specialist dealer.

The necessary information and documents for repair and maintenance will be made available to the specialist dealer on request.

If you have any questions, please contact the Bischoff & Bischoff customer service.

You must check your wheelchair regularly according to the maintenance schedule and find it to be in perfect working order.

### 9.2.1. Maintenance advice for the specialist trade

Maintenance must be carried out in accordance with the Bischoff & Bischoff maintenance plan. A maintenance schedule will be provided upon request. A training programme for the specialised trade is also available.

You must check your wheelchair regularly according to the maintenance schedule and find it to be in perfect working order.

### 10. Repair

Repairs and maintenance work may only be carried out by authorised partners (medical supply stores or specialist dealers). Please call our customer service department to find a partner in your area. The customer service contact details can be found on the back of these operating instructions.

### 11. Parking

Do not park your power wheelchair in damp or humid conditions, as this will lead to mould growth and rapid deterioration of the upholstery parts. Also ensure that the drives are engaged (see chapter 6.15).

### 12. Lifting and holding points

Due to its weight, it is preferable to use ramps for loading the electric wheelchair. If you need to lift the wheelchair nevertheless, only pick it up by the frame (fig. 53). Do not grip the panel, electronics or other attachment parts. You can secure the electric wheelchair for transport with the four lugs (fig. 54 and fig. 55) in the front and back under the seat.

### 13. Transport

The wheelchair must not be used as a vehicle seat. It must not be used as such under any circumstances. It does not fulfil the requirements of ISO 7176-19.

can cause extremely serious injuries and death in the event of an accident.

You can recognise that your wheelchair is not approved as a seat in a motor vehicle for the transport of persons with reduced mobility (KMP) by the symbols on the type plate of your product (chapter 3.2).

**Wheelchair not suitable as a seat in the KMP (Fig. 56).**

**Wheelchair suitable as a seat in the KMP (Fig. 57).**

This Taiga wheelchair is not approved for the transport of disabled persons in accordance with ISO 10542-2 and must not be used as a vehicle seat.

### 14. Passing on the wheelchair

The wheelchair is suitable for reuse. Before it is passed on to another user, it must be serviced and hygienically prepared by a specialist dealer. When passing on the wheelchair, please remember to give this user manual and the address of your specialist dealer to the new user.

### 15. Storage and dispatch

If the wheelchair is to be stored or shipped, all plug-in and non-fixed parts, except the drive wheels, must be removed and packed individually in suitable cardboard boxes. The electrical circuit must be opened for longer storage periods. The individually packed parts can then be packed together in a larger box. It is advisable to keep the original packaging and store it for this purpose so that it is available when needed. Your wheelchair is then optimally protected against environmental influences during storage or transport.

When storing your wheel chair, please observe the information in the technical data. To prevent mould growth and damage to the upholstered parts, the storage location should be as dry as possible and not exposed to direct sunlight.

## 16. Disposal

If the wheelchair is no longer used and you wish to dispose of it, please contact your mobility dealer.

If you would like to dispose of the waste yourself, enquire with local waste disposal companies about the disposal regulations for your place of residence.

## 17. Warranty conditions

1. For products supplied by Bischoff & Bischoff, the period of warranty is 24 months from purchase date. The products are delivered free from production and material defects. In case of evidence of material or manufacturing defect, defective parts will be replaced free of charge.
2. This warranty excludes claims, if a repair or replacement of a product or a part is required for the following reasons:
  - a) normal wear and tear, which include but is not limited to the following parts where fitted; Batteries, motor brushes, handles, armrests, upholstery, tyres, brakes, caps, etc.
  - b) overloading the product such as, for example, by exceeding the maximum user weight or load.
  - c) the product or part has not been maintained or serviced in accordance with the Instruction manual, the care and hygiene instructions or the manufacturer's recommendations as shown in the maintenance instructions.
  - d) Accessories have been used which are not original accessories.
  - e) The product or part has been damaged through negligence, accident or incorrect use.
  - f) changes/modifications have been made to the product or parts, which deviate from the manufacturer's specifications.

g) Repairs have been carried out, before our customer service has been informed of the circumstances.

3. To enforce the warranty, please contact your dealer with the exact details of the nature of the difficulty. The product must be repaired by a dealer authorised by Bischoff & Bischoff.

4. For parts, which have been repaired or exchanged within the scope of the warranty, the warranty will be extended to the remaining warranty period for the product in accordance with point 1).

5. For original spare parts which have been fitted at the customer's expense, these will have a 12 months guarantee, following the fitting.

6. If the repair fails after a reasonable time, the customer can, at his discretion, request a reduction, compensation or the termination of the purchase contract.

7. The guarantee is subject to the law of the country in which the product was purchased from Bischoff & Bischoff & Bischoff.

Please also observe the Bischoff & Bischoff care, warranty, hygiene and maintenance instructions. Your dealer will be happy to provide these to you.

For replacement devices and repairs, the warranty period is one year.

Non-observance of the operating instructions as well as improper use, maintenance jobs and technical changes or adaptations without the permission of the Bischoff & Bischoff company, lead to guarantee loss and the loss of product liability in general.

Warranty-Certificate

**Model: Taiga**

Serial Number\*

\* (to be filled in by authorized dealer)

**Dealer:**

Date and stamp

(Please note our general business terms and conditions)

\*Information sign can be found on the side frame

**Note!**

No liability is accepted for printing errors, product and colour modifications, or changes to technical specifications. The colours of the product shown may vary slightly in reality. Figures may contain accessories. The current version of each of the operating manuals can be found in the download area of our website:

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