

Basic price: **Ejoy RD XL**

€ 5.871,-

HMV-No.: 18.50.04.3002 | Max. total capacity 170 kg

**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)

<b>Telescopic seating system</b>	€
Seat width mm <b>420 mm</b> 470 mm 520 mm effective seat width   +120 mm extendable side panels <sup>1</sup>	<input type="radio"/>
Seat cushion 50 mm high, max user weight 170kg   Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 170kg   black   for seat width   420 mm / Order No. 844092000   470 mm / Order No. 844092100   520 mm / Order No. 844092200	303,45
Seat depth mm 400 425 <b>450</b> 475 500 525 550	<input type="radio"/>
Seat height rear mm   without seat cushion 430 450 <b>470</b> 490 510 530	<input type="radio"/>
Seat angle 0° – 20°   manually adjustable	<input checked="" type="radio"/>
Seat angle 0°– 20°   electrically adjustable   only with control module JSM-A-L or JSM-L-LED   Order No. 8530961000	1.402,20
Adaptable back height 500 mm	<input checked="" type="radio"/>
Contour backrest cushion   black   for seat width   420 mm / Order No. 8440923000   470 mm / Order No. 8440924000   520 mm / Order No. 8440925000	303,45
Back support angle 90°– 120°   manually adjustable, 4 steps	<input checked="" type="radio"/>
Back support angle 100°– 142°   electrically adjustable <sup>2</sup>   only with control module JSM-A-L oder JSM-L-LED   Order No. 8440960000	1.402,20
<b>Side panel</b>	
Adjustable in height and lateral* Side panels extendable each side   +60 mm	<input checked="" type="radio"/>


<sup>1</sup>Effective seat width by extendable side panels.

<sup>2</sup>Headrest is mandatory by powered backrest recline 100°–142°.

▼ Please mark (bold = default setting)

<b>Legrest</b>	€
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable   for seat width 420 mm / Order No. 8440386000 470 mm / Order No. 8440387000 520 mm / Order No. 8440388000 (Length of footboard has to be the same size as seat width)	216,70
Manual elevating legrest   left / right   unit complete left or right	123,10
for seat width   420 mm / right (Order No. 1352350001)   470 mm / right (Order No. 1352350002)   520 mm / right (Order No. 1352350003)   420 mm / left (Order No. 1352351001)   470 mm / left (Order No. 1352351002)   520 mm / left (Order No. 1352351003)	
Powered elevating legrest left / right   unit	870,-
only with control unit JSM-A-L or JSM-L-LED   left / Order No. 8355361000   right / Order No. 8355360000	
Lower leg length: <sup>3</sup> <b>370</b> 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest   left / right   unit	179,60
left / Order No. 1352371001   right / Order No. 1352370001	

<sup>3</sup>The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.



# ORDER FORM

## EJOY RD XL Electrically powered wheelchair

**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting) €

<b>Motor und Battery</b>	
Motor 6 km/h, 350 W   HMV-No. 18.50.04.3002	<input checked="" type="checkbox"/>
Motor 10 km/h, 350 W   HMV-No.: none   all-in price left and right   Order No. l. 8440703000, r. 8440702000	827,-
Battery 2 x 12 V/56 Ah	<input checked="" type="checkbox"/>
Battery 2 x 12 V/80 Ah   upcharge per piece   Order No. 8530700000	288,30
Charger 24 V/5 A	<input checked="" type="checkbox"/>
Charger 24 V/8 A <sup>4</sup>   upcharge   Order No. 8130845000	103,-
<b>Lights</b>	
Light system incl. indicators   according to German street regulations	<input checked="" type="checkbox"/>
<b>Remote</b>	
Swing-away length adjustable remote bracket	<input checked="" type="checkbox"/>
VR2 JSM-L	<input checked="" type="checkbox"/>
VR2 JSM-A-L   up to 2 powered seat adjustments   Order No. 8355805000	364,-
VR2 attendant remote   incl. bracket   Order No. 8401829001	1.035,30
R-net JSM-L-LED   up to 4 powered seat adjustments   Order No. 8540846000	1.090,-
R-Net attendant remote   incl. bracket   Order No. 8401829002	1.423,10
Remote mounting position left <b>right</b>	<input type="radio"/>
Attendant remote mounting position left <b>right</b>	<input type="radio"/>
<b>Tyres</b>	
Drive wheels suspensions	<input checked="" type="checkbox"/>
Drive wheel 14", PU   abrasion resistant, black   pair   Order No. 8440500000	294,70
Drive wheel 14", air   abrasion resistant, black   pair   Order No. 8440505000	<input checked="" type="checkbox"/>

▼ Please mark (bold = default setting) €

Castor 9", PU   abrasion resistant, black   pair   Order No. 8440400000	315,70
Castor 9", air   abrasion resistant, black   pair   Order No. 8440405000	<input checked="" type="checkbox"/>
Anti-tipper	<input checked="" type="checkbox"/>
Fender   front and rear	<input checked="" type="checkbox"/>
<b>Accessories</b>	
Lap belt	<input checked="" type="checkbox"/>
Headrest   multidimensional adjustments incl. bracket   Order No. 8440917000	215,-
Calf strap   420 mm / Order No. 8440320000   470 mm / Order No. 8440321000   520 mm / Order No. 8440322000	47,80
Rear view mirror left / right   unit incl. bracket   left / Order No. mirror 8401889000, holder 8401809000   right / Order No. mirror 8401884000, holder 8401809000	140,30
Cane holder for 2 canes left right   unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder   Order No. 8166750000	43,60
<b>Frame colour</b>	
White	<input checked="" type="checkbox"/>

<sup>4</sup>When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.