



Company:

Commission:

Date:

Basic price: **Ejoy RD XXL / 200 kg** € 8.310,-
 With VR2 remote | max. total capacity 200 kg
 HMV-No.: 18.50.04.5007

Basic price: **Ejoy RD XXL / 250 kg** € 14.654,-
 With R-Net remote | max. total capacity 250 kg
 HMV-No.: 18.50.04.5007

▼ Please mark (bold = default setting)

Telescopic seating system	€
Seat width mm 580 mm 630 mm 680 mm effective seat width +120 mm extendable side panels ¹	<input type="radio"/>
Seat cushion 50 mm high, max user weight 250kg Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 250kg black for seat width 580 mm / Order No. 8540920000 630 mm / Order No. 8540921000 680 mm / Order No. 8540922000	303,45
Seat depth mm 400 425 450 475 500 525 550	<input type="radio"/>
Seat height rear mm without seat cushion 430 450 470 490 510 530	<input type="radio"/>
Seat angle 0° – 20° manually adjustable	<input checked="" type="radio"/>
Seat support angle 0° – 20° electrically adjustable only with control module JSM-A-L or JSM-L-LED Order No. 8530961000	1.402,20
Contour backrest cushion black for seat width 580 mm / Order No. 8540923000 630 mm / Order No. 8540924000 680 mm / Order No. 8540925000	303,45
Back support angle 90° – 120° manually adjustable, 4 steps	<input checked="" type="radio"/>
Back support angle 100° – 142° electrically adjustable ² only with remote module JSM-A-L oder JSM-L-LED Order No. 8440960000	1.402,20
Side panel	
Adjustable in height and lateral Side panels extendable each side +60 mm	<input checked="" type="radio"/>

▼ Please mark (bold = default setting)

Legrest	€
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable for seat width 580 mm / Order No. 8540386000 630 mm / Order No. 8540387000 680 mm / Order No. 8540388000 (Length of footboard has to be the same size as seat width)	<input checked="" type="radio"/>
Manual elevating legrest left / right unit complete left or right	123,10
right (Order No. 1352350003) left (Order No. 1352351003)	
Powered elevating legrest left / right unit	870,-
only with control unit JSM-A-L or JSM-L-LED left / Order No. 8355361000 right / Order No. 8355360000	
Lower leg length: ³ 370 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest left / right unit	179,60
left / Order No. 1352371001 right / Order No. 1352370001	
Motor und Battery	
Motor 6 km/h, 450 W HMV-No. 18.50.04.5007	<input checked="" type="radio"/>
Motor 10 km/h, 450 W HMV-No.: none all-in price left and right	827,-
Order No. l. 8440703000, r. 8440702000	
Battery 2 x 12 V/80 Ah	<input checked="" type="radio"/>
Charger 24 V/8 A	<input checked="" type="radio"/>

¹Effective seat width by extendable side panels.

²Headrest is mandatory by powered backrest recline 100°–142°.

³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

