

Basic price: **Ejoy RD**

€ 5.460,-

HMV-No.: 18.50.04.0222 | Max. total capacity 140 kg

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

Telescopic seating system	€
Seat width mm 420 mm 470 mm 520 mm effective seat width +50 mm extendable side panels ¹	<input type="radio"/>
Seat cushion 50 mm high, max user weight 140kg Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 140kg black for seat width 420 mm / Order No. 844092000 470 mm / Order No. 844092100 520 mm / Order No. 844092200	303,45
Seat depth mm 400 425 450 475 500 525 550	<input type="radio"/>
Seat height rear mm without seat cushion 430 450 470 490 510 530	<input type="radio"/>
Seat angle 0° – 20° manually adjustable	<input checked="" type="radio"/>
Seat angle 0°– 20° electrically adjustable only with control module JSM-A-L or JSM-L-LED Order No. 8530961000	1.402,20
Adaptable back height 500 mm	<input checked="" type="radio"/>
Contour backrest cushion black for seat width 420 mm / Order No. 8440923000 470 mm / Order No. 8440924000 520 mm / Order No. 8440925000	303,45
Back support angle 90°– 120° manually adjustable, 4 steps	<input checked="" type="radio"/>
Back support angle 100°– 142° electrically adjustable ² only with control module JSM-A-L oder JSM-L-LED Order No. 8440960000	1.402,20
Side panel	
Adjustable in height and lateral* Side panels extendable each side +25 mm	<input checked="" type="radio"/>

¹Effective seat width by extendable side panels.

²Headrest is mandatory by powered backrest recline 100°–142°.

▼ Please mark (bold = default setting)

Legrest	€
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable for seat width 420 mm / Order No. 8440386000 470 mm / Order No. 8440387000 520 mm / Order No. 8440388000 (Length of footboard has to be the same size as seat width)	216,70
Manual elevating legrest left / right unit complete left or right	123,10
for seat width 420 mm / right (Order No. 1352350001) 470 mm / right (Order No. 1352350002) 520 mm / right (Order No. 1352350003) 420 mm / left (Order No. 1352351001) 470 mm / left (Order No. 1352351002) 520 mm / left (Order No. 1352351003)	
Powered elevating legrest left / right unit	870,-
only with control unit JSM-A-L or JSM-L-LED left / Order No. 8355361000 right / Order No. 8355360000	
Lower leg length: ³ 370 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest left / right unit	179,60
left / Order No. 1352371001 right / Order No. 1352370001	


³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

Standard equipment

Free of charge option

1 v. 2



Company:

Commission:

Date:

▼ Please mark (bold = default setting) €

Motor und Battery	
Motor 6 km/h, 350 W HMV-No. 18.50.04.0222	●
Motor 10 km/h, 350 W HMV-No.: none all-in price left and right	827,-
Order No. l. 8440703000, r. 8440702000	
Battery 2 x 12 V/56 Ah	●
Battery 2 x 12 V/80 Ah upcharge per piece	288,30
Order No. 8530700000	
Charger 24 V/5 A	●
Charger 24 V/8 A ⁴ upcharge	103,-
Order No. 8130845000	
Lights	
Light system incl. indicators according to German street regulations	●
Remote	
Swing-away length adjustable remote bracket	●
VR2 JSM-L	●
VR2 JSM-A-L up to 2 powered seat adjustments	364,-
Order No. 8355805000	
VR2 attendant remote incl. bracket	1.035,30
Order No. 8401829001	
R-net JSM-L-LED up to 4 powered seat adjustments	1.090,-
Order No. 8540846000	
R-Net attendant remote incl. bracket	1.423,10
Order No. 8401829002	
Remote mounting position left right	○
Attendant remote mounting position left right	○
Tyres	
Drive wheels suspensions	●
Drive wheel 14", PU abrasion resistant, black pair	294,70
Order No. 8440500000	
Drive wheel 14", air abrasion resistant, black pair	●
Order No. 8440505000	

▼ Please mark (bold = default setting) €

Castor 9", PU abrasion resistant, black pair Order No. 8440400000	315,70
Castor 9", air abrasion resistant, black pair Order No. 8440405000	●
Anti-tipper	●
Fender front and rear	●
Accessories	
Lap belt	●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-
Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80
Rear view mirror left / right unit incl. bracket	140,30
left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror. 8401884000, holder 8401809000	
Cane holder for 2 canes left right unit	64,-
for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	
Licence plate holder Order No. 8166750000	43,60
Frame colour	
White	●

⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.