



# ORDER FORM

## PYRO NEXT Light weight wheelchair

**Company:**

Basic price: **Pyro Next** € 747.-

**Customer:**

Basic price: **Pyro Next with hub brake** € 935.-

**Date:**

Max. capacity 125 kg

▼ Please mark!

Seat	€
Seat width in mm 400 430 460 490	<input type="radio"/>
Seat height front in mm 470 490 510	<input type="radio"/>
Seat depth 425 mm	<input checked="" type="radio"/>
Seat cushion   Seat widths 400 - 490 mm	24,96
Seat upholstery, Nylon	<input checked="" type="radio"/>
<b>Back support</b>	
Back support height 430 mm   with lumbar curve	<input checked="" type="radio"/>
Back support upholstery, Nylon	<input checked="" type="radio"/>
Handgrip elevating, removable*   Pair	158,50
Head support*   incl. fixation	195,52
<b>Sidepanels</b>	
Swing-away sidepanels   with arm support	<input checked="" type="radio"/>
<b>Leg support</b>	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	<input checked="" type="radio"/>
Leg support elevating left right	111,90
Footplate with heel strap	<input checked="" type="radio"/>
<b>Tyres</b>	
Rear wheel 24" x 1 3/8", PU   incl. quick-release axle	<input checked="" type="radio"/>
Front wheel 8" x 1 1/4", PU	<input checked="" type="radio"/>
PU-Tyres   solid rubber	<input checked="" type="radio"/>
Adapter block	<input checked="" type="radio"/>
Spoke protection, transparent 24"	150,80

▼ Please mark!

Wheel locks	€
Knee lever wheel lock   push to lock	<input checked="" type="radio"/>
Brake lever extension	23,92
<b>Accessories</b>	
Togo power push aid incl. mounting set	2.544,46
Mounting aid for adaptation	197,08
Anti-tip device   Pair	91.-
Tilting aid	<input checked="" type="radio"/>
Safety belt	27,87
Crutch holder	24,65
Therapy tray	223,60
Oxygen/infusion holder and bracket	175,76
APINO Shopper	18,93
HERCULES padlock	11,28
Passive light	<input checked="" type="radio"/>
<b>Frame color</b>	
Silver Blue	<input type="radio"/>

\* Head support and elevating handgrips cannot be mounted together

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment     Optional, no extra charge