



Company:

Basic price: **S-Eco 2 without drum brake** € 724.-

Basic price: **S-Eco 2 with drum brake** € 983,-

Customer:

Max. capacity 125 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 370 400 430 460 490	<input type="radio"/>
520	45,50
Seat height 500 mm	<input checked="" type="radio"/>
Seat depth 420 mm	<input checked="" type="radio"/>
Seat cushion Seat width 370 - 490 mm	26,40
Seat cushion Seat cushion 520 mm	37,40
Seat upholstery, Nylon	<input checked="" type="radio"/>
Sidepanels	
Duo-Sidepanel long swing-away	<input checked="" type="radio"/>
Sidepanel elevating left right	119,10
Sidepanel elevating with foldable arm pads left right	135,90
Leg support	
Swing-away leg support	<input checked="" type="radio"/>
Lower leg support elevating left right piece	118,40
Amputee support left right piece	172,70
Footplate with heel strap	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8 "	<input checked="" type="radio"/>
Front wheels 7" x 30 mm	<input checked="" type="radio"/>
PU-tyres solid rubber	<input checked="" type="radio"/>
Spoke protection, transparent 24" pair	159,50
Wheel back pair	66,90
One arm drive only without drum brake left right	357,50
Hand rim rubber-coated pair	173,80

▼ Please mark!

Back	€
Back support height 400 mm	<input checked="" type="radio"/>
Handgrip elevating, removable pair *	167,70
Head support incl. fixation *	206,80
Back support upholstery, Nylon	<input checked="" type="radio"/>
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension piece	25,30
Wheel lock for one-arm drive mechanical left right	97,50
Accessories	
Togo power push aid incl. mounting set	2.676.-
Mounting aid for Togo adaptation	208,50
Anti-tip device pair	96,30
Tilting aid	<input checked="" type="radio"/>
Safety belt	29,50
Crutch holder	26,10
Therapy tray	236,50
Oxygen/infusion holder and bracket	185,90
APINO Shopper	35,50
Retrofit kit Transit	231,-
Passive Beleuchtung	<input checked="" type="radio"/>
Frame colour	
Anthracite metallic	<input checked="" type="radio"/>

* Head support and elevating handgrips cannot be mounted together

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge