



**Company:**

**Customer:**

**Date:**

Base Price: **Neo XXL** | VR2 € 9.676,-  
 Base Price: **Neo XXL** | R-net € 9.936,-  
 Base Price: **Neo XXL** | extended seat frame | VR2 € 11.184,-  
 Base Price: **Neo XXL** | extended seat frame | R-net € 11.444,-  
 Max. capacity 200 kg

▼ Please mark!

Seat	€
Seat width 450 mm   effective 400 - 510 mm	●
Seat width adjustable 520 - 600 mm	234,60
Seat height 500 mm	●
Seat depth adjustable 420 - 510 mm	●
Seat cushion	●
Seat angle adjustable (-2° bis 19°) electrically mandatory with additional option R-Net control	1.637,-
Seat angle adjustable (-2° bis 19°) mechanically	●
Standard seat   black	●
Comfort seat   mandatory with extended seating frame	642,30
Recaro seat unit   seat and back support not electrically adjustable	3.430,-
Back	
Back support height 440 mm	●
Back support angle adjustment (-5° bis 75°) electrically mandatory with additional option R-Net control	2.276,-
Back support angle adjustment (-5° bis 75°) mechanically	●
Back support adjustable by strap system	175,60
Head support   incl. fixation	164,70
Sidepanels	
Adjustable in height and lateral	●
Leg support	
Leg support adjustable mechanically	●
Leg support elevating with length adjustment, mechanically left right   mandatory with calf supportn   piece	554,80
Calf support for horizontally elevating leg support   piece	159,-
Electrically elevating leg support   mandatory with additional option R-Net control left right   piece	1.475,-
Footplate angle adjustable	●
Motor and Battery	
Max. range 40 km	●
Max. Speed 6 km/h	●
Battery 2x 12V / 60 Ah	●
Charger 24V / 8 A	●

▼ Please mark!

Lighting	€
Front and rear lights, flashing lights LED	●
Control device	
mounted on right side	●
mounted on left side	49,60
swing-away control device, adjustable in length mounted on left side mounted on right side	131,50
VR-2 control	●
VR-2 control device for assistant   complete with fixation	1.035,-
R-Net control device for assistant   complete with fixation	1.423,-
Software kit: Software, cable (VR-2)	572,-
Software kit: Software, cable (R-net)	577,40
Tyres	
Front and rear wheel single suspension	●
Rear wheels 3.00 - 8   Front wheels 3.00 - 4   pneumatic tyres	●
conversion solid tyres   RolkoFill	606,30
conversion pneumatic tyres   block profile	255,-
conversion solid tyres   block profile   RolkoFill	1.048,-
Rear mudguard   pair not in combination with block profile	341,70
Accessories	
Anti-tip device	●
Transport attachment	●
Safety belt	81,40
Crutch holder	64,-
Rear mirror   can only be mounted with the bracket	58,80
Bracket for rear mirror	81,40
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge