



ORDER FORM

S-ECO 300 S-ECO 300 XL Standard wheelchair

Company:

Customer:

Date:

Basic price: **S-Eco 300 without hub brake** € 599,-

Basic price: **S-Eco 300 with hub brake** € 798,-

Max. capacity 125 kg

Basic price: **S-Eco 300 XL without hub brake** € 890,-

Basic price: **S-Eco 300 XL with hub brake** € 1.089,-

Max. capacity 170 kg

▼ Please mark!

Seat	€
Seat width in mm	
370 400 430 460 490	<input type="radio"/>
520	38,-
Seat width in mm XL-Version	<input type="radio"/>
520 550 580	
Seat height front in mm	
510	<input checked="" type="radio"/>
470 490	45,50
Seat depth 420 mm	<input checked="" type="radio"/>
Seat cushion Seat widths 370 – 460 mm	24,-
Seat cushion Seat width 490 mm	24,-
Seat cushion Seat widths 520 – 550 mm	34,-
Seat cushion Seat width 580 mm	40,-
Seat upholstery, Nylon	<input checked="" type="radio"/>
Back support	
Back support height 400 mm	<input checked="" type="radio"/>
Handgrip elevating, removable* Pair	152,40
Head support* incl. fixation not for seat width 580 mm	188,-
Back support upholstery, Nylon	<input checked="" type="radio"/>
Sidepanels	
Sidepanel removable with arm support, swing-away, elevating	<input checked="" type="radio"/>
Sidepanel elevating left right	108,20
Sidepanel elevating with foldable arm pads left right	123,50
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	<input checked="" type="radio"/>
Lower leg support elevating left right	107,60
Amputee support left right	157,-
Footplate with heel strap	<input checked="" type="radio"/>

* Head support and elevating handgrips cannot be mounted together

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 200 x 50 mm	<input checked="" type="radio"/>
PU-Tyres	<input checked="" type="radio"/>
Wheel base extension	<input checked="" type="radio"/>
Adapter block	<input checked="" type="radio"/>
Spoke protection, transparent 24"	145,-
One arm drive only without drum brake left right	325,-
Hand rim rubber-coated Pair 22" 24"	110,-
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension Piece	23,-
Knee lever wheel lock for one-arm drive left right	88,60
Accessories	
Togo power push aid incl. mounting set only up to SW 550 mm	2.446,60
Mounting aid for adaptation	189,50
Long telescopic tubes Togo assembly seat width 550 mm**	
Anti-tip device Pair	87,50
Tilting aid	<input checked="" type="radio"/>
Safety belt	26,80
Crutch holder	23,70
Therapy tray only up to SW 550 mm	215,-
Oxygen/infusion holder and bracket	169,-
HERCULES padlock	10,85
Passive light	<input checked="" type="radio"/>
Frame color	
Silver	<input checked="" type="radio"/>

** For seat width 550 mm these additional telescopic tubes are required.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge