



Company:

Customer No.:

Consignment:

Date:

Offer

Price: **VIA**

€ 3.456,-

Model No.: **8600**

Order

Max. payload: **120 kg**

Seat unit	€
Seat Width 430 mm	●
Seat Depth 430 mm	●
Seat Height 520 mm	●
Backrest	
Backrest Height 560 mm	●
Push Bar	●
Armrest	
Armrest incl. arm-pad (320 x 55 mm L x W)	●
Wheels	
Castor 8" PU	●
Drive Wheel 12,5" pneumatic	●
Leg-support	
Flip-up foot plate	●
Lower-leg length 440mm (incl. cushion 490 mm)	●

Electronic	€
PG - VSI 50 A	●
Holder for attendant control*	51,50
Best. Nr. 5000017396	
Max. Speed	
6 km/h	●
Propulsion	
2x 180 W	●
Battery	
24 V / 20 Ah Lithium-Ion, 450 Wh	●
Charger 2,5A	●
Accessories	
Crutch holder	●
Frame color	
White	●

*The existing control unit can be attached to the push bar using this bracket.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge