



ORDER FORM

CENTURO S1 4-WHEEL SCOOTER

Company:

Basic price: **Centuro S1 | 6 km/h**
 HMV: 18.51.05.1044 | Max. payload 160 kg

€ 2.693.-

Customer:

Date:

▼ Please Mark !

Technical data	€
Range approx. 35 km	●
Ground clearance 80 mm	●
Obstacle height 50 mm	●
Climbing ability 12° 21%	●
Reversing width 2050 mm	●
Tyres 10" pneumatic	●
Sprung suspension front and rear	●
Brakes electromagnetic	●
Mechanical parking brake	●
Handlebars ergonomic	●
Adjustable steering column	●
Motor and electrical equipment	
Motor 6 km/h, 550 W	●
Battery 2 x 12 V/35 Ah	●
Charger 3 A	●
Cockpit with LED display	●
Horn	●
Lighting according to German road traffic regulations	●
Authorized for public transportation	●

▼ Bitte ankreuzen

Seat and back support	€
Kapitänssitz drehbar, höheneinstellbar - Rückenlehne abklappbar und winkeleinstellbar - Armlehnen abklappbar und winkeleinstellbar - Sitzschienen zum Einstellen des Abstands zur Lenksäule	●
Zubehör	
Shopping basket removable	●
Floormat	●
Rear mirror right	●
Rear mirror holder right	●
Rear mirror left* Order No. 8162813000	20,50
Rear mirror holder left Order No. 8162815000	3,-
Crutch holder for 2 Crutches** Mounting at the seat frame Order No. 8162881000	79,00
Rain cover Order No. 8162886000	105,70
Frame color	
Silver	●

*just mountable with rear mirror holder

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge