



ORDER FORM

CENTURO S4 4-WHEEL SCOOTER

Company:

Basic price: **Centuro S4** | 6 km/h

€ 3.568.-

Customer:

Basic price: **Centuro S4** | 10 km/h

€ 4.311.-

Date:

Max. payload 160 kg

▼ Please Mark !

Technical Data	€
Range approx. 35 km	●
Ground clearance 115 mm	●
Obstacle height 50 mm	●
Climbing ability 12° 21%	●
Reversing width 1690 mm	●
Tyres 12" pneumatic	●
Sprung suspension front and rear	●
Brakes electromagnetic	●
Mechanical parking brake	●
Handlebars ergonomic	●
Adjustable steering column	●
Motor and electrical equipment	
Motor 6 km/h, 550 W	●
Motor 10 km/h, 700 W	●
Battery 2 x 12 V/56 Ah	●
Battery 2 x 12 V/80 Ah Order No. 2 x 8530700000	(2 pcs.) 288,30
Battery charger 5 A	●
USB-Port	●
Cockpit with LED display	●
Horn	●
Lighting according to German road traffic regulations	●
Electronic speed reduction when cornering (only for 10 km/h)	●
Seat and back support	
Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	●
XL Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	165,-

▼ Please Mark !

Accessories	€
Shopping basket removable	●
Floormat	●
Rear mirror right	●
Rear mirror holder right	●
Rear mirror left* Order No. 8162813000	20,50
Rear mirror holder left Order No. 8162815000	3,-
Crutch holder for 1 Crutch Mounting on the multifunction holder Order No. 8162882000	140,50
Crutch holder for 2 Crutches Mounting at the seat frame Order No. 8162881000	79,00
Holder for cross-foldable rollators Mounting on the multifunction holder Order No. 8166761000	247,10
Holder for oxygen bottle Mounting on the multifunction holder Order No. 8164888000	53,10
Holder for transport box Mounting on the multifunction holder Order No. 8150982000	55,64
Transport box with lock Order No. 8150987000	227,60
License plate holder required for > 6 km/h ** Order No. 8166730000	27,70
Rain cover Order No. 8162886000	105,70
Frame color	
White Anthracite	○

*only mountable with rear mirror holder

** In Germany electric vehicles with a maximum speed of more than 6 km/h have to own a liability insurance. To prove this you have to put an insurance license plate on the rear side of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge