



# ORDER FORM

# CENTURO S7 | S7 MAX | S7 XXL

## 4-Wheel Scooter

**Company:**

**Customer:**

**Date:**

Basic price: **Centuro S7** € 4.714,-  
 Max. payload 160 kg

Basic price: **Centuro S7 MAX** € 4.849,-  
 Max. payload 160 kg

Basic price: **Centuro S7 XXL** € 5.247,-  
 Max. payload 250 kg

▼ Please Mark !

	S7	S7 MAX	S7 XXL
<b>Technical Data</b>			
Motor 15 km /h	●	X	X
Motor 10 km /h	X	●	●
Range approx. 45 km	●	●	X
Range approx. 35 km	X	X	●
Battery 2 x 12 V / 80 Ah		●	
Battery charger 8 A		●	
Tyres 13"   pneumatic		●	
Ground clearance 130 mm		●	
Obstacle height 100 mm		●	
Climbing ability 12°   21 %		●	
Reversing width 2300 mm		●	
<b>Standard equipment</b>			
Sprung suspension   front and rear		●	
Brakes electromagnetic		●	
Mechanical parking brake		●	
Handlebars ergonomic		●	
Adjustable steering column		●	
USB-Port		●	
Cockpit with LCD-Display		●	
Horn		●	
Rear mirror   right and left		●	
Lighting according to German road traffic regulations		●	
Electronic speed reduction when cornering		●	
<b>Seat and back support</b>			
Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	●	●	X

▼ Please Mark !

	S7	S7 MAX	S7 XXL
XL Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	€ 171,-*	€ 171,-*	●
<b>Accessories</b>			
Shopping basket   detachable		●	
Floor mat		●	
Cup holder		●	
Crutch holder for 1 Crutch Mounting on the multifunction holder Order no. 8162882000		€ 140,50	
Crutch holder for 2 Crutches Mounting at the seat frame Order no. 8162881000		€ 79,00	
Holder for cross-foldable rollators Mounting on the multifunction holder Order no. 8166761000		€ 247,10	
Holder for oxygen bottle Mounting on the multifunction holder Order no. 8164888000		€ 53,50	
Holder for transport box Mounting on the multifunction holder Order no. 8150982000		€ 55,60	
Transport box with lock Order no. 8150987000		€ 227,60	
License plate holder   required for >6 km/h ** Order no. 8166730000		€ 27,70	
Rain cover Order no. 8162886000		€ 105,70	
<b>Frame color</b>			
White    Anthracite		○	

\* Surcharge / Option price less than price of standard equipment.  
 This price differs from the unit price for the seat unit.

\*\* In Germany electric vehicles with a maximum speed of more than 6 km / h have to own a liability insurance . To prove this you have to put an insurance license plate on the back of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment    ○ Optional, no extra charge    X not an Option