



Company:

Basic price: **Protego without drum brake** € 2.587.-

Basic price: **Protego with drum brake** € 2.856.-

Customer:

Max. capacity 130 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 390 - 415 440 - 465 490 - 515	<input type="radio"/>
Seat height front 430 mm without cushion	<input checked="" type="radio"/>
Seat depth adjustable from 420 - 480 mm	<input checked="" type="radio"/>
Seat angle adjustable 0 - 15°	<input checked="" type="radio"/>
Medial support	106.-
Seat upholstery Dartex, night grey	<input checked="" type="radio"/>
Back	
Back support ergonomic shape with lateral support	<input checked="" type="radio"/>
Back support height 580 - 640 mm	<input checked="" type="radio"/>
Back support angle 0 - 44°	<input checked="" type="radio"/>
Trunk support swing away left right	108.40
Tool holder lateral trunk support	15.90
Back support adjustable by strap system	<input checked="" type="radio"/>
Back support upholstery Dartex, night grey	<input checked="" type="radio"/>
Handgrip elevating, removable	<input checked="" type="radio"/>
Sidepanels	
Sidepanels adjustable in height and depth, removable with arm support upholsteries Dartex upholstery	<input checked="" type="radio"/>
Hemiplegic arm support left right	341.70
Leg support	
Leg support elevating and angle adjustable with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" x 45 mm PU	<input checked="" type="radio"/>
Spoke protection, transparent 24"	165.80
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension	26.30
Accessories	
Togo power push aid incl. mounting set	2.783.-
Mounting aid for adaptation	216.80
Tilting aid	<input checked="" type="radio"/>
Safety belt	30.60
Therapy tray	245.90
Additional package Therapy tray with lateral trunk support	361.80
Oxygen/infusion holder 110 mm 160 mm	193.30
Passive light	<input checked="" type="radio"/>
Frame colour	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge