

# ORDER FORM

## PYRO LIGHT Light weight wheelchair

**Company:**

Basic price: **Pyro Light without drum brake** € 1.131.-

Basic price: **Pyro Light with drum brake** € 1.400.-

**Customer:**

Max. capacity 125 kg

**Date:**

▼ Please mark!

| Seat   | €                                |
|--|----------------------------------|
| Seat width in mm<br>380    410    430    450    480  | <input type="radio"/>            |
| 520  | 47.30                            |
| Seat height front in mm<br>490   | <input checked="" type="radio"/> |
| 420    440    450    470    520    540   | 52.10                            |
| Seat depth   adjustable from 420 - 440 mm  | <input checked="" type="radio"/> |
| Seat cushion   Seat width 380 - 450 mm<br>Sitzbreite 380 mm   Order no. 9800080000<br>Sitzbreite 410 mm   Order no. 9800081000<br>Sitzbreite 430 mm   Order no. 9800082000<br>Sitzbreite 450 mm   Order no. 9800083000 | 27.40                            |
| Seat cushion  <br>Seat width 480 mm   Order no. 9800084000   | 29.70                            |
| Seat cushion<br>Seat width 520 mm   Order no. 9800085000   | 38.90                            |
| Seat upholstery, Vinyl   | <input checked="" type="radio"/> |
| <b>Back support</b>  |                                  |
| Back support height 420 mm   adjustable 400 - 440 mm   | <input checked="" type="radio"/> |
| Back support with lumbar curve   | <input checked="" type="radio"/> |
| Back support adjustable by strap-system<br>Seat width 380 - 520 mm   | <input checked="" type="radio"/> |
| Back support without lumbar curve<br>Order no.   1330932000  | 50.10                            |
| Handgrip elevating*, removable   Pair<br>Without drum brake Order no.   1270806000<br>With drum brake Order no.   1270807000   | 174.40                           |
| Head support*   incl. fixation<br>für SB 380 – 430 mm Order no.   9802900000<br>für SB 450 – 520 mm Order no.   9802901000   | 215.-                            |
| Back support upholstery, Vinyl   | <input checked="" type="radio"/> |


\* Head support and handgrips elevating cannot be mounted together

▼ Please mark!

| Sidepanels  | €                                |
|---|----------------------------------|
| Combi-sidepanel   with arm support, swing-away, elevating   | <input checked="" type="radio"/> |
| Sidepanel elevating    left    right<br>links Order no.   1352805000<br>rechts Order no.   1352806000   | 123.80                           |
| <b>Leg support</b>  |                                  |
| Swing-away leg support with individual footplates   | <input checked="" type="radio"/> |
| Leg support elongation 100 mm   Pair<br>Order no.   13300333000   | 69.60                            |
| Lower leg support elevating    left    right<br>links SB 370 – 400 Order no.   1352351001<br>links SB 430 – 460 Order no.   1352351002<br>links SB 490 – 520 Order no.   1352351003 | 123.10                           |
| rechts SB 37-40 Order no.   1352350001<br>rechts SB 43-46 Order no.   1352350002<br>rechts SB 49-52 Order no.   1352350003  |                                  |
| Amputee support    left    right<br>left Order no.   1352371001<br>right Order no.   1352370001   | 179.60                           |
| Footplates with heel straps   | <input checked="" type="radio"/> |
| <b>Tyres</b>  |                                  |
| Rear wheels 24" x 1 3/8" PU   incl. quick-release axle  | <input checked="" type="radio"/> |
| Front wheels 6" x 1 1/4" Soft rubber  | <input checked="" type="radio"/> |
| Hand rim rubber-coated 24"   pair<br>Order no.   980006600  | 125.80                           |
| Spoke protection, transparent 24"<br>Order no.   980006400  | 165.80                           |
| Wheel base extension  | <input checked="" type="radio"/> |
| Front wheel boom<br>Order no.   1260842000  | 76.60                            |
| One arm drive   only without drum brake    left    right<br>links für SB 37-40 Order no.   5000017646<br>rechts für SB 43-52 Order no.   5000017645                                 | 649.-                            |
| continuation page 2   |                                  |

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment     Optional, no extra charge



# ORDER FORM

**PYRO LIGHT**  
Light weight wheelchair

**Company:**

**Customer:**

**Date:**

▼ Please mark!

▼ Please mark!

| <b>Wheel locks</b>   |                                  |
|--|----------------------------------|
| Knee lever wheel lock   push to lock   | <input checked="" type="radio"/> |
| Wheel lock lever extension<br>Order no.   0009677000   | 26.30                            |
| Knee lever wheel lock for one-arm drive      left<br>right<br>links Order no.   9000611000<br>rechts Order no.   9000610000  | 101.40                           |
| <b>Accessories</b>   |                                  |
| Togo power push aid incl. mounting set (not available at the moment)   | 2.783,-                          |
| Mounting aid for adaptation<br>Order no.   8010199000  | 216.80                           |
| Anti-tip device   Pair<br>Order no.   1270845000   | 100.10                           |
| Tilting aid  | <input checked="" type="radio"/> |
| Lap belt<br>Order no.   9800012000   | 30.60                            |
| Crutch holder<br>Order no.   9800030002  | 27.10                            |
| Therapy tray<br>S (380-430)      L (430-480)      XL (520)<br>S (380-430) Order no.   9800006000<br>L (430-480) Order no.   9800007000<br>XL (520)      Order no.   9800008000 | 245.90                           |
| Oxygen/infusion holder and bracket<br>Order no.   1010884000   | 193.30                           |
| Passive light  | <input checked="" type="radio"/> |
| <b>Frame colour</b>  |                                  |
| Silver      Blue   | <input type="radio"/>            |

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge