

Basic price: **Ejoy RD XXL** / 250 kg  
 | max. total capacity 250 kg  
 HMV-No.: 18.50.04.5007

€ 14.654,-

**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)


<b>Telescopic seating system</b>	€
Seat width mm <b>580 mm</b> 630 mm 680 mm effective seat width   +120 mm extendable side panels <sup>1</sup>	<input type="radio"/>
Seat cushion 50 mm high, max user weight 250kg   Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 250kg   black   for seat width   580 mm / Order No. 8540920000   630 mm / Order No. 8540921000   680 mm / Order No. 8540922000	303,45
Seat depth mm 400 425 <b>450</b> 475 500 525 550	<input type="radio"/>
Seat height rear mm   without seat cushion 430 450 <b>470</b> 490 510 530	<input type="radio"/>
Seat angle 0° – 20°   manually adjustable	<input checked="" type="radio"/>
Contour backrest cushion   black   for seat width   580 mm / Order No. 8540923000   630 mm / Order No. 8540924000   680 mm / Order No. 8540925000	303,45
Back support angle 90°– 120°   manually adjustable, 4 steps	<input checked="" type="radio"/>
<b>Side panel</b>	
Adjustable in height and lateral Side panels extendable each side   +60 mm	<input checked="" type="radio"/>

<sup>1</sup>Effective seat width by extendable side panels.

▼ Please mark (bold = default setting)

<b>Legrest</b>	€
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable   for seat width 580 mm / Order No. 8540386000 630 mm / Order No. 8540387000 680 mm / Order No. 8540388000 (Length of footboard has to be the same size as seat width)	<input checked="" type="radio"/>
Manual elevating legrest   left / right   unit complete left or right	123,10
right (Order No. 1352350003)   left (Order No. 1352351003)	
Lower leg length: <sup>2</sup> <b>370</b> 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest   left / right   unit	179,60
left / Order No. 1352371001   right / Order No. 1352370001	
<b>Motor und Battery</b>	
Motor 6 km/h, 450 W   HMV-No. 18.50.04.5007	<input checked="" type="radio"/>
Battery 2 x 12 V/80 Ah	<input checked="" type="radio"/>
Charger 24 V/8 A	<input checked="" type="radio"/>

<sup>2</sup>The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.



**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)	€
<b>Lights</b>	
Light system incl. indicators   according to German street regulations	<input checked="" type="radio"/>
<b>Remote</b>	
Swing-away length adjustable remote bracket	<input checked="" type="radio"/>
<b>Electrical adjustments options</b>	
Without electrical adjustments options	<input checked="" type="radio"/>
With 4 electrical adjustment options	1.090,-
Order No. 8445812000	
R-Net attendant remote   incl. bracket	1.423,10
Order No. 8401829002	
Seat support angle 0°– 20°   electrically adjustable	
only with control module JSM-A-L or JSM-L-LED	1.402,20
Order No. 8530961000	
Back support angle 100°– 142°   electrically adjustable <sup>3</sup>	
only with remote module JSM-A-L oder JSM-L-LED	1.402,20
Order No. 8440960000	
Powered elevating legrest left / right   unit	
only with control unit JSM-A-L or JSM-L-LED	870,-
left / Order No. 8355361000	
right / Order No. 8355360000	
Remote mounting position left <b>right</b>	<input type="radio"/>
Attendant remote mounting position left <b>right</b>	<input type="radio"/>
<b>Tyres</b>	
Drive wheels suspensions	<input checked="" type="radio"/>
Drive wheel 14", PU   abrasion resistant, black   pair	294,70
Order No. 8440500000	
Drive wheel 14", air   abrasion resistant, black   pair	<input checked="" type="radio"/>
Order No. 8440505000	

▼ Please mark (bold = default setting)	€
Castor 9", PU   abrasion resistant, black   pair	315,70
Order No. 8440400000	
Castor 9", air   abrasion resistant, black   pair	<input checked="" type="radio"/>
Order No. 8440405000	
Fender   front and rear	<input checked="" type="radio"/>
Anti-tipper	<input checked="" type="radio"/>
<b>Accessories</b>	
Lap belt	<input checked="" type="radio"/>
Headrest   multidimensional adjustments incl. bracket	215,-
Order No. 8440917000	
Calf strap	
580mm / Order No. 8540320000	
630 mm / Order No. 8540321000	47,80
680 mm / Order No. 8540322000	
Rear view mirror left / right	
unit incl. bracket	140,30
left / Order No. mirror. 8401889000, holder 8401809000	
right / Order No. Rücksp. 8401884000, holder 8401809000	
Cane holder for 2 canes left right   unit	64,-
for mechanical back support, Order No. 5000016048	
for electrical back support, Order No. 5000016070	
Licence plate holder	43,60
Order No. 8166750000	
<b>Frame colour</b>	
White	<input checked="" type="radio"/>

<sup>3</sup>Headrest is mandatory by powered backrest recline 100°–142°.