

Basic price: **Ejoy RD XL | 6 km/h** € 5.871,-
 HMV-No.: 18.50.04.3002 | Max. total capacity 170 kg

Basispreis: **Ejoy RD | 10 km/h³** € 6.698,-
 / Max. Belastbarkeit 170 kg

Company:

Commission:

Date:

▼ Please mark (bold = default setting) €

Telescopic seating system	
Seat width mm 420 mm 470 mm 520 mm	○
effective seat width + 50 mm extendable side panels ¹	
Seat cushion 50 mm high, max user weight 170kg Nylon	●
Contoured seat cushion max. user weight 170kg black for seat width 420 mm / Order No. 8440920000 470 mm / Order No. 8440921000 520 mm / Order No. 8440922000	303,45
Seat depth mm 400 425 450 475 500 525 550	○
Seat height rear mm without seat cushion 430 450 470 490 510 530	○
Seat angle 0° – 20° manually adjustable	●
Adaptable back height 500 mm	●
Contour backrest cushion black for seat width 420 mm / Order No. 8440923000 470 mm / Order No. 8440924000 520 mm / Order No. 8440925000	303,45
Back support angle 90°– 120° manually adjustable, 4 steps	●
Side panel	
Adjustable in height and lateral* Side panels extendable each side +25 mm	●
Legrest	
Swing away legrest with angle adjustable foot plate	●
Continuous footboard, swing away, angle adjustable for seat width 420 mm / Order No. 8440386000 470 mm / Order No. 8440387000 520 mm / Order No. 8440388000 (Length of footboard has to be the same size as seat width)	216,70

¹Effective seat width by extendable side panels.


▼ Please mark (bold = default setting) €

Manual elevating legrest left / right unit complete left or right	123,10
for seat width 420 mm / right (Order No. 1352350001) 470 mm / right (Order No. 1352350002) 520 mm / right (Order No. 1352350003) 420 mm / left (Order No. 1352351001) 470 mm / left (Order No. 1352351002) 520 mm / left (Order No. 1352351003)	
Lower leg length: ² 370 390 410 430 450 470 490 510 530	○
the lower leg length must not be greater than the seat height	
Amputee legrest left / right left / Order No. 1352371001 right / Order No. 1352370001	179,60
Battery	
Battery 2 x 12 V/56 Ah	●
Battery 2 x 12 V/80 Ah upcharge for 2 pieces Order No. 2 x 8530700000	288,30
Charger 24 V/5 A Charger 24 V/8 A ⁴ upcharge Order No. 8130845000	103,-
Lights	
Light system incl. indicators according to German street regulations	●
continue next page	

²The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

³In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

**Company:****Commission:****Date:**

▼ Please mark (bold = default setting) €

Remote	
Swing-away length adjustable remote bracket	●
Control without electrical adjustments	●
VR2 Control up to 2 powered adjustments Order No. 8445810000	364,-
VR2 Attendant remote incl. bracket Order No. 8401829001	1.035,30
VR2 Control up to 4 powered adjustments Order No. 8445812000	1.090,-
VR2 Attendant remote incl. bracket Order No. 8401829002	1.423,10
Seat angle 0°– 20° electrically adjustable only with control module JSM-A-L or JSM-L-LED Order No. 8530961000	1.402,20
Back support angle 100°– 142° electrically adjustable ⁵ only with control module JSM-A-L oder JSM-L-LED Order No. 8440960000	1.402,20
Powered elevating legrest left / right unit only with control unit JSM-A-L or JSM-L-LED left / Order No. 8355361000 right / Order No. 8355360000	870,-
Remote mounting position left right	○
Attendant remote mounting position left right	○
Tyres	
Drive wheels suspensions	●
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000	●
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70
Castor 9", PU abrasion resistant, black pair Order No. 8440400000	●
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50
Anti-tipper	●
Fender front and rear	●

▼ Please mark (bold = default setting) €

Accessories	
Lap belt	●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-
Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80
Rear view mirror left / right unit incl. bracket left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	140,30
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder ⁶ Order No. 8166750000	43,60
Frame colour	
White	●

⁶In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

⁵Headrest is mandatory by powered backrest recline 100°–142°.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment ○ Free of charge option

2 v. 2