

Basic price: **Ejoy RD | 6 km/h** € 5.460,-
 HMV-No.: 18.50.04.0222 | Max. total capacity 140 kg

Basic price: **Ejoy RD | 10 km/h³** € 6.287,-
 / Max. total capacity 140 kg

Company:

Commission:

Date:

▼ Please mark (bold = default setting) €

Telescopic seating system	
Seat width mm	
420 mm / Standard	
470 mm / Order No. seat 8440911000, backrest 8440931000	<input type="radio"/>
520 mm / Order No. seat 8440912000, backrest 8440932000	
effective seat width +50 mm extendable side panels ¹	
Seat cushion 50 mm high, max user weight 140kg Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 140kg black	
for seat width	
420 mm / Order No. 8440920000	303,45
470 mm / Order No. 8440921000	
520 mm / Order No. 8440922000	
Seat depth mm	
400 425 450 475 500 525 550	<input type="radio"/>
Seat height rear mm without seat cushion	
430 450 470 490 510 530	<input type="radio"/>
Seat angle 0° – 20° manually adjustable	<input checked="" type="radio"/>
Contour backrest cushion black	
for seat width	
420 mm / Order No. 8440923000	303,45
470 mm / Order No. 8440924000	
520 mm / Order No. 8440925000	
Back support angle 90° – 120°	<input checked="" type="radio"/>
manually adjustable, 4 steps	
Side panel	
Adjustable in height and lateral	<input checked="" type="radio"/>
Side panels extendable each side +25 mm	
Legrest	
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable	
for seat width	
420 mm / Order No. 8440386000	
470 mm / Order No. 8440387000	216,70
520 mm / Order No. 8440388000	
(Length of footboard has to be the same size as seat width)	

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting) €

Manual elevating legrest left / right	
unit complete left or right	123,10
for seat width	
420 mm / right (Order No. 1352350001)	
470 mm / right (Order No. 1352350002)	
520 mm / right (Order No. 1352350003)	
420 mm / left (Order No. 1352351001)	
470 mm / left (Order No. 1352351002)	
520 mm / left (Order No. 1352351003)	
Lower leg length: ²	
370 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest left / right	unit 179,60
left / Order No. 1352371001	
right / Order No. 1352370001	
Battery	
Battery 2 x 12 V/56 Ah	<input checked="" type="radio"/>
Battery 2 x 12 V/80 Ah	upcharge for 2 pieces 288,30
Order No. 2 x 8530700000	
Charger 24 V/5 A	<input checked="" type="radio"/>
Charger 24 V/8 A ⁴	upcharge 103,-
Order No. 8130845000	
Lights	
Light system incl. indicators according to German street regulations	<input checked="" type="radio"/>

²The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.


³In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

Standard equipment Free of charge option

1 v. 2



Company:

Commission:

Date:

▼ Please mark (bold = default setting) €

Remote			
Swing-away length adjustable remote bracket			●
Control without electrical adjustments			●
VR2 Control up to 2 powered adjustments		364,-	
Order No. 8442810000			
VR2 Attendant remote incl. bracket		1.035,30	
Order No. 8401829001			
R-netControl up to 4 powered adjustments		1.090,-	
Order No. 8445812000			
R-net Attendant remote incl. bracket		1.423,10	
Order No. 8401829002			
Seat angle 0°– 20° electrically adjustable		1.402,20	
Order No. 8530961000			
Back support angle 100°– 142° electrically adjustable ⁵		1.402,20	
Order No. 8440960000			
Powered elevating legrest	left / right	unit	870,-
left / Order No. 8355361000			
right / Order No. 8355360000			
Remote mounting position	left	right	○
Attendant remote mounting position	left	right	○
Tyres			
Drive wheels suspensions			●
Drive wheel 14", PU abrasion resistant, black pair			●
Order No. 8440500000			
Drive wheel 14", air abrasion resistant, black pair		294,70	
Order No. 8440505000			
Castor 9", PU abrasion resistant, black pair			●
Order No. 8440400000			
Castor 9", air abrasion resistant, black pair		216,50	
Order No. 8440405000			
Anti-tipper			●
Fender front and rear			●

▼ Please mark (bold = default setting) €

Accessories			
Lap belt			●
Headrest multidimensional adjustments incl. bracket		215,-	
Order No. 8440917000			
Calf strap			
420 mm / Order No. 8440320000			
470 mm / Order No. 8440321000			
520 mm / Order No. 8440322000			
Rear view mirror	left / right	unit incl. bracket	140,30
left / Order No. mirror 8401889000, holder 8401809000			
right / Order No. mirror. 8401884000, holder 8401809000			
Cane holder for 2 canes	left	right unit	64,-
for mechanical back support, Order No. 5000016048			
for electrical back support, Order No. 5000016070			
Licence plate holder			43,60
Order No. 8166750000			
Frame colour			
White			●

⁵Headrest is mandatory by powered backrest recline 100°–142°.

● Standard equipment ○ Free of charge option