

# ORDER FORM

## PYRO LIGHT OPTIMA XL

Light weight wheelchair

**Company:**

Basic price: **Pyro Light Optima XL without drum brake** € 1.746.-

Basic price: **Pyro Light Optima XL with drum brake** € 2.015.-

**Customer:**

Max. capacity 170 kg

**Date:**

▼ Please mark!

Seat	€
Seat width in mm 520    550*    580	○
Seat height front in mm 500	●
400   included front wheel 6", rear wheel 22"	52.10
425    450    475	52.10
Seat depth 450 mm (Standard) adjustable from 450 mm up to 500 mm in 25 mm steps	●
Seat cushion height 7 cm	
Sitzbreite 520 mm   Order no. 9800089000	34,-
Sitzbreite 550 mm   Order no. 9800097000	46,-
Sitzbreite 580 mm   Order no. 9800098000	34,-
Seat cushion height 5 cm	
Sitzbreite 520 mm   Order no. 9800085000	38,90
Sitzbreite 550 mm   Order no. 9800087000	46,-
Sitzbreite 580 mm   Order no. 9800088000	46,-
Seat cushion   Seat width 580 mm	
Seat upholstery, Nylon	●
<b>Back</b>	
Back support height 420 mm	●
Back support height adjustable 400 - 460 mm	●
Back support with lumbar curve	●
Back support adjustable by strap-system	●
Head support   incl. fixation   Not available for sw 580	
SB 520 mm    SB 550 mm   Order no. 9802901000	215,-
Handgrip elevating, removable	●
Back support upholstery, Nylon	●
<b>Sidepanels</b>	
Sidepanels adjustable in depth and height	●


▼ Please mark!

Leg support	€
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elongation 100 mm   piece   Order no. 1330333000	34.80
Lower leg support elevating	
left   Order no. 1352351003	123.10
right   Order no. 1352350003	
Amputee support                      left              right	179.60
Footplates adjustable in angle and in width	●
<b>Tyres</b>	
Rear wheels 24" x 1 3/8" PU   incl. quick-release axle	●
Front wheel 8"   200 x 50 mm, PU	●
Hand rim rubber-coated 24"   Pair   Order no. 9800066000	180.80
Spoke protection, transparent 24"   Order no. 9800054000	165.80
Wheel base extension   Order no. 1270841000	69.50
One arm drive   only without drum brake	
left für SB 520   Order no. 5000017646	
right für SB 550   Order no. 5000017645	649.-
for SB 580: not available	
<b>Wheel locks</b>	
Knee lever wheel lock   push to lock	●
Wheel lock lever extension   Order no. 0009677000	26.30
Knee lever wheel lock for one-arm drive	
left   Order no. 9000611000	101.40
right   Order no. 9000610000	
Accessories next page	

\* For seat width 550 mm these additional telescopic tubes are required.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment    ○ Optional, no extra charge



# ORDER FORM

## PYRO LIGHT OPTIMA XL

Light weight wheelchair

**Company:**

**Customer:**

**Date:**

▼ Please mark!

<b>Accessories</b>	
Togo power push aid incl. mounting set (not available at the moment)	...
only up to SW 550 mm	
Mounting aid for Togo adaptation	216.80
Order no. 8010199000	
Long telescopic tubes   Togo assembly seat width 550 mm*	203.-
Order no. 1x8010115000; 1x8010114000; 2x8010116000	
Tilting aid	<input checked="" type="radio"/>
Anti-tip device   Pair   Order no. 1270845000	100.10
Lap belt XL   Order no. 9800017000	35.30
Cane holder   Order no. 9800030002	27.10
Therapy tray	
SB 520 SB 550   Order no. 9800008000	245.90
SB 580: nicht verfügbar	
Passive light	<input checked="" type="radio"/>
<b>Frame colour</b>	
Silver	<input checked="" type="radio"/>

\* For seat width 550 mm these additional telescopic tubes are required.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge