



# ORDER FORM

## CENTURO S2 4-WHEEL SCOOTER

**Company:**

Basic price: **Centuro S2** | 6 km/h

€ 2.904,-

Basic price: **Centuro S2** | 10 km/h

€ 3.507,-

**Customer:**

Max. payload 160 kg

**Date:**

▼ Please Mark !

Technical Data	€
Range approx. 35 km	●
Ground clearance 80 mm	●
Obstacle height 50 mm	●
Climbing ability 6°   10%	●
Reversing width 2050 mm	●
Tyres 10"   pneumatic	●
Sprung suspension   front and rear	●
Brakes electromagnetic	●
Mechanical parking brake	●
Handlebars ergonomic	●
Adjustable steering column	●
<b>Motor and electrical equipment</b>	
Motor 6 km/h, 550 W	●
Motor 10 km/h, 600 W	●
Battery 2 x 12 V/55 Ah	●
Battery charger 5 A	●
USB-Port	●
Cockpit with LED display	●
Horn	●
Lighting according to German road traffic regulations	●
Authorized for public transportation	●
Electronic speed reduction when cornering (only for 10 km/h)	●
<b>Seat and back support</b>	
Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	●

▼ Please Mark !

Accessories	€
Shopping basket   removable	●
Floormat	●
Rear mirror   right	●
Rear mirror holder   right	●
Rear mirror   left*	20,50
Order No. 8162813000	
Rear mirror holder   left	3,-
Order No. 8162815000	
Crutch holder for 2 Crutches**	79,00
Mounting at the seat frame	
Order No. 8162881000	
License plate holder   required for > 6 km/h ***	27,70
Order No. 8166730000	
Rain cover	105,70
Order No. 8162886000	
<b>Frame color</b>	
White      Anthracite	○

\*just mountable with rear mirror holder

\*\* The overall length of the scooter is extended by attaching the crutch holder.

\*\*\*In Germany electric vehicles with a maximum speed of more than 6 km / h have to own a liability insurance . To prove this you have to put an insurance license plate on the back of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge