

**Company:**

**Commission:**

**Date:**

Basic price: **Ejoy FD XXL / 6 km/h** € 15.550,-  
With VR2 remote, inclusive 2 electrical adjustments

| max. total capacity 200 kg

HMV-No.: 18.50.04.5008

Basic price: **Ejoy FD XXL / 6 km/h** € 16.790,-  
With R-Net remote, incl. 4 electrical adjustments

| max. total capacity 250 kg

HMV-No.: 18.50.04.5008

▼ Please mark (bold = default setting)

Telescopic seating system	€
Seat width mm <b>580 mm</b> 630 mm 680 mm effective seat width   +120 mm extendable side panels <sup>1</sup>	<input type="radio"/>
Seat cushion 50 mm high, max user weight 250kg   Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 250kg   black   for seat width   580 mm / Order No. 854092000   630 mm / Order No. 8540921000   680 mm / Order No. 8540922000	303,45
Seat depth mm 400 425 <b>450</b> 475 500 525 550	<input type="radio"/>
Seat height rear mm   without seat cushion 430 450 <b>470</b> 490 510 530	<input type="radio"/>
Seat angle 0° – 20°   manually adjustable	<input checked="" type="radio"/>
Seat support angle 0° – 20°   electrically adjustable   only with control module JSM-A-L or JSM-L-LED   Order No. 8530961000	1.402,20
Contour backrest cushion   black   for seat width   580 mm / Order No. 8540923000   630 mm / Order No. 8540924000   680 mm / Order No. 8540925000	303,45
Back support angle 90° – 120°   manually adjustable, 4 steps	<input checked="" type="radio"/>
Back support angle 100° – 142°   electrically adjustable <sup>2</sup>   only with remote module JSM-A-L oder JSM-L-LED   Order No. 8440960000	1.402,20
<b>Side panel</b>	
Adjustable in height and lateral Side panels extendable each side   +60 mm	<input checked="" type="radio"/>

<sup>1</sup>Effective seat width by extendable side panels.


<sup>2</sup>Headrest is mandatory by powered backrest recline 100°–142°.

▼ Please mark (bold = default setting)

Legrest	€
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable   for seat width 580 mm / Order No. 8540386000 630 mm / Order No. 8540387000 680 mm / Order No. 8540388000 (Length of footboard has to be the same size as seat width)	<input checked="" type="radio"/>
Manual elevating legrest   left / right   unit complete left or right	123,10
right (Order No. 1352350003)   left (Order No. 1352351003)	
Powered elevating legrest left / right   unit	870,-
only with control unit JSM-A-L or JSM-L-LED   left / Order No. 8355361000   right / Order No. 8355360000	
Lower leg length: <sup>3</sup> <b>370</b> 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest   left / right   unit	179,60
left / Order No. 1352371001   right / Order No. 1352370001	
<b>Battery</b>	
Motor 10 km/h, 450 W   HMV-No.: none   all-in price left and right	827,-
Order No. l. 8440703000, r. 8440702000	
Battery 2 x 12 V/80 Ah	<input checked="" type="radio"/>
Charger 24 V/8 A	<input checked="" type="radio"/>

<sup>3</sup>The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

<sup>4</sup>When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.



**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)

<b>Lights</b>		€
Light system incl. indicators   according to German street regulations		<input checked="" type="radio"/>
<b>Remote</b>		
Swing-away length adjustable remote bracket		<input checked="" type="radio"/>
VR2 attendant remote   incl. bracket   Order No. 8401829001	1.035,30	
R-net Control with 4 electrical adjustment options   included as standard for 250 kg model   Order No. 8445812000	1.090,-	
Attendant remote   incl. bracket   Order No. 8401829002	1.423,10	
Remote mounting position	left	<b>right</b> <input type="radio"/>
Attendant remote mounting position	left	<b>right</b> <input type="radio"/>
<b>Tyres</b>		
Drive wheels suspensions		<input checked="" type="radio"/>
Drive wheel 14", PU   abrasion resistant, black   pair   Order No. 8440500000	294,70	
Drive wheel 14", air   abrasion resistant, black   pair   Order No. 8440505000		<input checked="" type="radio"/>
Castor 9", PU   abrasion resistant, black   pair   Order No. 8440400000	315,70	
Castor 9", air   abrasion resistant, black   pair   Order No. 8440405000		<input checked="" type="radio"/>
Fender   front and rear		<input checked="" type="radio"/>
Anti-tipper		<input checked="" type="radio"/>

▼ Please mark (bold = default setting)

<b>Accessories</b>		€
Lap belt		<input checked="" type="radio"/>
Headrest   multidimensional adjustments incl. bracket   Order No. 8440917000	215,-	
Calf strap   580mm / Order No. 8540320000   630 mm / Order No. 8540321000   680 mm / Order No. 8540322000	47,80	
Rear view mirror	left / right	
	unit incl. bracket	140,30
	left / Order No. mirror. 8401889000, holder 8401809000   right / Order No. Rücksp. 8401884000, holder 8401809000	
Cane holder for 2 canes	left right   unit	64,-
	for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	
Licence plate holder   Order No. 8166750000		43,60
<b>Frame colour</b>		
White		<input checked="" type="radio"/>