

Basic price: **Ejoy FD | 6 km/h** € 5.875,-  
 HMV-No.: 18.50.04.0223 | Max. total capacity 140 kg

Basic price: **Ejoy FD | 10 km/h<sup>2</sup>** € 6.665,-  
 Max. total capacity 140 kg

**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)

<b>Telescopic seating system</b>	€
Seat width mm <b>420 mm</b> 470 mm 520 mm	<input type="radio"/>
effective seat width   +50 mm extendable side panels <sup>1</sup>	
Seat cushion 50 mm high, max user weight 140kg   Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 140kg   black   for seat width   420 mm / Order No. 8440920000   470 mm / Order No. 8440921000   520 mm / Order No. 8440922000	303,45
Seat depth mm 400 425 <b>450</b> 475 500 525 550	<input type="radio"/>
Seat height rear mm   without seat cushion 430 450 <b>470</b> 490 510 530	<input type="radio"/>
Seat angle 0° – 20°   manually adjustable	<input checked="" type="radio"/>
Adaptable back height 500 mm	<input checked="" type="radio"/>
Contour backrest cushion   black   for seat width   420 mm / Order No. 8440923000   470 mm / Order No. 8440924000   520 mm / Order No. 8440925000	303,45
Back support angle 90°– 120°   manually adjustable, 4 steps	<input checked="" type="radio"/>
<b>Side panel</b>	
Adjustable in height and lateral* Side panels extendable each side   +25 mm	<input checked="" type="radio"/>
<b>Legrest</b>	
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable   for seat width 420 mm / Order No. 8440386000 470 mm / Order No. 8440387000 520 mm / Order No. 8440388000 (Length of footboard has to be the same size as seat width)	216,70

<sup>1</sup>Effective seat width by extendable side panels.


▼ Please mark (bold = default setting)

Manual elevating legrest   left / right   unit complete left or right	123,10
for seat width   420 mm / right (Order No. 1352350001)   470 mm / right (Order No. 1352350002)   520 mm / right (Order No. 1352350003)   420 mm / left (Order No. 1352351001)   470 mm / left (Order No. 1352351002)   520 mm / left (Order No. 1352351003)	
Lower leg length: <sup>3</sup> <b>370</b> 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest   left / right   unit	179,60
left / Order No. 1352371001   right / Order No. 1352370001	
<b>Battery</b>	
Battery 2 x 12 V/56 Ah	<input checked="" type="radio"/>
Battery 2 x 12 V/80 Ah   upcharge per 2 pieces	288,30
Order No. 2 x 8530700000	
Charger 24 V/5 A	<input checked="" type="radio"/>
Charger 24 V/8 A <sup>4</sup>   upcharge	103,-
Order No. 8130845000	
<b>Lights</b>	
Light system incl. indicators   according to German street regulations	<input checked="" type="radio"/>

<sup>2</sup>In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

<sup>3</sup>The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

<sup>4</sup>When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.



**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting) €

<b>Remote</b>			
Swing-away length adjustable remote bracket			●
Control without electrical adjustments			●
VR2 Control up to 2 powered adjustments		364,-	
Order No. 8355805000			
VR2 Attendant remote   incl. bracket		1.035,30	
Order No. 8401829001			
R-net Control up to 4 powered adjustments		1.090,-	
Order No. 8540846000			
R-net Attendant remote   incl. bracket		1.423,10	
Order No. 8401829002			
Seat angle 0°– 20°   electrically adjustable		1.402,20	
Order No. 8530961000			
Back support angle 100°– 142°   electrically adjustable <sup>5</sup>		1.402,20	
Order No. 8440960000			
Powered elevating legrest	left / right	unit	870,-
left / Order No. 8355361000			
right / Order No. 8355360000			
Remote mounting position	left	<b>right</b>	○
Attendant remote mounting position	left	<b>right</b>	○
<b>Tyres</b>			
Drive wheels suspensions			●
Drive wheel 14", PU   abrasion resistant, black   pair			●
Order No. 8440500000			
Drive wheel 14", air   abrasion resistant, black   pair		294,70	
Order No. 8440505000			
Castor 9", PU   abrasion resistant, black   pair			●
Order No. 8440400000			
Castor 9", air   abrasion resistant, black   pair		216,50	
Order No. 8440405000			
Anti-tipper			●
Fender   front and rear			●

<sup>5</sup>Headrest is mandatory by powered backrest recline 100°–142°.

▼ Please mark (bold = default setting) €

<b>Accessories</b>			
Lap belt			●
Headrest   multidimensional adjustments incl. bracket		215,-	
Order No. 8440917000			
Calf strap			
420 mm / Order No. 8440320000			
470 mm / Order No. 8440321000			
520 mm / Order No. 8440322000			
Rear view mirror	left / right	unit incl. bracket	140,30
left / Order No. mirror 8401889000, holder 8401809000			
right / Order No. mirror. 8401884000, holder 8401809000			
Cane holder for 2 canes	left	right   unit	64,-
for mechanical back support, Order No. 5000016048			
for electrical back support, Order No. 5000016070			
Licence plate holder <sup>6</sup>			43,60
Order No. 8166750000			
<b>Frame colour</b>			
White			●

<sup>6</sup>In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

● Standard equipment    ○ Free of charge option

Version 2025.03.03 Prices plus VAT at the legally applicable rate, subject to technical changes.