

ORDER FORM

EJOY MD XL

Electrically powered wheelchair

Basic price: **Ejoy MD XL | 6 km/h** € 6.750,-
Max. total capacity 170 kg

Basic price: **Ejoy MD XL | 10 km/h²** € 7.600,-
Max. total capacity 170 kg

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

Telescopic seating system	€
Seat width mm 420 mm 470 mm 520 mm	<input type="radio"/>
effective seat width + 50 mm extendable side panels ¹	
Seat cushion 50 mm high, max user weight 170kg Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 170kg black for seat width 420 mm / Order No. 8440920000	303,45
470 mm / Order No. 8440921000	
520 mm / Order No. 8440922000	
Seat depth mm 400 425 450 475 500 525 550	<input type="radio"/>
Seat height rear mm without seat cushion 430 450 470 490 510 530	<input type="radio"/>
Seat angle 0° – 20° manually adjustable	<input checked="" type="radio"/>
Adaptable back height 500 mm	<input checked="" type="radio"/>
Contour backrest cushion black for seat width 420 mm / Order No. 8440923000	303,45
470 mm / Order No. 8440924000	
520 mm / Order No. 8440925000	
Back support angle 90°– 120° manually adjustable, 4 steps	<input checked="" type="radio"/>
Side panel	
Adjustable in height and lateral* Side panels extendable each side +25 mm	<input checked="" type="radio"/>
Legrest	€
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable for seat width 420 mm / Order No. 8440386000	216,70
470 mm / Order No. 8440387000	
520 mm / Order No. 8440388000 (Length of footboard has to be the same size as seat width)	

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting)

Manual elevating legrest left / right unit complete left or right	123,10
for seat width 420 mm / right (Order No. 1352350001)	
470 mm / right (Order No. 1352350002)	
520 mm / right (Order No. 1352350003)	
420 mm / left (Order No. 1352351001)	
470 mm / left (Order No. 1352351002)	
520 mm / left (Order No. 1352351003)	
Lower leg length: ³ 370 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest left / right unit	179,60
left / Order No. 1352371001	
right / Order No. 1352370001	
Battery	
Battery 2 x 12 V/56 Ah	<input checked="" type="radio"/>
Battery 2 x 12 V/80 Ah upcharge for 2 pieces	288,30
Order No. 2 x 8530700000	
Charger 24 V/5 A	<input checked="" type="radio"/>
Charger 24 V/8 A ⁴ upcharge	103,-
Order No. 8130845000	
Lights	
Light system incl. indicators according to German street regulations	<input checked="" type="radio"/>

²In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.


⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

Standard equipment

Free of charge option

1 v. 2



Company:

Commission:

Date:

▼ Please mark (bold = default setting) €

Remote	
Swing-away length adjustable remote bracket	●
Control without electrical adjustments	●
VR2 Control up to 2 powered adjustments Order No. 8355805000	364,-
VR2 Attendant remote incl. bracket Order No. 8401829001	1.035,30
R-net Control up to 4 powered adjustments Order No. 8540846000	1.090,-
R-net Attendant remote incl. bracket Order No. 8401829002	1.423,10
Seat angle 0° – 20° electrically adjustable Order No. 8530961000	1.402,20
Back support angle 100° – 142° electrically adjustable ⁵ Order No. 8440960000	1.402,20
Powered elevating legrest left / right unit left / Order No. 8355361000 right / Order No. 8355360000	870,-
Remote mounting position left right	○
Attendant remote mounting position left right	○
Tyres	
Drive wheels suspensions	●
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000	●
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70
Castor 9", PU abrasion resistant, black pair Order No. 8440400000	●
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50
Anti-tipper	●
Fender front and rear	●

▼ Please mark (bold = default setting) €

Accessories	
Lap belt	●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-
Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80
Rear view mirror left / right unit incl. bracket left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	140,30
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder ⁶ Order No. 8166750000	43,60
Frame colour	
White	●

⁶In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

⁵Headrest is mandatory by powered backrest recline 100°–142°.