

Basic price: **Ejoy RD XL | 6 km/h** € 5.875,-  
 HMV-No.: 18.50.04.3002 | Max. total capacity 170 kg

Basic price: **Ejoy RD XL | 10 km/h<sup>2</sup>** € 6.725,-  
 / Max. total capacity 170 kg

**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)	€
<b>Telescopic seating system</b>	
Seat width mm <b>420 mm</b> 470 mm 520 mm	<input type="radio"/>
effective seat width   + 50 mm extendable side panels <sup>1</sup>	
Seat cushion 50 mm high, max user weight 170kg   Nylon	<input checked="" type="radio"/>
Contoured seat cushion max. user weight 170kg   black   for seat width   420 mm / Order No. 8440920000	303,45
470 mm / Order No. 8440921000	
520 mm / Order No. 8440922000	
Seat depth mm 400 425 <b>450</b> 475 500 525 550	<input type="radio"/>
Seat height rear mm   without seat cushion 430 450 <b>470</b> 490 510 530	<input type="radio"/>
Seat angle 0° – 20°   manually adjustable	<input checked="" type="radio"/>
Adaptable back height 500 mm	<input checked="" type="radio"/>
Contour backrest cushion   black   for seat width   420 mm / Order No. 8440923000	303,45
470 mm / Order No. 8440924000	
520 mm / Order No. 8440925000	
Back support angle 90°– 120°   manually adjustable, 4 steps	<input checked="" type="radio"/>
<b>Side panel</b>	
Adjustable in height and lateral* Side panels extendable each side   +25 mm	<input checked="" type="radio"/>
<b>Legrest</b>	
Swing away legrest with angle adjustable foot plate	<input checked="" type="radio"/>
Continuous footboard, swing away, angle adjustable   for seat width 420 mm / Order No. 8440386000	216,70
470 mm / Order No. 8440387000	
520 mm / Order No. 8440388000 (Length of footboard has to be the same size as seat width)	

<sup>1</sup>Effective seat width by extendable side panels.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.


▼ Please mark (bold = default setting)	€
Manual elevating legrest   left / right   unit complete left or right	123,10
for seat width   420 mm / right (Order No. 1352350001)	
470 mm / right (Order No. 1352350002)	
520 mm / right (Order No. 1352350003)	
420 mm / left (Order No. 1352351001)	
470 mm / left (Order No. 1352351002)	
520 mm / left (Order No. 1352351003)	
Lower leg length: <sup>3</sup> <b>370</b> 390 410 430 450 470 490 510 530	<input type="radio"/>
the lower leg length must not be greater than the seat height	
Amputee legrest   left / right   unit	179,60
left / Order No. 1352371001	
right / Order No. 1352370001	
<b>Battery</b>	
Battery 2 x 12 V/56 Ah	<input checked="" type="radio"/>
Battery 2 x 12 V/80 Ah   upcharge for 2 pieces	288,30
Order No. 2 x 8530700000	
Charger 24 V/5 A	<input checked="" type="radio"/>
Charger 24 V/8 A <sup>4</sup>   upcharge	103,-
Order No. 8130845000	
<b>Lights</b>	
Light system incl. indicators   according to German street regulations	<input checked="" type="radio"/>
continue next page	

<sup>2</sup>In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

<sup>3</sup>The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

<sup>4</sup>When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

Standard equipment  Free of charge option



**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting) €

<b>Remote</b>	
Swing-away length adjustable remote bracket	●
Control without electrical adjustments	●
Control up to 2 powered adjustments   Order No. 8445810000	364,-
Attendant remote   incl. bracket   Order No. 8401829001	1.035,30
Control up to 4 powered adjustments   Order No. 8445812000	1.090,-
Attendant remote   incl. bracket   Order No. 8401829002	1.423,10
Seat angle 0°– 20°   electrically adjustable   Order No. 8530961000	1.402,20
Back support angle 100°– 142°   electrically adjustable <sup>5</sup>   Order No. 8440960000	1.402,20
Powered elevating legrest left / right   unit   left / Order No. 8355361000   right / Order No. 8355360000	870,-
Remote mounting position left <b>right</b>	○
Attendant remote mounting position left <b>right</b>	○
<b>Tyres</b>	
Drive wheels suspensions	●
Drive wheel 14", PU   abrasion resistant, black   pair   Order No. 8440500000	●
Drive wheel 14", air   abrasion resistant, black   pair   Order No. 8440505000	294,70
Castor 9", PU   abrasion resistant, black   pair   Order No. 8440400000	●
Castor 9", air   abrasion resistant, black   pair   Order No. 8440405000	216,50
Anti-tipper	●
Fender   front and rear	●
<b>Accessories</b>	
Lap belt	●

▼ Please mark (bold = default setting) €

Headrest   multidimensional adjustments incl. bracket   Order No. 8440917000	215,-
Calf strap   420 mm / Order No. 8440320000   470 mm / Order No. 8440321000   520 mm / Order No. 8440322000	47,80
Rear view mirror left / right   unit incl. bracket	140,30
left / Order No. mirror 8401889000, holder 8401809000   right / Order No. mirror 8401884000, holder 8401809000	
Cane holder for 2 canes left right   unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder <sup>6</sup>   Order No. 8166750000	43,60
<b>Frame colour</b>	
White	●

<sup>6</sup>In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

<sup>5</sup>Headrest is mandatory by powered backrest recline 100°–142°.