

ORDER FORM

EJOY RD XXL

Electrically powered wheelchair

Company:

Commission:

Date:

Basic price: **Ejoy RD XXL** | 6 km/h € 14.650,-
 With VR2 Control, incl. 2 electrical adjustments
 | max. total capacity 200 kg
 HMV-No.: 18.50.04.5007

Basic price: **Ejoy RD XXL** | 6 km/h € 16.790,-
 With R-Net Control, incl. 4 electrical adjustments
 | max. total capacity 250 kg
 HMV-Nr.: 18.50.04.5007

▼ Please mark (bold = default setting) €

Telescopic seating system	
Seat width mm 580 mm 630 mm 680 mm ○	
effective seat width +120 mm extendable side panels ¹	
Seat cushion 50 mm high, max user weight 250kg Nylon ●	
Contoured seat cushion max. user weight 250 kg black for seat width	
580 mm / Order No. 8540920000 303,45	
630 mm / Order No. 8540921000	
680 mm / Order No. 8540922000	
Seat depth mm ○	
400 425 450 475 500 525 550	
Seat height rear mm without seat cushion ○	
430 450 470 490 510 530	
Seat angle 0° – 20° manually adjustable ●	
Contour backrest cushion black	
for seat width	
580 mm / Order No. 8540923000 303,45	
630 mm / Order No. 8540924000	
680 mm / Order No. 8540925000	
Back support angle 90°– 120° ●	
manually adjustable, 4 steps	
Side panel	
Adjustable in height and lateral ●	
Side panels extendable each side +60 mm	
Legrest	
Swing away legrest with angle adjustable foot plate ●	
Continuous footboard, swing away, angle adjustable ●	
for seat width	
580 mm / Order No. 8540386000	
630 mm / Order No. 8540387000	
680 mm / Order No. 8540388000	
(Length of footboard has to be the same size as seat width)	
Manual elevating legrest left / right 123,10	
right (Order No. 1352350003)	
left (Order No. 1352351003)	

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting) €

Lower leg length: ² 370 390 410 430 450 470 490 510 530 ○	
the lower leg length must not be greater than the seat height	
Amputee legrest left / right unit 179,60	
left / Order No. 1352371001	
right / Order No. 1352370001	
Battery	
Battery 2 x 12 V/80 Ah ●	
Charger 24 V/8 A ●	
Lights	
Light system incl. indicators according to German street regulations ●	
Remote	
Swing-away length adjustable remote bracket ●	
Electrical adjustments options	
Attendant remote incl. bracket 1.035,30	
Bestellnr. 8401829001	
R-net Control with 4 electrical adjustment options 1.090,-	
included as standard for 250 kg model	
Order No. 8445812000 1.423,10	
Attendant remote incl. bracket	
Order No. 8401829002	
Seat support angle 0° – 20° electrically adjustable 1.402,20	
Order No. 8530961000	
Back support angle 100°– 142° electrically adjustable ³ 1.402,20	
Order No. 8440960000	
Powered elevating legrest left / right unit 870,-	
left / Order No. 8355361000	
right / Order No. 8355360000	
Remote mounting position left right ○	
Attendant remote mounting position left right ○	
continue next page	


²The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

³Headrest is mandatory by powered backrest recline 100°–142°.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment ○ Free of charge option

1 v. 2



ORDER FORM

EJOY RD XXL Electrically powered wheelchair

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

€

Tyres	
Drive wheels suspensions	<input checked="" type="radio"/>
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000	<input checked="" type="radio"/>
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70
Castor 9", PU abrasion resistant, black pair Order No. 8440400000	<input checked="" type="radio"/>
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50
Fender front and rear	<input checked="" type="radio"/>
Anti-tipper	<input checked="" type="radio"/>
Accessories	
Lap belt	<input checked="" type="radio"/>
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-
Calf strap 580mm / Order No. 8540320000 630 mm / Order No. 8540321000 680 mm / Order No. 8540322000	47,80
Rear view mirror left / right unit incl. bracket	140,30
left / Order No. mirror. 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder Order No. 8166750000	43,60
Frame colour	
White	<input checked="" type="radio"/>

Standard equipment

Free of charge option

2 v. 2

Version: 2025.03.03 Prices plus VAT at the legally applicable rate, subject to technical changes.