



Company:

Basic price: **Protego without drum brake** € 2.595.-

Basic price: **Protego with drum brake** € 2.865.-

Customer:

Max. capacity 130 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 390 - 415 440 - 465 490 - 515	<input type="radio"/>
Seat height front 430 mm without cushion	<input checked="" type="radio"/>
Seat depth adjustable from 420 - 480 mm	<input checked="" type="radio"/>
Seat angle adjustable 0 - 15°	<input checked="" type="radio"/>
Medial support mit Halter Order No. 1531915000	106.-
Seat upholstery Dartex, night grey	<input checked="" type="radio"/>
Back	
Back support ergonomic shape with lateral support	<input checked="" type="radio"/>
Back support height 580 – 640 mm	<input checked="" type="radio"/>
Back support angle 0 – 44°	<input checked="" type="radio"/>
Trunk support swing away Stück Left 1531911000	108.50
Right 1531910000	
Tool holder lateral trunk support Stück Order No. 1531914000	15.90
Back support adjustable by strap system	<input checked="" type="radio"/>
Back support upholstery Dartex, night grey	<input checked="" type="radio"/>
Handgrip elevating, removable	<input checked="" type="radio"/>
Sidepanels	
Sidepanels adjustable in height and depth, removable with arm support upholsteries Dartex upholstery	<input checked="" type="radio"/>
Hemiplegic arm support left right Stück Order No. 1520803000	341.70
Leg support	
Leg support elevating and angle adjustable with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" x 45 mm PU	<input checked="" type="radio"/>
Spoke protection, transparent 24" Pair Order No. 9800054000	165.90
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension Order No. 0009677000 Stück	26.30
Accessories	
Togo power push aid incl. mounting set Max. capacity 125 kg Order no. 8010000000	2.770.-
Mounting aid for adaptation Order No. 8010199000	216.80
Tilting aid	<input checked="" type="radio"/>
Safety belt Order No. 9800012000	30.60
Therapy tray	
Seat width (S) 370 – 430 Order No. 9800006000	255.-
Seat width (L) 430 – 490 Order No. 9800007000	255.-
Seat width (XL) 490 – 560 Order No. 9800008000	285.-
Therapy tray Hemi Order No. 9800002000	315.-
Passive light	<input checked="" type="radio"/>
Frame colour	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge