



Company:

Basic price: **Protego SL without drum brake** € 3.740,-

Basic price: **Protego SL with drum brake** € 4.060,-

Customer:

Max. capacity 150 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 440 - 465 490 - 515	<input type="radio"/>
Seat height 520 - 600 mm with seat cushion	<input checked="" type="radio"/>
Seat depth 480 - 520 mm	<input checked="" type="radio"/>
Contour seat cushion	<input checked="" type="radio"/>
Seat angle adjustable 0 - 15°	<input checked="" type="radio"/>
Medial support Order no. 1531915000	106.-
Seat upholstery Dartex, night grey	<input checked="" type="radio"/>
Back	
Back support ergonomic shape with lateral support	<input checked="" type="radio"/>
Back support height 540 - 650 mm with seat cushion	<input checked="" type="radio"/>
Back support angle 0 - 44°	<input checked="" type="radio"/>
Trunk support swing away left Order no. 1531911000 right Order no. 1531910000	108.40
Clamp holder for side support piece Order no. 1531914000	15.90
Head support incl. fixation	<input checked="" type="radio"/>
Back support upholstery Dartex, night grey	<input checked="" type="radio"/>
Handgrip elevating, removable	<input checked="" type="radio"/>
Sidepanels	
Sidepanels adjustable in height and depth, removable, with Dartex arm support	<input checked="" type="radio"/>
Hemiplegic arm support left right Order no. 1531810000	341.70

▼ Please mark!

Leg support	€
Leg support elevating and angle adjustable with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" x 45 mm PU	<input checked="" type="radio"/>
Spoke protection, transparent 24" Order no. 9800054000	165.80
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension Order no. 1270877000	26.30
Accessories	
Togo power push aid incl. mounting set Max. capacity 125 kg Order no. 8010000000	2.770.-
Mounting aid for adaptation Order no. 8010199000	216.80
Tilting aid	<input checked="" type="radio"/>
Safety belt Order no. 9800012000	30.60
Therapy tray Seat width (L) 430 - 490 Order No. 9800007000 Seat width (XL) 490 - 560 Order No. 9800008000 Therapy tray Hemi Order No. 9800002000	255.- 285.- 315.-
Passive light	<input checked="" type="radio"/>
Frame colour	
Anthracite	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge