

# ORDER FORM

## PYRO LIGHT OPTIMA

Light weight wheelchair

**Company:**

Basic price: **Pyro Light Optima without drum brake** € 1.650,-

Basic price: **Pyro Light Optima with drum brake** € 1.850,-

Max. capacity 125 kg

**Customer:**

**Date:**

▼ Please mark!

<b>Seat</b>	€
Seat width in mm 370    400    430    460    490	<input type="radio"/>
520	47.30
Seat height front in mm 500	<input checked="" type="radio"/>
400   included front wheel 6". rear wheel 22"	52.10
425    450    475	52.10
Seat depth 425 mm (Standard) adjustable from 400 mm up to 500 mm in 25 mm steps	<input checked="" type="radio"/>
Seat cushion height 5 cm	
Seat width 370 mm   Order no. 9800080000	27.50
Seat width 400 mm   Order no. 9800081000	27.50
Seat width 430 mm   Order no. 9800082000	27.50
Seat width 460 mm   Order no. 9800083000	27.50
Seat width 490 mm   Order no. 9800084000	29.70
Seat width 520 mm   Order no. 9800085000	38.90
Seat upholstery. Nylon	<input checked="" type="radio"/>
<b>Back</b>	
Back support height 420 mm	<input checked="" type="radio"/>
Back support height adjustable 400 – 460 mm	<input checked="" type="radio"/>
Back support with lumbar curve	<input checked="" type="radio"/>
Back support adjustable by strap-system	<input checked="" type="radio"/>
Head support*   incl. fixation für SB 380 – 430 mm   Order no. 9802900000	215.-
für SB 450 – 520 mm   Order no. 9802901000	
Handgrip elevating removable	<input checked="" type="radio"/>
Back support upholstery Nylon	<input checked="" type="radio"/>
<b>Sidepanels</b>	€
Sidepanels adjustable in depth and height	<input checked="" type="radio"/>


▼ Please mark!

<b>Leg support</b>	
Swing-away leg support with individual footplates plastic. removable. adjustable in height	<input checked="" type="radio"/>
Leg support elongation 100 mm   piece   Order no. 1330333000	34.80
Lower leg support elevating   left / right left SB 370 – 400   Order no. 1352351001 left SB 430 – 460   Order no. 1352351002 left SB 490 – 520   Order no. 1352351003	123.10
right SB 370-400   Order no. 1352350001 right SB 430-460   Order no. 1352350002 right SB 490-520   Order no. 1352350003	
Amputee support   left / right left Order no.   1352371001 right Order no.   1352370001	179.60
Footplates adjustable in angle and in width	<input checked="" type="radio"/>
<b>Tyres</b>	
Rear wheels 24" x 1 3/8" PU   incl. quick-release axle Front wheel 8" x 1.25" PU	<input checked="" type="radio"/>
Front wheel 200 x 50 mm   PU solid rubber   Order no. 9500132000	31.60
Front wheel 6" complete   PU solid rubber   Order no. 9500114000	42.70
Hand rim rubber-coated 24"   Pair   Order no. 9800066000	180.80
Spoke protection. transparent 24"   Order no. 9800054000	165.80
One arm drive   only without drum brake left / right left für SB 370 – 400   Order no. 5000017646 right für SB 430 – 520   Order no. 5000017645	649.-
continue next page	

\* Height-adjustable headrest and push tubes cannot be fitted as a combination

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment     Optional, no extra charge



# ORDER FORM

## PYRO LIGHT OPTIMA

Light weight wheelchair

**Company:**

**Customer:**

**Date:**

▼ Please mark!

<b>Wheel locks</b>	
Knee lever wheel lock   push to lock	<input checked="" type="radio"/>
Wheel lock lever extension   piece	26.30
Order no.   0009677000	
Knee lever wheel lock for one-arm drive left / right	
left Order no.   9000611000	101.40
right Order no.   9000610000	
<b>Accessories</b>	
Togo power push aid incl. mounting set	2.770.-
Max. capacity 125 kg	
Order no. 8010000000	
Mounting aid for adaptation   Order no. 8010199000	216.80
Tilting aid	<input checked="" type="radio"/>
Anti-tip device   Order no. 1270845000   pair	100.10
Lap belt   Order no. 9800012000	30.60
Lap belt XL   Order no. 9800017000	35.30
Cane holder   Order no. 9800030002	27.10
Therapy tray	
Seat width (S) 370 – 430   Order No. 9800006000	255.-
Seat width (L) 430 – 490   Order No. 9800007000	255.-
Seat width (XL) 490 – 560   Order No. 9800008000	285.-
Therapy tray Hemi   Order No. 9800002000	315.-
Passive light	<input checked="" type="radio"/>
<b>Frame colour</b>	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge