



ORDER FORM

PYRO LIGHT XL Light weight wheelchair

Company:

Basic price: **Pyro Light XL without hub brake** € 1.950,-

Basic price: **Pyro Light XL with hub brake** € 2.100,-

Customer:

Max. capacity 170 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 510 560	<input type="radio"/>
Seat height front in mm 500	<input checked="" type="radio"/>
420 440 450 470 520 540	52.10
Seat depth adjustable 440 - 460 mm	<input checked="" type="radio"/>
Seat cushion Seat width 510 mm Order no. 9800086000	36.60
Seat cushion Seat width 560 mm Order no. 9800087000	46.-
Seat upholstery, Vinyl	<input checked="" type="radio"/>
Back support	
Back support height 420 mm	<input checked="" type="radio"/>
Adaption kit for height adjustment 400 - 440 mm Order no. 1270875000	126.-
Back support adjustable by strap-system 51cm: Order no. 1320806000	175.70
56cm: Order no. 1320807000	
Handgrip elevating, removable* Pair Without hub brake Order no. 1270806000	174.40
With hub brake Order no. 1270807000	
Head support* incl. fixation SB 510 mm Order no. 9802901000	215.-
SB 560: not available	
Back support upholstery, Vinyl	<input checked="" type="radio"/>
Sidepanels	
Desk, short version, swing away	<input type="radio"/>
Long version, swing away	
Sidepanel elevating left Order no. 1352805000	123.80
right Order no. 1352806000	


▼ Please mark!

Leg support	€
Swing-away leg support with individual footplates plastic, removable, adjustable in height	<input checked="" type="radio"/>
Leg support elongation 100 mm Pair Order no. 1330333000	69.60
Lower leg support elevating left Order no. 1352351003	123.10
right Order no. 1352350003	
Amputee support left Order no. 1352371001	179.60
right Order no. 1352370001	
Footplate with heel strap	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 200 x 50 mm PU	<input checked="" type="radio"/>
Hand rim rubber-coated 24" Pair Order no. 9800066000	180.80
Spoke protection, transparent 24" Pair Order no. 9800054000	165.80
Front wheel boom piece Order no. 1260842000	40.-
One arm drive only without drum brake SB 510 left right Order no. 5000017645	371.80
SB 560: not available	
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension Piece Order no. 0009677000	26.30
Knee lever wheel lock for one-arm drive left Order no. 9000611000	101.40
right Order no. 9000610000	
Accessories next page	

* Head support and elevating handgrips cannot be mounted together

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge



ORDER FORM

PYRO LIGHT XL Light weight wheelchair

Company:

Customer:

Date:

▼ Please mark!

Accessories	
Togo power push aid** incl. mounting set Max. capacity 125 kg Order no. 8010000000	2.770.-
Telescopic tubes for seat width 560 mm** Order no. 1x8010115000; 1x8010114000; 2x8010116000	224.20
Mounting aid for Togo adaptation Order no. 8010199000	216.80
Anti-tip device Pair Order no. 1270845000	100.10
Tilting aid	<input checked="" type="radio"/>
Lap belt Order no. 9800015000	35.30
Crutch holder Order no. 9800030002	27.10
Therapy tray	
Seat width (S) 370 – 430 Order No. 9800006000	255.-
Seat width (L) 430 – 490 Order No. 9800007000	255.-
Seat width (XL) 490 – 560 Order No. 9800008000	285.-
Therapy tray Hemi Order No. 9800002000	315.-
Passive light	<input checked="" type="radio"/>
Frame Colour	
Silver	<input checked="" type="radio"/>

** For sw 560 mm additional telescopic tubes are required

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge

Version 2025.02.24 Prices plus VAT at the legally applicable rate, subject to technical changes.