

ORDER FORM

PYRO START PLUS

Light weight wheelchair

Company:

Customer:

Date:

Basic price: **Pyro Start Plus without drum brake** € 1.650,-

Basic price: **Pyro Start Plus with drum brake** € 1.730,-

Basic price: **Pyro Start Plus with drum brake and HA** € 1.850,-

| with height adjustment for höhenverstellbaren side panels and adjustable back

Max. capacity 125 kg

▼ Please mark!

€

Seat		
Seat width in mm	370 400 430 460 490	<input type="radio"/>
	520	47,30
Seat height front in mm	500	<input checked="" type="radio"/>
	450 475	52,10
Seat depth adjustable from 400 - 440 mm		<input checked="" type="radio"/>
Seat cushion height 5 cm Seat width 370 - 520 mm		
Seat width 370 mm Order No. 9800080000		27,50
Seat width 400 mm Order No. 9800081000		27,50
Seat width 430 mm Order No. 9800082000		27,50
Seat width 460 mm Order No. 9800083000		27,50
Seat width 490 mm Order No. 9800084000		29,70
Seat width 520 mm Order No. 9800085000		38,90
Seat upholstery, Nylon		<input checked="" type="radio"/>
Back support		
Back support with lumbar curve		<input checked="" type="radio"/>
Back support adjustable 420 - 440 mm		<input checked="" type="radio"/>
Back support adjustable by strap-system		175,70
Seat width 370 mm Order No. 1352935000		
Seat width 400 mm Order No. 1352936000		
Seat width 430 mm Order No. 1352937000		
Seat width 460 mm Order No. 1352938000		
Seat width 490 mm Order No. 1352939000		
Seat width 520 mm Order No. 1352940000		
Handgrip elevating*, removable Pair		174,40
Head support* incl. fixation		
Seat width 370 - 440 mm Order No. 9802900000		215,10
Seat width 460 - 560 mm Order No. 9802901000		
Back support upholstery, Nylon		<input checked="" type="radio"/>

* Head support and handgrips elevating cannot be mounted together

▼ Please mark!


€

Sidepanels		
Sidepanel removable with arm support, swing-away, elevating		<input checked="" type="radio"/>
Sidepanel elevating left / right		
left Order No. 1352805000		123,90
right Order No. 1352806000		
Sidepanel elevating with foldable arm pads left / right		
Left Order No. 1352811000		141,30
Right Order No. 1352812000		
Leg support		
Swing-away leg support with individual footplates plastic, removable, adjustable in height		<input checked="" type="radio"/>
Lower leg support elevating left / right piece		
left SB 370 - 400 Order No. 1352351001		
left SB 430 - 460 Order No. 1352351002		
left SB 490 - 520 Order No. 1352351003		123,10
right SB 370-400 Order No. 1352350001		
right SB 430-460 Order No. 1352350002		
right SB 490-520 Order No. 1352350003		
Leg support elongation 100 mm Pair		69,60
Amputee support left / right		
left Order No. 1352371001		179,60
right Order No. 1352370001		
Footplates with heel straps		<input checked="" type="radio"/>
Wheel locks		
Knee lever wheel lock push to lock		<input checked="" type="radio"/>
Wheel lock lever extension		26,30
Order No. 0009677000		
Knee lever wheel lock for one-arm drive left / right		
left Order No. 9000611000		101,40
right Order No. 9000610000		

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge



ORDER FORM

PYRO START PLUS

Light weight wheelchair

Company:

Customer:

Date:

▼ Please mark!

€

Tyres		
Rear wheels 24" x 1 3/8" PU incl. quick-release axle		<input checked="" type="radio"/>
Front wheels 8" x 1 1/4", PU		<input checked="" type="radio"/>
Front wheels 200 x 50 mm. PU Pair	69,60	
Order No. 9500132000		
PU tyres		<input checked="" type="radio"/>
Wheel base extension		<input checked="" type="radio"/>
Adapter block		<input checked="" type="radio"/>
Spoke protection, transparent 24" pair	165,90	
Order No. 9800054000		
One arm drive only without drum brake left / right		
left for seat width 370 – 400 Order No. 5000017646	649.-	
right for seat width 430 – 520 Order No. 5000017645		
Hand rim rubber-coated 24" Pair	180,80	
Order No. 9800066000		
Accessories		
Togo power push aid incl. mounting set		
Max. capacity 125 kg	2.770,-	
Order no. 8010000000		
Mounting aid for adaptation	216,80	
Order No. 8010199000		
Anti-tip device Pair	100,20	
Order No. 1270845000		
Tilting aid		<input checked="" type="radio"/>
Lap belt	30,70	
Order No. 9800012000		
Crutch holder	27,10	
Order No. 9800030002		
Therapy tray		
Seat width (S) 370 – 430 Order No. 9800006000	255.-	
Seat width (L) 430 – 490 Order No. 9800007000	255.-	
Seat width (XL) 490 – 560 Order No. 9800008000	285.-	
Therapy tray Hemi Order No. 9800002000	315.-	
Passive light		<input checked="" type="radio"/>
Frame colour		
Silver	Blue	<input type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge