



ORDER FORM

PYRO START PLUS SL

Light weight wheelchair

Company:

Basic price: **Pyro Start Plus SL without drum brake** € 1.750.-

Basic price: **Pyro Start Plus SL with drum brake** € 1.950.-

Customer:

Max. capacity 150 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 430 460 490 520	<input type="radio"/>
Seat height front 525 mm adjustable from 525 - 600 mm in 25 mm steps	<input checked="" type="radio"/>
Seat depth 460 mm adjustable from 460 - 520 mm in 20mm steps.	<input checked="" type="radio"/>
Seat cushion Seat width 460 - 480 mm	33,-
Seat cushion Seat width 500 - 520 mm	37,40
Seat upholstery, Nylon	<input checked="" type="radio"/>
Back support	
Back support with lumbar curve	<input checked="" type="radio"/>
Back support height 460 mm adjustable 460 - 500 mm in 20 mm steps	<input checked="" type="radio"/>
Back support adjustable by strap-system	<input checked="" type="radio"/>
Handgrip elevating*, removable pair	167,70
Head support* incl. fixation	206,80
Back support upholstery, Nylon	<input checked="" type="radio"/>
Sidepanels	
Sidepanel removable with arm support, swing-away, elevating	<input checked="" type="radio"/>
Sidepanel elevating left right piece	119,10
Sidepanel elevating with foldable arm pads left right piece	135,90
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	<input checked="" type="radio"/>
Leg support elongation 100 mm pair	67,-
Lower leg support elevating left right piece	118,40
Amputee support left right piece	172,70
Footplates with heel straps	<input checked="" type="radio"/>

* Head support and handgrips elevating cannot be mounted together

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 8" x 1 1/4", PU	<input checked="" type="radio"/>
Front wheels 200 x 50 mm PU pair	67,-
PU-tyres	<input checked="" type="radio"/>
Wheel base extension	<input checked="" type="radio"/>
Adapter block	<input checked="" type="radio"/>
Spoke protection, transparent 24" pair	159,50
One arm drive only without drum brake left right	357,50
Hand rim rubber-coated 24" pair	173,80
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension piece	25,30
Knee lever wheel lock for one-arm drive left right	97,50
Accessories	
Togo power push aid incl. mounting set Max. capacity 125 kg Order no. 8010000000	2.770.-
Mounting aid for Togo adaptation	216,80
Retrofit kit Transit	231.-
Anti-tip device pair	96,30
Tilting aid	<input checked="" type="radio"/>
Safety belt	29,50
Crutch holder	26,10
Therapy tray	
Seat width (S) 370 – 430 Order No. 9800006000	255.-
Seat width (L) 430 – 490 Order No. 9800007000	255.-
Seat width (XL) 490 – 560 Order No. 9800008000	285.-
Therapy tray Hemi Order No. 9800002000	315.-
Passive light	<input checked="" type="radio"/>
Frame colour	
Anthracite	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge