



**Company:**

Basic price: **S - VR without drum brake**

€ 1.850,-

Basic price: **S - VR with drum brake**

€ 2.050,-

**Customer:**

Max. capacity 120 kg

**Date:**

▼ Please mark!

<b>Seat</b>	€
Seat width in mm 400    440    480	<input type="radio"/>
Seat height 480 mm	<input checked="" type="radio"/>
Seat depth 420 mm	<input checked="" type="radio"/>
Seat angle adjustable 4°	<input checked="" type="radio"/>
Seat cushion   50 mm height	24,-
<b>Back</b>	
Back support height 430 / 850 mm	<input checked="" type="radio"/>
Back support angle -4° - +30° / +82° adjusted mechanically	<input checked="" type="radio"/>
Neck holder adjustable in height, upholstered, removable	<input checked="" type="radio"/>
Back support upholstery, vinyl black	<input checked="" type="radio"/>
<b>Sidepanels</b>	
Sidepanels removable	<input checked="" type="radio"/>
Arm support elongated	<input checked="" type="radio"/>

▼ Please mark!

<b>Leg support</b>	€
Swing-away leg support with individual footplates   plastic, removable, adjustable in height	<input checked="" type="radio"/>
<b>Tyres</b>	
Rear wheels 24" x 1 3/8" PU	<input checked="" type="radio"/>
Front wheels 200 x 50 mm PU	<input checked="" type="radio"/>
<b>Wheel locks</b>	
Knee lever wheel lock   push to lock	<input checked="" type="radio"/>
<b>Accessories</b>	
Anti-tip device   Pair	<input checked="" type="radio"/>
Lap belt	26,-
Therapy tray	
Seat width (S) 370 – 430   Order No. 9800006000	255.-
Seat width (L) 430 – 490   Order No. 9800007000	255.-
Therapy tray Hemi   Order No. 9800002000	315.-
Fold blocking   anti-theft device	<input checked="" type="radio"/>
<b>Frame colour</b>	
Anthracite	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge