



ORDER FORM

TRITON Multifunctional wheelchair

Company:

Basic price: **Triton without drum brake** € 2.499,-

Basic price: **Triton with drum brake** € 2.758,-

Customer:

Max. capacity 125 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm	
Variation I: 360 - 430	<input type="radio"/>
Variation II: 430 - 500	
Seat height 415 mm with seat cushion 520 mm	<input checked="" type="radio"/>
Seat depth 420 - 480 mm	<input checked="" type="radio"/>
Seat angle adjustable 0° - 15°	<input checked="" type="radio"/>
Medial support with holder	102,-
Seat upholstery Dartex night grey	<input checked="" type="radio"/>
Back	
Comfort back support, shaped ergonomically	<input checked="" type="radio"/>
Back support height adjustable	<input checked="" type="radio"/>
Back support reclining from -3° up to +43°	<input checked="" type="radio"/>
Lateral trunk support swing away	<input checked="" type="radio"/>
Handgrip elevating, removable 85 - 102 cm	<input checked="" type="radio"/>
Head support adjustable in 3 directions	<input checked="" type="radio"/>
Back support upholstery Dartex night grey	<input checked="" type="radio"/>
Head support extension	72,60
Sidepanels	
Sidepanels adjustable in height and depth, removable with upholstered arm support	<input checked="" type="radio"/>
Hemiplegic arm support	left right 328,60

▼ Please mark!

Leg support	€
Leg support elevating horizontally, with length adjustment	<input checked="" type="radio"/>
Anterior lower leg support	<input checked="" type="radio"/>
Posterior lower leg support	<input checked="" type="radio"/>
Individual footplates angle adjustable	<input checked="" type="radio"/>
Tyres	
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheels 7" PU	<input checked="" type="radio"/>
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Accessories	
Anti-tip device pair	96,30
Lap belt	29,50
Therapy tray	
Seat width (S) 370 - 430 Order No. 9800006000	255,-
Seat width (L) 430 - 490 Order No. 9800007000	255,-
Seat width (XL) 490 - 560 Order No. 9800008000	285,-
Therapy tray Hemi Order No. 9800002000	315,-
Passive light	<input checked="" type="radio"/>
Frame colour	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment Optional, no extra charge