

BX 11 (with continuous footplate)



BX 11 S (with detachable legrests)

Cost proposal

Order


Invoice address/Customer number	
Cust. no.	
Komission	
Name	
Street	
Postcode/Area	
E-mail	
Tel.	

Shipping address	
Cust. no.	
Komission	
Name	
Street	
Postcode/Area	
E-mail	
Tel.	

The order takes place subject to the selected combination being checked for technical feasibility.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Version 2025.02.24 Prices plus VAT at the legally applicable rate, subject to technical changes.



General information										
User weight 135 kg Material aluminium Quick-release axle Depth-adjustable seat sling Adjustable back upholstery, breathable, black Clothes guard, aluminium Knee lever brake (B+B active pressure brake) Aluminium handrims, black, powder-coated With the drum brake option: handrims and rims in silver, PU tyres										
Front frame									Base price*	
BX 11 closed front frame with continuous footplate HMV-No.: 18.50.03.0200									2.100,-	
BX 11 S two-piece front frame for removable legrests HMV-No.: 18.50.03.0200									2.150,-	
Legrests										
One piece, with aluminium footplate Standard with BX 11 angle-adjustable, flip-up ●										
Two-piece, removable with plastic footplates Standard with BX 11 S pair in frame colour footplates angle-adjustable, can be flipped up ●										
Legrest horizontal adjustable in black only with BX 11 S										
left									118,40	
right									118,40	
Amputee support in black only with BX 11 S										
left									172,70	
right									172,70	
Frame colours										
Black	White	Blue	Red							<input type="radio"/>
User-related dimensions										
Front seat height (FSH) in mm <i>see matrix below</i>										
490	500	510	520	530	540	550	560	<input type="radio"/>		
Rear seat height (RSH) in mm <i>see matrix below</i> with 24" drive wheels, with 22" approx. 25 mm less										
420	440	460	480	500	520	<input type="radio"/>				
Attention: from a difference of 60 mm between FSH and RSH, the angle-adjustable back must be installed.										
Backrest height (RH) in mm										
275	300	325	350	375	400	<input type="radio"/>				
Lower leg length BX 11 (UL) in mm <i>see matrix below</i> minimum gap from footplate to floor is 40 mm										
430	445	460	475	490	505	520	535	<input type="radio"/>		
Lower leg length BX 11 S (UL) in mm minimum gap from footplate to floor is 40 mm										
430	450	470	490	510	<input type="radio"/>					
Seat width (SW) & Seat depth (SD) in mm <i>combination only inside the table possible</i>										
SW		SD		SW		SD				
360	380			420	440	480	500	Please select only from <u>one</u> of the three tables ! <input type="radio"/>		
380	400			440	460	500	520			
400	420**			460	480**	520	540**			

* Price plus statutory VAT

**Seat depth not possible with BX 11 S variants

		UL							
		430	445	460	475	490	505	520	535
FSH ¹	490	x	x						
	500	x	x	x	x				
	510	x	x	x	x				
	520	x	x	x	x	x	x		
	530	x	x	x	x	x	x	x	
	540	x	x	x	x	x	x	x	
	550	x	x	x	x	x	x	x	x
	560	x	x	x	x	x	x	x	x

¹Front Seat Height


		Castor Size			
		4"	5"	6"	7"
FSH ¹	490	x	x		
	500	x	x	x	
	510	x	x	x	x
	520	x	x	x	x
	530	x	x	x	x
	540	x	x	x	x
	550		x	x	x
	560			x	x

		Castor Size			
		4"	5"	6"	7"
RSH ²	420	x	x	x	x
	440	x	x	x	x
	460	x	x	x	x
	480	x	x	x	x
	500	x	x	x	x
	520	x	x	x	x

With 7" front wheels, only centres of gravity 30 - 54 mm possible.

²Rear Seat Height

● Basic configuration ○ optional without additional cost



Frame accessories		
Tipping aids	left	78,30
	right	78,30
Anti-tip wheels pair		225,50
Seat cushion thickness 50mm		
Desired dimension in mm:		
Backrest system and sliding tubes		
Straight backrest, adjustable sliding tubes		●
Angle-adjustable backrest, height-adjustable sliding tubes, removable		255,-
Sideguards		
Aluminium mudguard Splashguard lip pair		●
Detachable, height-adjustable with armrest pair		204,-



One piece legrests, with footplates made of aluminium



Two-piece legrests, removable with footplates, plastic



Legrest horizontal adjustable



Amputee support



Angle-adjustable backrest



Height-adjustable detachable



Aluminium clothes guard (Splashguard lip)



Sideguard, detachable, height-adjustable




Anti-tip wheels




Sopportes pisador

● Basic configuration ○ optional without additional cost



Front wheels <i>Please note the interdependency between the lower leg length and front seat height (matrix page 2).</i>		
With 7" front wheels, only centres of gravity 30 - 54 mm possible.		
4" plastic rim, soft rubber, black		<input type="radio"/>
5" plastic rim, soft rubber, black		<input type="radio"/>
6" plastic rim, soft rubber, black		<input type="radio"/>
6" plastic rim, PU grey		<input type="radio"/>
7" plastic rim, soft rubber, black		<input type="radio"/>
Drive wheels <i>Please note the interdependency on the rear seat height (matrix page 2)</i>		
24" PU on aluminium hollow rim pair		<input type="radio"/>
22" pneumatic tyres with aluminium hollow rim pair		<input type="radio"/>
22" Schwalbe Marathon plus on aluminium hollow rim pair		173,80
24" pneumatic tyres with aluminium hollow rim pair		<input type="radio"/>
24" Schwalbe Marathon plus on aluminium hollow rim pair		173,80
Rear wheel camber		
0° Standard 1° 2°		<input type="radio"/>
Centre of gravity in mm <i>measured from the front edge of the backrest tube, see sketch below (66 and 78 not possible on SW 360 380 400)</i>		
30 42 54 Standard 66 78 30 = passive, 78 = active		<input type="radio"/>
Fitting handrims		
20 mm gap Standard		<input type="radio"/>
30 mm gap <i>not possible with drum brake and PU rear wheels</i>		
Drive wheel accessories		
Transparent 22" spokeguard pair		159,50
Transparent 24" spokeguard pair		159,50
Brake		
Drum brake option <i>including drum brake wheels 24" PU</i>		316,80
Brake lever extension left		25,30
Brake lever extension right		25,30
Seat accessories		
Lap belt		29,50

Notes



Centre of gravity, side view

● Basic configuration ○ optional without additional cost