

ORDER FORM

PYRO LIGHT OPTIMA XL

Light weight wheelchair

Company:

Basic price: **Pyro Light Optima XL without drum brake** € 2.150.-

Customer:

Basic price: **Pyro Light Optima XL with drum brake** € 2.350.-

Date:

Max. capacity 170 kg

▼ Please mark!

Seat	€
Seat width in mm 520 550* 580	○
Seat height front in mm 500	●
400 included front wheel 6", rear wheel 22"	52.10
425 450 475	52.10
Seat depth 450 mm (Standard) adjustable from 450 mm up to 500 mm in 25 mm steps	●
Seat cushion height 7 cm Sitzbreite 520 mm Order no. 9800089000	34,-
Sitzbreite 550 mm Order no. 9800097000	46,-
Sitzbreite 580 mm Order no. 9800098000	34,-
Seat cushion height 5 cm Sitzbreite 520 mm Order no. 9800085000	38,90
Sitzbreite 550 mm Order no. 9800087000	46,-
Sitzbreite 580 mm Order no. 9800088000	46,-
Seat cushion Seat width 580 mm	
Seat upholstery, Nylon	●
Back	
Back support height 420 mm	●
Back support height adjustable 400 - 460 mm	●
Back support with lumbar curve	●
Back support adjustable by strap-system	●
Head support incl. fixation Not available for sw 580 SB 520 mm SB 550 mm Order no. 9802901000	215,-
Handgrip elevating, removable	●
Back support upholstery, Nylon	●
Sidepanels	
Sidepanels adjustable in depth and height	●

▼ Please mark!


Leg support	€
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elongation 100 mm piece Order no. 1330333000	34.80
Lower leg support elevating left Order no. 1352351003	123.10
right Order no. 1352350003	
Amputee support left right	179.60
Footplates adjustable in angle and in width	●
Tyres	
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	●
Front wheel 8" 200 x 50 mm, PU	●
Hand rim rubber-coated 24" Pair Order no. 9800066000	180.80
Spoke protection, transparent 24" Order no. 9800054000	165.80
One arm drive only without drum brake left für SB 520 Order no. 5000017646	649.-
right für SB 550 Order no. 5000017645 for SB 580: not available	
Wheel locks	
Knee lever wheel lock push to lock	●
Wheel lock lever extension Order no. 0009677000	26.30
Knee lever wheel lock for one-arm drive left Order no. 9000611000	101.40
right Order no. 9000610000	
Accessories next page	

* For seat width 550 mm these additional telescopic tubes are required.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge



ORDER FORM

PYRO LIGHT OPTIMA XL

Light weight wheelchair

Company:

Customer:

Date:

▼ Please mark!

Accessories		
Togo power push aid incl. mounting set		
only up to SW 550 mm		2.770,-
Max. capacity 125 kg		
Order no. 8010000000		216.80
Mounting aid for Togo adaptation		
Order no. 8010199000		
Long telescopic tubes Togo assembly seat width 550 mm*		203.-
Order no. 1x8010115000; 1x8010114000; 2x8010116000		
Tilting aid	●	
Anti-tip device Pair Order no. 1270845000		100.10
Lap belt XL Order no. 9800017000		35.30
Cane holder Order no. 9800030002		27.10
Therapy tray		
Seat width (S) 370 – 430 Order No. 9800006000		255.-
Seat width (L) 430 – 490 Order No. 9800007000		255.-
Seat width (XL) 490 – 560 Order No. 9800008000		285.-
SB 580: nicht verfügbar		
Therapy tray Hemi Order No. 9800002000		315.-

▼ Please mark!

Passive light	●
Wheelchair accessory set consisting of:	
Universal cane holder for 2 walking sticks, 1 pc	
"Press" standard safety belt, 1 pc	
Seat cushion, black, 1 pc	
1 pair of anti-tip wheels	
Accessory-Set S (SW 37-40 cm) Order. no. 9802928000	229.-
Accessory-Set M (SW 43-46 cm) Order. no. 9802929000	229.-
Accessory-Set L (SW 48 cm) Order. no. 9802930000	229.-
Accessory-Set XL (SW 51 cm) Order. no. 9802931000	249.-
Accessory-Set XXL (SW 56 cm) Order. no. 9802932000	249.-
Frame colour	
Silver	●

* For seat width 550 mm these additional telescopic tubes are required.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge