

ORDER FORM

PYRO START
Light weight wheelchair

Company:

Customer:

Date:

Basic price: **Pyro Start** € 1.490,-Basic price: **Pyro Start with hub brake** € 1.570,-

Max. capacity 125 kg

Basic price: **Togo Bundle Pyro Start** € 4.230,-

Pushing and braking aid with Pyro Start fully assembled

Max. capacity 125 kg

▼ Please mark!

Seat	€
Seat width in mm	
370 400 430 460 490	○
520	45.50
Seat height front in mm	
510	●
470 490	50.10
Seat depth 420 mm	●
Seat cushion height 5 cm	
Sitzbreite 370 mm Order no. 9800080000	27.50
Sitzbreite 400 mm Order no. 9800081000	27.50
Sitzbreite 430 mm Order no. 9800082000	27.50
Sitzbreite 460 mm Order no. 9800083000	27.50
Sitzbreite 490 mm Order no. 9800084000	29.70
Sitzbreite 520 mm Order no. 9800085000	38.90
Seat upholstery. Nylon	●
Back support	
Back support height 400 mm	●
Back support upholstery. Vinyl	●
Head support* incl. fixation	215.10
Handgrip elevating. removable* Pair	
without hub brake Order no. 1270806000	174.40
with hub brake Order no. 1270807000	
Sidepanels	
Sidepanel removable	●
with arm support. swing-away. elevating	
Sidepanel elevating left right pcs.	
left Order no. 1352805000	123.90
right Order no. 1352806000	

* Head support and elevating handgrips cannot be mounted together


▼ Please mark!

Leg support	€
Swing-away leg support with individual footplates plastic, removable, adjustable in height	●
Leg support elevating left / right	123.10
left SB 370 – 400 Order no. 1352351001	
left SB 430 – 460 Order no. 1352351002	
left SB 490 – 520 Order no. 1352351003	
right SB 370-400 Order no. 1352350001	
right SB 430-460 Order no. 1352350002	
right SB 490-520 Order no. 1352350003	
Leg support elongation 100 mm piece Order no. 1330333000	34.80
Amputee support left / right	
left Order no. 1352371001	179.60
right Order no. 1352370001	
Footplate with heel strap	●
Wheel locks	
Knee lever wheel lock push to lock	●
Brake lever extension	26.30
Order no. 0009677000	
Knee lever wheel lock for one-arm drive left / right	
left Order no. 9000611000	101.40
right Order no. 9000610000	
Tyres	
Rear wheel 24" x 1 3/8 " PU incl. quick-release axle	●
Front wheel 200 x 50 mm PU	●
PU-Tyres solid rubber	●
Wheel base extension	●
Adapter block	●
continue next page	

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge



ORDER FORM

PYRO START
Light weight wheelchair

Company:

Customer:

Date:

▼ Please mark!

Spoke protection, transparent 24" Order no. 9800054000	165.90
One arm drive only without drum brake left / right left für SB 370 – 400 Order no. 5000017646 right für SB 430 – 520 Order no. 5000017645	649,-
Hand rim rubber-coated Pair Order no. 9800066000	180.80
Accessories	
Anti-tip device Pair Order no. 1270845000	100.20
Tilting aid	●
Lap belt Order no. 9800012000	30.70
Crutch holder Order no. 9800030002	27.10
Therapy tray	
Seat width (S) 370 – 430 Order No. 9800006000	255.-
Seat width (L) 430 – 490 Order No. 9800007000	255.-
Seat width (XL) 490 – 560 Order No. 9800008000	285.-
Therapy tray Hemi Order No. 9800002000	315.-

▼ Please mark!

Passive light	●
Wheelchair accessory set consisting of: Universal cane holder for 2 walking sticks, 1 pc "Press" standard safety belt, 1 pc Seat cushion, black, 1 pc 1 pair of anti-tip wheels	
Accessory-Set S (SW 37-40 cm) Order. no. 9802928000	229.-
Accessory-Set M (SW 43-46 cm) Order. no. 9802929000	229.-
Accessory-Set L (SW 48 cm) Order. no. 9802930000	229.-
Accessory-Set XL (SW 51 cm) Order. no. 9802931000	249.-
Accessory-Set XXL (SW 56 cm) Order. no. 9802932000	249.-
Frame color	
Silver	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge