

ORDER FORM

S-ECO 2

Standard wheelchair

Company:

Basic price: **S-Eco 2 without drum brake**

€ 979,-

Basic price: **S-Eco 2 with drum brake**

€ 1059,-

Customer:

Max. capacity 125 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm	
370 400 430 460 490	○
520	45,50
Seat height 500 mm	●
Seat depth 420 mm	●
Sitzkissen Höhe 5 cm	
Seat width 370 mm Order No. 9800080000	27,50
Seat width 400 mm Order No. 9800081000	27,50
Seat width 430 mm Order No. 9800082000	27,50
Seat width 460 mm Order No. 9800083000	27,50
Seat width 490 mm Order No. 9800084000	29,70
Seat width 520 mm Order No. 9800085000	38,90
Seat upholstery, Nylon	●
Sidepanels	
Duo-Sidepanel long SWing-away	●
Sidepanel elevating left / right	
links Best. Nr. 1352805000	123,90
rechts Best. Nr. 1352806000	
Sidepanel elevating with foldable arm pads left / right	
links Best. Nr. 1352811000	141,30
rechts Best. Nr. Rechts 1352812000	
Leg support	
Swing-away leg support	●
Lower leg support elevating left / right piece	
left SW 370 – 400 Order No. 1352351001	
left SW 430 – 460 Order No. 1352351002	
left SW 490 – 520 Order No. 1352351003	123,10
right SW 370 – 400 Order No. 1352350001	
right SW 430 – 460 Order No. 1352350002	
right SW 490 – 520 Order No. 1352350003	

▼ Please mark!


Amputee support left / right	piece	
left Order No. 1352371001		179,60
right Order No. 1352370001		
Footplate with heel strap		●
Tyres		
Rear wheels 24" x 1 3/8 "		●
Front wheels 7" x 30 mm		●
PU-tyres solid rubber		●
Spoke protection, transparent 24" pair		165,90
Order No. 9800054000		
One arm drive only without drum brake left / right		
left für SW 370 – 400 Order No. 5000017646		649,-
right für SW 430 – 520 Order No. 5000017645		
Hand rim rubber-coated 24"	pair	180,80
Order No. 9800066000		
Back		
Back support height 400 mm		●
Handgrip elevating, removable pair *		
without DB Order No. 1270806000		174,40
with DB Order No. Order No. 1270807000		
Head support incl. fixation *		
für SW 380 – 430 mm Order No. 9802900000		206,80
für SW 450 – 520 mm Order No. 9802901000		
Back support upholstery, Nylon		●
More next page		

* Head support and elevating handgrips cannot be mounted together

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge



ORDER FORM

S-ECO 2
Standard wheelchair

Company:

Customer:

Date:

▼ Please mark!

Wheel locks		
Knee lever wheel lock push to lock		●
Wheel lock lever extension piece Order No. 0009677000	26,30	
Wheel lock for one-arm drive mechanical left / right left Order No. 9000611000 right Order No. 9000610000	101,40	
Accessories		
Togo power push aid incl. mounting set Max. capacity 125 kg Order no. 8010000000	2.770,-	
Mounting aid for Togo adaptation Order No. 8010199000	216,80	
Anti-tip device pair Order No. 1270845000	100,20	
Tilting aid		●
Lap belt Order No. 9800012000	30,70	
Crutch holder Order No. 9800030002	27,10	
Therapy tray		
Seat width (S) 370 – 430 Order No. 9800006000	255.-	
Seat width (L) 430 – 490 Order No. 9800007000	255.-	
Seat width (XL) 490 – 560 Order No. 9800008000	285.-	
Therapy tray Hemi Order No. 9800002000	315.-	

▼ Please mark!

Passive Beleuchtung	●
Wheelchair accessory set consisting of: Universal cane holder for 2 walking sticks, 1 pc "Press" standard safety belt, 1 pc Seat cushion, black, 1 pc 1 pair of anti-tip wheels	
Accessory-Set S (SW 37-40 cm) Order. no. 9802928000	229.-
Accessory-Set M (SW 43-46 cm) Order. no. 9802929000	229.-
Accessory-Set L (SW 48 cm) Order. no. 9802930000	229.-
Accessory-Set XL (SW 51 cm) Order. no. 9802931000	249.-
Accessory-Set XXL (SW 56 cm) Order. no. 9802932000	249.-
Frame colour	
Anthracite metallic	●

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge