

ORDER FORM

EJOY FD XL

Electrically powered wheelchair

Company:

Commission:

Date:

Basic price: **Ejoy FD XL | 6 km/h** € 6.350,-
HMV-No.: 18.50.04.3003 | Max. total capacity 170 kg

Basic price: **Ejoy FD XL | 10 km/h** € 7.200,-
Max. total capacity 170 kg

▼ Please mark (bold = default setting)

Telescopic seating system		€
Seat width mm		
420 mm 470 mm 520 mm		○
effective seat width + 120 mm extendable side panels ¹		
Seat cushion 50 mm high, max user weight 170kg Nylon		●
Contoured seat cushion max. user weight 170kg black		
for seat width		
420 mm / Order No. 8440920000	303,45	
470 mm / Order No. 8440921000		
520 mm / Order No. 8440922000		
Seat depth mm		○
400 425 450 475 500 525 550		
Seat height rear mm without seat cushion		○
430 450 470 490 510 530		
Seat angle 0°– 20° manually adjustable		●
Adaptable back height 500 mm		●
Contour backrest cushion black		
for seat width		
420 mm / Order No. 8440923000	303,45	
470 mm / Order No. 8440924000		
520 mm / Order No. 8440925000		
Back support angle 90°– 120°		●
manually adjustable, 4 steps		
Side panel		
Adjustable in height and lateral*		●
Side panels extendable each side +25 mm		
Legrest		
Swing away legrest with angle adjustable foot plate		●
Continuous footboard, swing away, angle adjustable		
for seat width		
420 mm / Order No. 8440386000	216,70	
470 mm / Order No. 8440387000		
520 mm / Order No. 8440388000		
(Length of footboard has to be the same size as seat width)		

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting)

Manual elevating legrest	left / right	
	unit complete left or right	123,10
for seat width		
420 mm / right (Order No. 1352350001)		
470 mm / right (Order No. 1352350002)		
520 mm / right (Order No. 1352350003)		
420 mm / left (Order No. 1352351001)		
470 mm / left (Order No. 1352351002)		
520 mm / left (Order No. 1352351003)		
Lower leg length: ²		○
370 390 410 430 450 470 490 510 530		
the lower leg length must not be greater than the seat height		
Amputee legrest	left / right	unit 179,60
left / Order No. 1352371001		
right / Order No. 1352370001		
Battery		
Battery 2 x 12 V/56 Ah		●
Battery 2 x 12 V/80 Ah	upcharge per 2 pieces	288,30
Order No. 8530700000		
Charger 24 V/5 A		●
Charger 24 V/8 A ³ upcharge		103,-
Order No. 8130845000		
Lights		
Light system incl. indicators	according to German street regulations	●


²The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.³When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

1 v. 2



ORDER FORM

EJOY FD XL
Electrically powered wheelchair

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

€

Remote		
Swing-away length adjustable remote bracket		●
Control without electrical adjustments		●
Control up to 2 powered adjustments Order No. 8355805000	364,-	
Attendant remote incl. bracket Order No. 8401829001	1.035,30	
Control up to 4 powered adjustments Order No. 8540846000	1.090,-	
Attendant remote incl. bracket Order No. 8401829002	1.423,10	
Seat angle 0°– 20° electrically adjustable Order No. 8530961000	1.402,20	
Back support angle 100°– 142° electrically adjustable ⁴ Order No. 8440960000	1.402,20	
Powered elevating legrest left / right unit left / Order No. 8355361000 right / Order No. 8355360000	870,-	
Remote mounting position left right		○
Attendant remote mounting position left right		○
Tyres		
Drive wheels suspensions		●
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000		●
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70	
Castor 9", PU abrasion resistant, black pair Order No. 8440400000		●
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50	
Anti-tipper		●
Fender front and rear		●
Accessories		

⁴Headrest is mandatory by powered backrest recline 100°–142°.

▼ Please mark (bold = default setting)

€

Lap belt		●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-	
Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80	
Rear view mirror left / right unit incl. bracket left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	140,30	
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-	
Licence plate holder ⁵ Order No. 8166750000	43,60	
Frame colour		
White		●

⁵In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

2 v. 2