

ORDER FORM

EJOY RD XL

Electrically powered wheelchair

Company:

Commission:

Date:

Basic price: **Ejoy RD XL | 6 km/h** € 5.875,-
HMV-No.: 18.50.04.3002 | Max. total capacity 170 kg

Basic price: **Ejoy RD XL | 10 km/h²** € 6.725,-
/ Max. total capacity 170 kg

▼ Please mark (bold = default setting) €

Telescopic seating system

Seat width mm	
420 mm 470 mm 520 mm	○
effective seat width + 120 mm extendable side panels ¹	
Seat cushion 50 mm high, max user weight 170kg Nylon	●
Contoured seat cushion max. user weight 170kg black	
for seat width	
420 mm / Order No. 8440920000	303,45
470 mm / Order No. 8440921000	
520 mm / Order No. 8440922000	
Seat depth mm	
400 425 450 475 500 525 550	○
Seat height rear mm without seat cushion	
430 450 470 490 510 530	○
Seat angle 0° – 20° manually adjustable	●
Adaptable back height 500 mm	●
Contour backrest cushion black	
for seat width	
420 mm / Order No. 8440923000	303,45
470 mm / Order No. 8440924000	
520 mm / Order No. 8440925000	
Back support angle 90°– 120°	●
manually adjustable, 4 steps	
Side panel	
Adjustable in height and lateral*	●
Side panels extendable each side +25 mm	
Legrest	
Swing away legrest with angle adjustable foot plate	●
Continuous footboard, swing away, angle adjustable	
for seat width	
420 mm / Order No. 8440386000	216,70
470 mm / Order No. 8440387000	
520 mm / Order No. 8440388000	
(Length of footboard has to be the same size as seat width)	

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting) €

Manual elevating legrest	left / right	
	unit complete left or right	123,10
for seat width		
420 mm / right (Order No. 1352350001)		
470 mm / right (Order No. 1352350002)		
520 mm / right (Order No. 1352350003)		
420 mm / left (Order No. 1352351001)		
470 mm / left (Order No. 1352351002)		
520 mm / left (Order No. 1352351003)		
Lower leg length: ³		
370 390 410 430 450 470 490 510 530		○
the lower leg length must not be greater than the seat height		
Amputee legrest	left / right	unit 179,60
left / Order No. 1352371001		
right / Order No. 1352370001		
Battery		
Battery 2 x 12 V/56 Ah		●
Battery 2 x 12 V/80 Ah	upcharge for 2 pieces	288,30
Order No. 2 x 8530700000		
Charger 24 V/5 A		●
Charger 24 V/8 A ⁴ upcharge		103,-
Order No. 8130845000		
Lights		
Light system incl. indicators	according to German street regulations	●
continue next page		

²In Germany, electric scooters with a maximum speed of more than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.


⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

1 v. 2



ORDER FORM

EJOY RD XL
Electrically powered wheelchair

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

€

Remote			
Swing-away length adjustable remote bracket			●
Control without electrical adjustments			●
Control up to 2 powered adjustments		364,-	
Order No. 8445810000			
Attendant remote incl. bracket		1.035,30	
Order No. 8401829001			
Control up to 4 powered adjustments		1.090,-	
Order No. 8445812000			
Attendant remote incl. bracket		1.423,10	
Order No. 8401829002			
Seat angle 0°– 20° electrically adjustable		1.402,20	
Order No. 8530961000			
Back support angle 100°– 142° electrically adjustable ⁵		1.402,20	
Order No. 8440960000			
Powered elevating legrest	left / right	unit	870,-
left / Order No. 8355361000			
right / Order No. 8355360000			
Remote mounting position	left	right	○
Attendant remote mounting position	left	right	○
Tyres			
Drive wheels suspensions			●
Drive wheel 14", PU abrasion resistant, black pair			●
Order No. 8440500000			
Drive wheel 14", air abrasion resistant, black pair		294,70	
Order No. 8440505000			
Castor 9", PU abrasion resistant, black pair			●
Order No. 8440400000			
Castor 9", air abrasion resistant, black pair		216,50	
Order No. 8440405000			
Anti-tipper			●
Fender front and rear			●
Accessories			
Lap belt			●

▼ Please mark (bold = default setting)

€

Headrest	multidimensional adjustments incl. bracket			215,-
Order No. 8440917000				
Calf strap				
420 mm / Order No. 8440320000				47,80
470 mm / Order No. 8440321000				
520 mm / Order No. 8440322000				
Rear view mirror	left	right		
			unit incl. bracket	140,30
left / Order No. mirror 8401889000, holder 8401809000				
right / Order No. mirror 8401884000, holder 8401809000				
Cane holder for 2 canes	left	right	unit	64,-
for mechanical back support, Order No. 5000016048				
for electrical back support, Order No. 5000016070				
Licence plate holder ⁶				43,60
Order No. 8166750000				
Frame colour				
White				●

⁶In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

⁵Headrest is mandatory by powered backrest recline 100°–142°.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

2 v. 2