

ORDER FORM

Page 1/4

BX 11
Adaptive wheelchair



BX 11 (with continuous footplate)



BX 11 S (with detachable legrests)

Cost proposal


Order

Invoice address/Customer number	
Cust. no.	
Komission	
Name	
Street	
Postcode/Area	
E-mail	
Tel.	

Shipping address	
Cust. no.	
Komission	
Name	
Street	
Postcode/Area	
E-mail	
Tel.	

The order takes place subject to the selected combination being checked for technical feasibility.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm



ORDER FORM

Page 2/4

BX 11

BX 11 S

Adaptive wheelchair

General information

User weight 135 kg | Material aluminium | Quick-release axle | Depth-adjustable seat sling | Adjustable back upholstery, breathable, black | Clothes guard, aluminium | Knee lever brake (B+B active pressur brake) | Aluminium handrims, black, powder-coated |
With the drum brake option: handrims and rims in silver, PU tyres

Front frame

Base price*

BX 11 | closed front frame with continuous footplate | HMV-No.: 18.50.03.0200

2.100,-

BX 11 S | two-piece front frame for removable legrests | HMV-No.: 18.50.03.0200

2.150,-

Legrests

One piece, with aluminium footplate | Standard with BX 11 | angle-adjustable, flip-up

Two-piece, removable with plastic footplates | Standard with BX 11 S | pair in frame colour | footplates angle-adjustable, can be flipped up

Legrest horizontal adjustable | in black | only with BX 11 S

left

118,40

right

118,40

Amputee support | in black | only with BX 11 S

left

172,70

right

172,70

Frame colours

Black White Blue Red

☐

User-related dimensions

Front seat height (FSH) | in mm | *see matrix below*

490 500 510 520 530 540 550 560

☐Rear seat height (RSH) | in mm | *see matrix below* | with 24" drive wheels, with 22" approx. 25 mm less

420 440 460 480 500 520

☐**Attention: from a difference of 60 mm between FSH and RSH, the angle-adjustable back must be installed.**

Backrest height (RH) | in mm

275 300 325 350 375 400

☐Lower leg length BX 11 (UL) | in mm | *see matrix below* | minimum gap from footplate to floor is 40 mm

430 445 460 475 490 505 520 535

☐

Lower leg length BX 11 S (UL) | in mm | minimum gap from footplate to floor is 40 mm

430 450 470 490 510

☐Seat width (SW) & Seat depth (SD) | in mm | *combination only inside the table possible*

SW	SD	SW	SD	SW	SD
360	380	420	440	480	500
380	400	440	460	500	520
400	420**	460	480**	520	540**

Please select only from one
of the three tables !

☐

* Price plus statutory VAT

** Seat depth not possible with BX 11 S variants

		UL							
FSH ¹		430	445	460	475	490	505	520	535
	490	x	x						
	500	x	x	x	x				
	510	x	x	x	x				
	520	x	x	x	x	x	x		
	530	x	x	x	x	x	x	x	
	540	x	x	x	x	x	x	x	
	550	x	x	x	x	x	x	x	x
	560	x	x	x	x	x	x	x	x


¹Front Seat Height

		Castor Size				
FSH ¹		4"	5"	6"	7"	
	490	x	x			
	500	x	x	x		
	510	x	x	x	x	
	520	x	x	x	x	
	530	x	x	x	x	
	540	x	x	x	x	
	550		x	x	x	
	560			x	x	

		Castor Size				
RSH ²		4"	5"	6"	7"	
	420	x	x	x	x	
	440	x	x	x	x	
	460	x	x	x	x	
	480	x	x	x	x	
	500	x	x	x	x	
	520	x	x	x	x	

With 7" front wheels, only centres of gravity 30 - 54 mm possible.

²Rear Seat Height☒ Basic configuration☐ optional without additional cost



Frame accessories		
Tipping aids	left	78,30
	right	78,30
Anti-tip wheels	pair	225,50
Seat cushion thickness 50mm		
Desired dimension in mm:		
Backrest system and sliding tubes		
Straight backrest, adjustable sliding tubes		●
Angle-adjustable backrest, height-adjustable sliding tubes, removable		255.-
Sideguards		
Aluminium mudguard Splashguard lip pair		●
Detachable, height-adjustable with armrest pair		204.-



One piece legrests, with footplates made of aluminium



Two-piece legrests, removable with footplates, plastic



Legrest horizontal adjustable




Amputee support



Angle-adjustable backrest



Height-adjustable detachable



Aluminium clothes guard (Splashguard lip)



Sideguard, detachable, height-adjustable




Anti-tip wheels



Sopportes pisador

● Basic configuration ○ optional without additional cost

**Front wheels** | Please note the interdependency between the lower leg length and front seat height (matrix page 2).

With 7" front wheels, only centres of gravity 30 - 54 mm possible.

- 4" plastic rim, soft rubber, black
- 5" plastic rim, soft rubber, black
- 6" plastic rim, soft rubber, black
- 6" plastic rim, PU grey
- 7" plastic rim, soft rubber, black

☐

☐

☐

☐

☐

Drive wheels | Please note the interdependency on the rear seat height (matrix page 2)

- 24" PU on aluminium rim | pair
- 22" pneumatic tyres with aluminium rim | pair
- 22" Schwalbe Marathon plus on aluminium rim | pair
- 24" pneumatic tyres with aluminium rim | pair
- 24" Schwalbe Marathon plus on aluminium rim | pair

☐

☐

173,80

☐

173,80

Rear wheel camber

- 0° Standard | 1° 2°

☐

Centre of gravity | in mm | measured from the front edge of the backrest tube, see sketch below (66 and 78 not possible on SW 360 | 380 | 400)

- 30 42 54 | Standard 66 78 | 30 = passive, 78 = active

☐

Fitting handrims

- 20 mm gap | Standard
- 30 mm gap | not possible with drum brake and PU rear wheels

☐

Drive wheel accessories

- Transparent 22" spokeguard | pair
- Transparent 24" spokeguard | pair

159,50

159,50

Brake

- Drum brake option | including drum brake wheels 24" | PU
- Brake lever extension left
- right

316,80


25,30

25,30

Seat accessories

- Lap belt

29,50

Notes

Centre of gravity, side view

● Basic configuration ○ optional without additional cost