



ORDER FORM

ECON XXL
Standard wheelchair

Company:

Customer:

Date:

Basic price: Seat width 600 mm with hub brake	€ 3.369,-
Basic price: Seat width 650 mm with hub brake	€ 3.369,-
Basic price: Seat width 700 mm with hub brake	€ 3.369,-
Basispreis: Sitzbreite 600 mm without hub brake	€ 3.289,-
Basispreis: Sitzbreite 650 mm without hub brake	€ 3.289,-
Basispreis: Sitzbreite 700 mm without hub brake	€ 3.289,-

Max. capacity 250 kg

▼ Please mark!

Seat	€
Seat width in mm 600 650 700	☐
Seat height 475 – 525 mm	☒
Seat depth 480 mm	☒
Seat cushion strengthened	☒
Back support	
Back support height 450 mm w/o seat cushion	☒
Back support adjustable by strap system	☒
Sidepanels	
Sidepanels swing away and height adjustable	☒
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	☒
Footplate elevating	☒
Lower leg support elevating left Stk Order no. 1910361000	140.90
right Stk Order no. 1910360000	
Amputee support left Stk Best. Nr.: 1910365000	179.60
right Stk Best. Nr.: 1910364000	

▼ Please mark!

Tyres	€
Rear wheels 24" PU strengthened axis	☒
Front wheels 7" PU	☒
Wheel locks	
Knee lever wheel lock push to lock	☒
Drum brake	☒
Accessories	
Tilting aid	☒
Anti-tip device Pair Order no. 1910845000	100.10
Crutch holder	☒
Stabilization tube	☒
Back support with bag	☒
Calf strap	☒
Lap belt Order no. 9800019000	45.50
Passive light	☒
Frame colour	
Anthracite	☒

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

☒ Basic equipment ☐ Optional, no extra charge