

# ORDER FORM

## EJOY RD XXL

Electrically powered wheelchair

**Company:**

**Commission:**

**Date:**

Basic price: **Ejoy RD XXL** | 6 km/h € 14.650,-  
 With VR2 Control, incl. 2 electrical adjustments  
 | max. total capacity 200 kg  
 HMV-No.: 18.50.04.5007

Basic price: **Ejoy RD XXL** | 6 km/h € 16.790,-  
 With R-Net Control, incl. 4 electrical adjustments  
 | max. total capacity 250 kg  
 HMV-Nr.: 18.50.04.5007

▼ Please mark (bold = default setting) €

<b>Telescopic seating system</b>	
Seat width mm <b>580 mm</b> 630 mm 680 mm <span style="float: right;">○</span> effective seat width   +120 mm extendable side panels <sup>1</sup>	
Seat cushion 50 mm high, max user weight 250kg   Nylon <span style="float: right;">●</span>	
Contoured seat cushion max. user weight 250 kg   black   for seat width   580 mm / Order No. 8540920000 <span style="float: right;">303,45</span>   630 mm / Order No. 8540921000   680 mm / Order No. 8540922000	
Seat depth mm 400 425 <b>450</b> 475 500 525 550 <span style="float: right;">○</span>	
Seat height rear mm   without seat cushion 430 450 <b>470</b> 490 510 530 <span style="float: right;">○</span>	
Seat angle 0° – 20°   manually adjustable <span style="float: right;">●</span>	
Contour backrest cushion   black   for seat width   580 mm / Order No. 8540923000 <span style="float: right;">303,45</span>   630 mm / Order No. 8540924000   680 mm / Order No. 8540925000	
Back support angle 90°– 120° <span style="float: right;">●</span>   manually adjustable, 4 steps	
<b>Side panel</b>	
Adjustable in height and lateral Side panels extendable each side   +60 mm <span style="float: right;">●</span>	
<b>Legrest</b>	
Swing away legrest with angle adjustable foot plate <span style="float: right;">●</span>	
Continuous footboard, swing away, angle adjustable   for seat width 580 mm / Order No. 8540386000 <span style="float: right;">●</span> 630 mm / Order No. 8540387000 680 mm / Order No. 8540388000 (Length of footboard has to be the same size as seat width)	
Manual elevating legrest   left / right <span style="float: right;">123,10</span>   unit complete left or right   right (Order No. 1352350003)   left (Order No. 1352351003)	

<sup>1</sup>Effective seat width by extendable side panels.

▼ Please mark (bold = default setting) €

Lower leg length: <sup>2</sup> <b>370</b> 390 410 430 450 470 490 510 530 <span style="float: right;">○</span>   the lower leg length must not be greater than the seat height	
Amputee legrest   left / right <span style="float: right;">  unit 179,60</span>   left / Order No. 1352371001   right / Order No. 1352370001	
<b>Battery</b>	
Battery 2 x 12 V/90 Ah <span style="float: right;">●</span>	
Charger 24 V/10 A <span style="float: right;">●</span>	
<b>Lights</b>	
Light system incl. indicators   according to German street regulations <span style="float: right;">●</span>	
<b>Remote</b>	
Swing-away length adjustable remote bracket <span style="float: right;">●</span>	
<b>Electrical adjustments options</b>	
Attendant remote   incl. bracket <span style="float: right;">1.035,30</span>   Bestellnr. 8401829001	
R-net Control with 4 electrical adjustment options <span style="float: right;">1.090,-</span>   included as standard for 250 kg model   Order No. 8445812000	
Attendant remote   incl. bracket <span style="float: right;">1.423,10</span>   Order No. 8401829002	
Seat support angle 0° – 20°   electrically adjustable <span style="float: right;">1.402,20</span>   Order No. 8530961000	
Back support angle 100°– 142°   electrically adjustable <sup>3</sup> <span style="float: right;">1.402,20</span>   Order No. 8440960000	
Powered elevating legrest left / right <span style="float: right;">  unit 870,-</span>   left / Order No. 8355361000   right / Order No. 8355360000	
Remote mounting position left <b>right</b> <span style="float: right;">○</span>	
Attendant remote mounting position left <b>right</b> <span style="float: right;">○</span>	
continue next page	


<sup>2</sup>The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

<sup>3</sup>Headrest is mandatory by powered backrest recline 100°–142°.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment ○ Free of charge option

1 v. 2



# ORDER FORM

**EJOY RD XXL**  
Electrically powered wheelchair

**Company:**

**Commission:**

**Date:**

▼ Please mark (bold = default setting)

€

<b>Tyres</b>	
Drive wheels suspensions	<input checked="" type="radio"/>
Drive wheel 14", PU   abrasion resistant, black   pair   Order No. 8440500000	<input checked="" type="radio"/>
Drive wheel 14", air   abrasion resistant, black   pair   Order No. 8440505000	294,70
Castor 9", PU   abrasion resistant, black   pair   Order No. 8440400000	<input checked="" type="radio"/>
Castor 9", air   abrasion resistant, black   pair   Order No. 8440405000	216,50
Fender   front and rear	<input checked="" type="radio"/>
Anti-tipper	<input checked="" type="radio"/>
<b>Accessories</b>	
Lap belt	<input checked="" type="radio"/>
Headrest   multidimensional adjustments incl. bracket   Order No. 8440917000	215,-
Calf strap   580mm / Order No. 8540320000   630 mm / Order No. 8540321000   680 mm / Order No. 8540322000	47,80
Rear view mirror left / right     unit incl. bracket	140,30
left / Order No. mirror. 8401889000, holder 8401809000   right / Order No. mirror 8401884000, holder 8401809000	
Cane holder for 2 canes left right   unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder   Order No. 8166750000	43,60
<b>Frame colour</b>	
White	<input checked="" type="radio"/>

Standard equipment

Free of charge option

2 v. 2

Version: 2025.03.03 Prices plus VAT at the legally applicable rate, subject to technical changes.