

ORDER FORM

EJOY FD

Electrically powered wheelchair

Company:

Commission:

Date:

Basic price: **Ejoy FD** | **6 km/h** € 5.875,-
HMV-No.: 18.50.04.0223 | Max. total capacity 140 kgBasic price: **Ejoy FD** | **10 km/h**² € 6.665,-
Max. total capacity 140 kg

▼ Please mark (bold = default setting)

Telescopic seating system										€
Seat width mm										
420 mm	470 mm	520 mm								○
effective seat width +50 mm extendable side panels ¹										
Seat cushion 50 mm high, max user weight 140kg Nylon										●
Contured seat cushion max. user weight 140kg black										
for seat width										
420 mm / Order No. 8440920000										303,45
470 mm / Order No. 8440921000										
520 mm / Order No. 8440922000										
Seat depth mm										○
400	425	450	475	500	525	550				
Seat heigth rear mm without seat cushion										○
430	450	470	490	510	530					
Seat angle 0° – 20° manually adjustable										●
Adaptable back height 500 mm										●
Contour backrest cushion black										
for seat width										
420 mm / Order No. 8440923000										303,45
470 mm / Order No. 8440924000										
520 mm / Order No. 8440925000										
Back support angle 90°– 120°										●
manually adjustable, 4 steps										
Side panel										
Adjustable in height and lateral*										●
Side panels extendable each side +25 mm										
Legrest										
Swing away legrest with angle adjustable foot plate										●
Continuous footboard, swing away, angle adjustable										
for seat width										
420 mm / Order No. 8440386000										216,70
470 mm / Order No. 8440387000										
520 mm / Order No. 8440388000										
(Length of footboard has to be the same size as seat width)										

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting)

Manual elevating legrest	left / right	
	unit complete left or right	123,10
for seat width		
420 mm / right (Order No. 1352350001)		
470 mm / right (Order No. 1352350002)		
520 mm / right (Order No. 1352350003)		
420 mm / left (Order No. 1352351001)		
470 mm / left (Order No. 1352351002)		
520 mm / left (Order No. 1352351003)		
Lower leg length: ³		○
370 390 410 430 450 470 490 510 530		
the lower leg length must not be greater than the seat height		
Amputee legrest	left / right	unit 179,60
left / Order No. 1352371001		
right / Order No. 1352370001		
Battery		
Battery 2 x 12 V/55 Ah		●
Battery 2 x 12 V/80 Ah	upcharge per 2 pieces	288,30
Order No. 2 x 8530700000		
Charger 24 V/8 A		●
Lights		
Light system incl. indicators	according to German street regulations	●
Remote		


²In Germany, electric scooters with a maximum speed of more than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

1 v. 2



ORDER FORM

EJOY FD
Electrically powered wheelchair

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

€

Swing-away length adjustable remote bracket	●
Control without electrical adjustments	●
VR2 Control up to 2 powered adjustments Order No. 8355805000	364,-
VR2 Attendant remote incl. bracket Order No. 8401829001	1.035,30
R-net Control up to 4 powered adjustments Order No. 8540846000	1.090,-
R-net Attendant remote incl. bracket Order No. 8401829002	1.423,10
Seat angle 0°–20° electrically adjustable Order No. 8530961000	1.402,20
Back support angle 100°–142° electrically adjustable ⁵ Order No. 8440960000	1.402,20
Powered elevating legrest left / right unit left / Order No. 8355361000 right / Order No. 8355360000	870,-
Remote mounting position left right	○
Attendant remote mounting position left right	○
Tyres	
Drive wheels suspensions	●
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000	●
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70
Castor 9", PU abrasion resistant, black pair Order No. 8440400000	●
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50
Anti-tipper	●
Fender front and rear	●
Accessories	

▼ Please mark (bold = default setting)

€

Lap belt	●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-
Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80
Rear view mirror left / right unit incl. bracket left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	140,30
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder ⁶ Order No. 8166750000	43,60
Frame colour	
White	●

⁶In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

⁵Headrest is mandatory by powered backrest recline 100°–142°.

● Standard equipment ○ Free of charge option

2 v. 2