

ORDER FORM

EJOY RD XL

Electrically powered wheelchair

Company:

Commission:

Date:

Basic price: **Ejoy RD XL | 6 km/h** € 5.875,-
HMV-No.: 18.50.04.3002 | Max. total capacity 170 kg

Basic price: **Ejoy RD XL | 10 km/h²** € 6.725,-
/ Max. total capacity 170 kg

▼ Please mark (bold = default setting)

€

Telescopic seating system

Seat width mm									
420 mm	470 mm	520 mm							○
effective seat width + 50 mm extendable side panels ¹									
Seat cushion 50 mm high, max user weight 170kg Nylon									●
Contoured seat cushion max. user weight 170kg black									
for seat width									
420 mm / Order No. 8440920000									303,45
470 mm / Order No. 8440921000									
520 mm / Order No. 8440922000									
Seat depth mm									○
400	425	450	475	500	525	550			
Seat height rear mm without seat cushion									○
430	450	470	490	510	530				
Seat angle 0° – 20° manually adjustable									●
Adaptable back height 500 mm									●
Contour backrest cushion black									
for seat width									
420 mm / Order No. 8440923000									303,45
470 mm / Order No. 8440924000									
520 mm / Order No. 8440925000									
Back support angle 90° – 120°									●
manually adjustable, 4 steps									
Side panel									
Adjustable in height and lateral*									●
Side panels extendable each side +25 mm									
Legrest									
Swing away legrest with angle adjustable foot plate									●
Continuous footboard, swing away, angle adjustable									
for seat width									
420 mm / Order No. 8440386000									216,70
470 mm / Order No. 8440387000									
520 mm / Order No. 8440388000									
(Length of footboard has to be the same size as seat width)									

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting)

€

Manual elevating legrest	left / right		
	unit complete left or right		123,10
for seat width			
420 mm / right (Order No. 1352350001)			
470 mm / right (Order No. 1352350002)			
520 mm / right (Order No. 1352350003)			
420 mm / left (Order No. 1352351001)			
470 mm / left (Order No. 1352351002)			
520 mm / left (Order No. 1352351003)			
Lower leg length: ³			
370	390	410	430 450 470 490 510 530
the lower leg length must not be greater than the seat height			
Amputee legrest	left / right	unit	179,60
left / Order No. 1352371001			
right / Order No. 1352370001			
Battery			
Battery 2 x 12 V/55 Ah			●
Battery 2 x 12 V/80 Ah			288,30
upcharge for 2 pieces			
Order No. 2 x 8530700000			
Charger 24 V/8 A			●
Lights			
Light system incl. indicators		according to German street regulations	●
continue next page			
Remote			


²In Germany, electric scooters with a maximum speed of more than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

1 v. 2



ORDER FORM

EJOY RD XL Electrically powered wheelchair

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

€

Swing-away length adjustable remote bracket	●
Control without electrical adjustments	●
Control up to 2 powered adjustments Order No. 8445810000	364,-
Attendant remote incl. bracket Order No. 8401829001	1.035,30
Control up to 4 powered adjustments Order No. 8445812000	1.090,-
Attendant remote incl. bracket Order No. 8401829002	1.423,10
Seat angle 0°– 20° electrically adjustable Order No. 8530961000	1.402,20
Back support angle 100°– 142° electrically adjustable ⁵ Order No. 8440960000	1.402,20
Powered elevating legrest left / right unit left / Order No. 8355361000 right / Order No. 8355360000	870,-
Remote mounting position left right	○
Attendant remote mounting position left right	○
Tyres	
Drive wheels suspensions	●
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000	●
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70
Castor 9", PU abrasion resistant, black pair Order No. 8440400000	●
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50
Anti-tipper	●
Fender front and rear	●
Accessories	
Lap belt	●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-

⁵Headrest is mandatory by powered backrest recline 100°–142°.

▼ Please mark (bold = default setting)

€

Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80
Rear view mirror left / right unit incl. bracket left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	140,30
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder ⁶ Order No. 8166750000	43,60
Frame colour	
White	●

⁶In Germany, electric scooters with a maximum speed of than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

● Standard equipment

○ Free of charge option

2 v. 2