

ORDER FORM

EJOY RD

Electrically powered wheelchair

Company:

Commission:

Date:

Basic price: **Ejoy RD | 6 km/h** € 5.460,-
HMV-No.: 18.50.04.0222 | Max. total capacity 140 kg

Basic price: **Ejoy RD | 10 km/h²** € 6.310,-
/ Max. total capacity 140 kg

▼ Please mark (bold = default setting)

€

Telescopic seating system

Seat width mm	
420 mm / Standard	
470 mm / Order No. seat 8440911000, backrest 8440931000	☐
520 mm / Order No. seat 8440912000, backrest 8440932000	
effective seat width +50 mm extendable side panels ¹	
Seat cushion 50 mm high, max user weight 140kg Nylon	☒
Contoured seat cushion max. user weight 140kg black	
for seat width	
420 mm / Order No. 8440920000	303,45
470 mm / Order No. 8440921000	
520 mm / Order No. 8440922000	
Seat depth mm	
400 425 450 475 500 525 550	☐
Seat height rear mm without seat cushion	
430 450 470 490 510 530	☐
Seat angle 0° – 20° manually adjustable	☒
Contour backrest cushion black	
for seat width	
420 mm / Order No. 8440923000	303,45
470 mm / Order No. 8440924000	
520 mm / Order No. 8440925000	
Back support angle 90° – 120°	☒
manually adjustable, 4 steps	
Side panel	
Adjustable in height and lateral	
Side panels extendable each side +25 mm	☒
Legrest	
Swing away legrest with angle adjustable foot plate	☒
Continuous footboard, swing away, angle adjustable	
for seat width	
420 mm / Order No. 8440386000	216,70
470 mm / Order No. 8440387000	
520 mm / Order No. 8440388000	
(Length of footboard has to be the same size as seat width)	

¹Effective seat width by extendable side panels.

▼ Please mark (bold = default setting)

€

Manual elevating legrest	left / right	
	unit complete left or right	123,10
for seat width		
420 mm / right (Order No. 1352350001)		
470 mm / right (Order No. 1352350002)		
520 mm / right (Order No. 1352350003)		
420 mm / left (Order No. 1352351001)		
470 mm / left (Order No. 1352351002)		
520 mm / left (Order No. 1352351003)		
Lower leg length: ³		
370 390 410 430 450 470 490 510 530		☐
the lower leg length must not be greater than the seat height		
Amputee legrest	left / right	unit 179,60
left / Order No. 1352371001		
right / Order No. 1352370001		
Battery		
Battery 2 x 12 V/55 Ah		☒
Battery 2 x 12 V/80 Ah	upcharge for 2 pieces	288,30
Order No. 2 x 8530700000		
Charger 24 V/8 A		☒
Lights		
Light system incl. indicators	according to German street regulations	☒
Remote		

²In Germany, electric scooters with a maximum speed of more than 6 km/h require their own liability insurance. An insurance licence plate must be affixed to the rear of the scooter as proof of this.


³The lower leg length must not be greater than the seat height to ensure sufficient ground clearance.

⁴When using 80 Ah batteries, the use of the 8 A charger is recommended to avoid very long charging times.

The dimensions given are subject to manufacturing tolerances of +/- 10 mm.

☒ Standard equipment☐ Free of charge option

1 v. 2



ORDER FORM

EJOY RD
Electrically powered wheelchair

Company:

Commission:

Date:

▼ Please mark (bold = default setting)

€

Swing-away length adjustable remote bracket	●
Control without electrical adjustments	●
VR2 Control up to 2 powered adjustments Order No. 8442810000	364,-
VR2 Attendant remote incl. bracket Order No. 8401829001	1.035,30
R-netControl up to 4 powered adjustments Order No. 8445812000	1.090,-
R-net Attendant remote incl. bracket Order No. 8401829002	1.423,10
Seat angle 0°– 20° electrically adjustable Order No. 8530961000	1.402,20
Back support angle 100°– 142° electrically adjustable ⁵ Order No. 8440960000	1.402,20
Powered elevating legrest left / right unit left / Order No. 8355361000 right / Order No. 8355360000	870,-
Remote mounting position left right	○
Attendant remote mounting position left right	○
Tyres	
Drive wheels suspensions	●
Drive wheel 14", PU abrasion resistant, black pair Order No. 8440500000	●
Drive wheel 14", air abrasion resistant, black pair Order No. 8440505000	294,70
Castor 9", PU abrasion resistant, black pair Order No. 8440400000	●
Castor 9", air abrasion resistant, black pair Order No. 8440405000	216,50
Anti-tipper	●
Fender front and rear	●
Accessories	

▼ Please mark (bold = default setting)

€

Lap belt	●
Headrest multidimensional adjustments incl. bracket Order No. 8440917000	215,-
Calf strap 420 mm / Order No. 8440320000 470 mm / Order No. 8440321000 520 mm / Order No. 8440322000	47,80
Rear view mirror left / right unit incl. bracket left / Order No. mirror 8401889000, holder 8401809000 right / Order No. mirror 8401884000, holder 8401809000	140,30
Cane holder for 2 canes left right unit for mechanical back support, Order No. 5000016048 for electrical back support, Order No. 5000016070	64,-
Licence plate holder Order No. 8166750000	43,60
Frame colour	
White	●

⁵Headrest is mandatory by powered backrest recline 100°–142°.