



ORDER FORM

TOGO
Power pushing aid

Company:

Customer:

Date:

Basic price: **Togo** € 2.780,-
Max. total capacity 180 kg
(including weight of pushing aid and wheelchair)

▼ Please mark!

Technical data	€
Togo for seat width 37 - 52 cm	●
Speed forward 1,0 up to 5,5 km/h 6 levels	●
Speed backwards approx. 1,5 km/h	●
Max. range approx. 12 km	●
Upward slope / downward slope max. 10%	●
Battery 24 V 11,5 Ah Standard	●
Operating voltage: 24 VDC	●
Battery charger	●
Drive wheels 200 x 50 solid rubber tires	●
Wheel lock electromagnetic	●

▼ Please mark!

Standard features	€
Mounting aid for Togo wheelchair mounting only with tool holder set Item No. 8011122000	0,-
Tool holder for manual wheelchairs Item No. 8011123000	0,-
Tool holder for multifunctional wheelchair Protego Item No. 8011122400	0,-
Mounting aid for adaptation Item No. 8011199000	219,-
Frame color	
Silver	●
Accessories	
Anti tipper (pair) for all models Item No. 1270845000	102,21
(for every wheelchair equipped with electric wheelchair push aid Togo, we strongly recommend using anti tipper.	

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment ○ Optional, no extra charge