



ORDER FORM

CENTURO S2
4-WHEEL SCOOTER

Company:

Basic price: **Centuro S2** | 6 km/h

€ 2.904,-

Basic price: **Centuro S2** | 10 km/h

€ 3.507,-

Customer:

Max. payload 160 kg

Date:

▼ Please Mark !

Technical Data		€
Range approx. 35 km	●	
Ground clearance 80 mm	●	
Obstacle height 50 mm	●	
Climbing ability 6° 10%	●	
Reversing width 2050 mm	●	
Tyres 10" pneumatic	●	
Sprung suspension front and rear	●	
Brakes electromagnetic	●	
Mechanical parking brake	●	
Handlebars ergonomic	●	
Adjustable steering column	●	
Motor and electrical equipment		
Motor 6 km/h, 550 W	●	
Motor 10 km/h, 600 W	●	
Battery 2 x 12 V/55 Ah	●	
Battery charger 8 A	●	
USB-Port	●	
Cockpit with LED display	●	
Horn	●	
Lighting according to German road traffic regulations	●	
Authorized for public transportation	●	
Electronic speed reduction when cornering (only for 10 km/h)	●	
Seat and back support		
Swing-away seat, height-adjustable - Back support foldable and angle adjustable - Arm support foldable and angle adjustable - Seat rails for adjusting the distance to the steering column	●	

▼ Please Mark !

Accessories		€
Shopping basket removable	●	
Floormat	●	
Rear mirror right	●	
Rear mirror holder right	●	
Rear mirror left*		20,50
Order No. 8162813000		
Rear mirror holder left		3,-
Order No. 8162815000		
Crutch holder for 2 Crutches** Mounting at the seat frame		79,00
Order No. 8162881000		
License plate holder required for > 6 km/h ***		27,70
Order No. 8166730000		
Rain cover		105,70
Order No. 8162886000		
Frame color		
White	Anthracite	○

*just mountable with rear mirror holder

** The overall length of the scooter is extended by attaching the crutch holder.

*** In Germany electric vehicles with a maximum speed of more than 6 km / h have to own a liability insurance . To prove this you have to put an insurance license plate on the back of the scooter.

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

● Basic equipment

○ Optional, no extra charge