




Pos.	Artikelnummer	Bezeichnung	Description
1	4004000001	Korb	Basket
2	4004000003	Verbindungshebel Rahmenteile	Frame connection lever
3	4004000004	Abdeckkappe Castorbuchse	Plastic cap
4	4004000005	Führungsprofil Faltmechanismus	Sliding profile for folding system
5	4004000006	Sitzplatte gepolstert	Seat plate padded
6	4004000007	Tablett	Tablet
7	4004000008	Faltsicherung	Lock system
8	4004000009	Radgabel R/L	Front fork R/L
9	4004000010	Lenklager Castorbuchse	Steering bearing
10	4004000011	Rad 200x50mm	Wheel 200x50mm
11	4004000025	Griff-/Schiebestange links	Pushing tube left
12	4004000026	Griff-/Schiebestange rechts	Pushing tube right
13	4004000027	Handgriff anatomisch rechts	Hand grip right
14	4004000028	Handgriff anatomisch links	Hand grip left
15	4004000030	Brems- und Bowdenzug	Brake- and Bowden cable
16	4004000031	Schraube Griffrohrbefestigung	Screw for pushing tube
17	4004000032	Vorderradlager	Front fork bearing
18	4004000033	Achsschrauben chrom 8x80mm	Axle screw 8x80mm
19	4004000034	Bremse rechts komplett	Brake complete right
20	4004000035	Bremse links komplett	Brake complete left
21	4004000036	Bremshebel rechts komplett	Brake lever right
22	4004000037	Bremshebel links komplett	Brake lever left
23	4004000039	Stockhalter	Cane holder
24	4004000002	Auflageprofil für Korb	Protection profile for basket
o. Abb.	4004000012	Ankipphilfe rechts	Tilt assist right
	4004000013	Reflektoren	Reflectors
	4004000029	Sitzplatte Kunststoff	Seat plate plastic
	4004000042	Bremsklotz komplett	Brake block complete

